
DOC272.97.80490

SIP 10 for TU5200

10/2015, Edition 1

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[Maintenance](#) on page 10

[Troubleshooting](#) on page 16

[Replacement parts and accessories](#) on page 16

Specifications

Specifications are subject to change without notice.

Specification	Details
Dimensions (W x D x H)	Sipper unit: 120 x 85 x 200 mm (4.7 x 3.4 x 7.9 in.)
Enclosure	IP30
Weight	Sipper unit: 0.5 kg (1.1 lb)
Protection class	II
Pollution degree	2
Installation category	II
Power supply	Power is supplied by the turbidimeter through the USB cable, 530 mA, 5 V
Operating temperature	10 to 40 °C (50 to 104 °F)
Storage temperature	-10 to 60 °C (14 to 140 °F)
Humidity	80% relative humidity, non-condensing
Rinsing agent	120 mL minimum to remove sample from the tubing
Flow speed	1 mL/sec
Sample temperature	2 to 70 °C (35.6 to 158 °F)
Sample salt content	65 g/L maximum
Interface	USB
Certifications	CE
Warranty	1 year (EU: 2 years)

General information

In no event will the manufacturer be liable for direct, indirect, special, incidental or consequential damages resulting from any defect or omission in this manual. The manufacturer reserves the right to make changes in this manual and the products it describes at any time, without notice or obligation. Revised editions are found on the manufacturer's website.

Safety information

NOTICE

The manufacturer is not responsible for any damages due to misapplication or misuse of this product including, without limitation, direct, incidental and consequential damages, and disclaims such damages to the full extent permitted under applicable law. The user is solely responsible to identify critical application risks and install appropriate mechanisms to protect processes during a possible equipment malfunction.

Please read this entire manual before unpacking, setting up or operating this equipment. Pay attention to all danger and caution statements. Failure to do so could result in serious injury to the operator or damage to the equipment.

Make sure that the protection provided by this equipment is not impaired. Do not use or install this equipment in any manner other than that specified in this manual.

Use of hazard information

▲ DANGER

Indicates a potentially or imminently hazardous situation which, if not avoided, will result in death or serious injury.

▲ WARNING

Indicates a potentially or imminently hazardous situation which, if not avoided, could result in death or serious injury.

▲ CAUTION

Indicates a potentially hazardous situation that may result in minor or moderate injury.

NOTICE

Indicates a situation which, if not avoided, may cause damage to the instrument. Information that requires special emphasis.

Precautionary labels


Read all labels and tags attached to the instrument. Personal injury or damage to the instrument could occur if not observed. A symbol on the instrument is referenced in the manual with a precautionary statement.

	This symbol, if noted on the instrument, references the instruction manual for operation and/or safety information.
	Electrical equipment marked with this symbol may not be disposed of in European domestic or public disposal systems. Return old or end-of-life equipment to the manufacturer for disposal at no charge to the user.

Product overview

The SIP 10 sipper kit for the TU5200 turbidimeter is used to pull a selected quantity of a water sample through the flow-thru head of the TU5200 before each turbidity measurement. Refer to [Figure 1](#). If selected, the sipper pulls a selected quantity of a user-supplied rinse agent through the flow-thru head after each turbidity measurement.

Figure 1 Product overview


1 Drain	4 Cleaning lid	7 TU5200 turbidimeter
2 Peristaltic pump	5 Flow-thru head	8 SIP 10 sipper unit
3 USB cable	6 Sample source	

Product components

Make sure that all components have been received. Refer to [Figure 2](#). If any items are missing or damaged, contact the manufacturer or a sales representative immediately.

Figure 2 Product components

1 SIP 10 sipper unit	4 Flow-thru head
2 Tygon® tubing, 1.70 m (5.6 ft)	5 Vial replacement tool ¹
3 USB cable, type AB, 1 m (3.3 ft)	

Installation

Install the flow-thru head





Pre-requisite: Clean the process vial. Refer to [Clean the vial](#) on page 10.

NOTICE

Do not touch or scratch the glass of the process vial. Contamination or scratches on the glass can cause measurement errors.


Install the flow-thru head as shown in the illustrated steps that follow.

¹ Remove before installation.


Install the tubing

Install the tubing as shown in the illustrated steps that follow.

Connect the USB cable

Connect the USB cable as shown in the illustrated steps that follow.

Operation

Configure the settings

1. Connect the turbidimeter power cable to an electrical outlet.
2. Push the power button to set the turbidimeter to on.
3. Push **Sipper**.

4. Select **On** to enable the sipper operation.

Note: When **On** is selected, the sipper unit completes a sipper cycle before each measurement and a purge cycle (if enabled) after each measurement.

5. Select and configure each option.

Option	Description
Sip Time	Sets the period of time that the sipper unit pulls the sample through the process vial before each measurement. Options: 1 to 300 seconds (default: 120 seconds). The pump flow speed is 1 mL/sec.
Settle Time	Sets the period of time that the sample sits in the process vial before it is measured. Options: 1 to 300 seconds (default: 60 seconds). Use the Settle Time setting to let the air bubbles and turbulence in the sample decrease before each measurement.
Purge Time	Sets the period of time that the sipper unit pulls a user-supplied rinse agent through the process vial when a purge cycle is done. Options: Off (default) or 1 to 300 seconds (default: 120 seconds). Note: A minimum of 120 mL of rinse agent is necessary to remove sample from the tubing.
Purge Start	Sets the purge cycle to start automatically after a measurement or when the user pushes Purge. Options: Auto (default) or Manual. Note: The Purge Start button is only enabled when Purge Time is set to on.

6. Push **OK**.

To exit and not save the changes, push **Cancel**. To change the settings to the factory default settings, push **Default**.

Flush the vial and tubing

Before initial use of a new process vial or tubing, flush the process vial and tubing as follows:

1. Put the inlet tubing in approximately 400 mL of deionized water.
Make sure that the tube opening is at the bottom of the container.
2. Push **Sipper**, then push **Purge** to start a purge cycle.
3. Do step 2 again until approximately 360 mL of deionized water is pulled through the process vial and tubing.
4. Push **OK**.

Measure a sample

1. Do a purge as follows:
 - a. Put the inlet tubing in deionized water.
 - b. Push **Sipper**, then push **Purge** to start a purge cycle.
 - c. Push **OK**.
2. Put the inlet tubing in the sample.
3. Push **Read**. The sipper unit completes a sipper cycle, then the turbidimeter measures the sample.
If the Purge Start setting is set to Auto, a purge cycle is completed after the measurement.
If the Purge Start setting is set to Manual, the Read button changes to a Purge button.
4. Push **Purge** to start a purge cycle if applicable.
5. Do steps 3 and 4 again until the sample measurements are complete.
6. Do a purge. Refer to step 1.

Prepare the flow-thru head for storage

Before short-term storage, flush the process vial and tubing. Refer to [Flush the vial and tubing](#) on page 9.

Before long-term storage, do the steps that follow:

1. Remove the process vial from the flow-thru head. Refer to the illustrated steps in [Replace the vial](#) on page 12.
2. Fill the process vial three times with sterilized water. Discard the water.
3. Fill the process vial with water that is deionized and sterilized.
4. Put a plug in the process vial. Use the plug for a sample vial.

Maintenance

Maintenance schedule

[Table 1](#) shows the recommended schedule of maintenance tasks. Facility requirements and operating conditions may increase the frequency of some tasks.

Table 1 Maintenance schedule

Task	1 year	As necessary
Clean the instrument on page 10		X
Clean the vial on page 10		X
Replace the vial on page 12	X	
Replace the pump tubing on page 14		X

Clean the instrument

Clean the exterior of the instrument with a moist cloth and a mild soap solution and then wipe the instrument dry as necessary.

Clean the vial

▲ WARNING	
	Personal injury hazard. Remove power from the instrument before this procedure is started.
NOTICE	
Do not touch or scratch the glass of the process vial. Contamination or scratches on the glass can cause measurement errors.	

Use standard laboratory procedures to clean the process vial before initial use and as necessary to remove contamination from the glass.

As an alternative, clean the process vial with the optional vial wiper. Refer to [Clean the vial with the optional wiper](#) on page 11.

1. Remove the process vial from the flow-thru head. Refer to the illustrated steps in [Replace the vial](#) on page 12.
2. Use standard laboratory procedures to clean the process vial.
3. Install the process vial in the flow-thru head. Refer to the illustrated steps in [Replace the vial](#) on page 12.
4. Clean the exterior of the process vial with a no-lint cloth to remove dirt, fingerprints or particles from the glass.



Clean the vial with the optional wiper



NOTICE

Carefully put the vial wiper into the process vial so that no water spills out.

Clean the process vial with the optional vial wiper and deionized water as shown in the illustrated steps that follow.

Clean the exterior of the process vial with a no-lint cloth to remove dirt, fingerprints or particles from the glass.

Replace the vial

Pre-requisite: Clean the new process vial. Refer to [Clean the vial](#) on page 10.

⚠ WARNING





Personal injury hazard. Remove power from the instrument before this procedure is started.



NOTICE

Do not touch or scratch the glass of the process vial. Contamination or scratches on the glass can cause measurement errors.

Replace the process vial when it has many scratches or after 1 year. Refer to the illustrated steps that follow. Make sure that no water or particles go into the process vial compartment.

When the illustrated steps are complete, flush the process vial and tubing. Refer to [Flush the vial and tubing](#) on page 9.





Replace the pump tubing

When damage or wear is seen on the pump tubing, replace the pump tubing as shown in the illustrated steps that follow.

Item to collect: Pump tubing (Lagoprene®) pre-assembled with peristaltic pump cover and connections

After the tubing is replaced, flush the process vial and tubing. Refer to [Flush the vial and tubing](#) on page 9.

Troubleshooting

Error	Description	Solution
Sipper module not connected. Please check cable.	There is no USB connection between the sipper and the turbidimeter.	Examine the USB cable. Make sure that the cable length is not more than 1 m (3.3 ft). Make sure that the USB cable is connected to the turbidimeter and the sipper.
Please check sipper and tube.	The pump tubing is not correctly installed.	Loosen the peristaltic pump cover, then install it again. Refer to the illustrated steps in Replace the pump tubing on page 14. Operate the pump for 15 seconds. Push down the lever and make sure that the tube is correctly installed around the rollers.

Replacement parts and accessories

⚠ WARNING

Personal injury hazard. Use of non-approved parts may cause personal injury, damage to the instrument or equipment malfunction. The replacement parts in this section are approved by the manufacturer.

Note: Product and Article numbers may vary for some selling regions. Contact the appropriate distributor or refer to the company website for contact information.

Replacement parts

Description	Item no.
Cable, USB type AB, 1 m (3.3 ft)	LZQ104
Cleaning lid	LZY919
Cleaning lid, screws (3x)	LZY921
Flow-thru head, TU5200	LZV967
Seal, cleaning lid	LZY914
Seal, process vial	LZY918
SIP 10 sipper unit with LZV940	LQV157.99.40002
Tubing, pump (Lagoprene) pre-assembled with peristaltic pump cover and connections	LZV877
Tubing kit, includes Tygon tubing 1.70 m (5.6 ft), pump tubing (Lagoprene) pre-assembled with peristaltic pump cover and connections	LZV940
Vial, process	LZY834
Vial nut	LZY917
Vial replacement tool	LZY906

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P.O. Box 389, Loveland, CO 80539-0389 U.S.A.
(800) 227-4224 (U.S.A. only)
HACH LANGE Srl