
DOC272.97.80548

# SIP 10 kit for TL23xx

02/2017, Edition 1

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## Specifications

Specifications are subject to change without notice.

Specification	Details
Dimensions (W x D x H)	Sipper unit: 120 x 85 x 200 mm (4.7 x 3.4 x 7.9 in.)
Enclosure	IP30
Weight	Sipper unit: 0.5 kg (1.1 lb)
Protection class	II
Pollution degree	2
Installation category	I
Power supply	Power is supplied by the turbidimeter through the USB cable, 530 mA, 5 V
Operating temperature	10 to 40 °C (50 to 104 °F)
Storage temperature	-10 to 60 °C (14 to 140 °F)
Humidity	80% relative humidity, non-condensing
Rinsing agent	120 mL minimum to remove sample from the tubing
Flow speed	1 mL/sec
Sample temperature	2 to 70 °C (35.6 to 158 °F)
Sample salt content	65 g/L maximum
Interface	USB
Certifications	CE
Warranty	1 year (EU: 2 years)

## General information

In no event will the manufacturer be liable for direct, indirect, special, incidental or consequential damages resulting from any defect or omission in this manual. The manufacturer reserves the right to make changes in this manual and the products it describes at any time, without notice or obligation. Revised editions are found on the manufacturer's website.

## Safety information

### NOTICE

The manufacturer is not responsible for any damages due to misapplication or misuse of this product including, without limitation, direct, incidental and consequential damages, and disclaims such damages to the full extent permitted under applicable law. The user is solely responsible to identify critical application risks and install appropriate mechanisms to protect processes during a possible equipment malfunction.

Please read this entire manual before unpacking, setting up or operating this equipment. Pay attention to all danger and caution statements. Failure to do so could result in serious injury to the operator or damage to the equipment.

Make sure that the protection provided by this equipment is not impaired. Do not use or install this equipment in any manner other than that specified in this manual.

## Use of hazard information

### ▲ DANGER

Indicates a potentially or imminently hazardous situation which, if not avoided, will result in death or serious injury.

### ▲ WARNING

Indicates a potentially or imminently hazardous situation which, if not avoided, could result in death or serious injury.

### ▲ CAUTION

Indicates a potentially hazardous situation that may result in minor or moderate injury.

### NOTICE

Indicates a situation which, if not avoided, may cause damage to the instrument. Information that requires special emphasis.

## Precautionary labels

Read all labels and tags attached to the instrument. Personal injury or damage to the instrument could occur if not observed. A symbol on the instrument is referenced in the manual with a precautionary statement.

	This symbol, if noted on the instrument, references the instruction manual for operation and/or safety information.
	Electrical equipment marked with this symbol may not be disposed of in European domestic or public disposal systems. Return old or end-of-life equipment to the manufacturer for disposal at no charge to the user.

## Certification

### Canadian Radio Interference-Causing Equipment Regulation, IECS-003, Class A:

Supporting test records reside with the manufacturer.

This Class A digital apparatus meets all requirements of the Canadian Interference-Causing Equipment Regulations.

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### FCC Part 15, Class "A" Limits

Supporting test records reside with the manufacturer. The device complies with Part 15 of the FCC Rules. Operation is subject to the following conditions:

1. The equipment may not cause harmful interference.
2. The equipment must accept any interference received, including interference that may cause undesired operation.

Changes or modifications to this equipment not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses and can

radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference, in which case the user will be required to correct the interference at their expense. The following techniques can be used to reduce interference problems:

1. Disconnect the equipment from its power source to verify that it is or is not the source of the interference.
2. If the equipment is connected to the same outlet as the device experiencing interference, connect the equipment to a different outlet.
3. Move the equipment away from the device receiving the interference.
4. Reposition the receiving antenna for the device receiving the interference.
5. Try combinations of the above.

## Product overview

### **⚠ WARNING**

Fire hazard. This product is not designed for use with flammable liquids.


### **⚠ CAUTION**

Chemical exposure hazard. Dispose of chemicals and wastes in accordance with local, regional and national regulations.

The SIP 10 sipper kit for the TL23xx turbidimeter is used to pull a selected quantity of a water sample through the flow cell in the TL23xx before each turbidity measurement. Refer to [Figure 1](#). If selected, the sipper pulls a selected quantity of a user-supplied rinse agent through the flow cell after each turbidity measurement.

**Figure 1 Product overview**


<b>1</b> USB cable	<b>4</b> Flow cell	<b>7</b> Sample source
<b>2</b> TL23xx turbidimeter	<b>5</b> Flow cell assembly	<b>8</b> SIP 10 sipper unit
<b>3</b> Flow cell cover	<b>6</b> Drain	<b>9</b> Peristaltic pump

## Product components

Make sure that all components have been received. Refer to [Figure 2](#). If any items are missing or damaged, contact the manufacturer or a sales representative immediately.

**Figure 2 Product components**

1 USB cable, type AB, 1 m (3.3 ft)	7 Flow cells (2x)
2 Tygon® tubing, 1/8-in. ID, 1.70 m (5.6 ft)	8 Bottom cap with O-ring
3 Flow cell cover	9 Y fitting for 1/8-in. ID tubing
4 Flow cell assembly	10 Flow cell frame
5 Top cap with O-ring	11 SIP 10 sipper unit
6 Replacement O-rings (2x)	

## Installation

### Clean and install the flow cell assembly

#### NOTICE

Do not touch or scratch the flow cell. Contamination or scratches on the glass can cause measurement errors.

1. Clean the flow cell assembly. Do steps 2 to 4 of [Clean the flow cell assembly](#) on page 12.
2. Install the flow cell assembly as shown in the illustrated steps that follow.

**Note:** Keep the tubing as short as possible to minimize air locking and lag time of the sample flow.





At step 4 of the illustrated steps, apply silicone oil to the flow cell as follows:

- a. Apply a small drop of silicone oil from the top to the bottom of the flow cell.


**Note:** Use only the provided silicone oil. This silicone oil has the same refractive index as the flow cell glass and masks minor glass scratches.

- b. Use the oiling cloth to apply the oil equally to the surface of the flow cell. Remove most of the oil. Make sure that the vial is almost dry.


**Note:** Put the oiling cloth in a plastic storage bag to keep the cloth clean.


**5**


**6**


**7**

## Connect the USB cable

Connect the USB cable. Refer to the illustration that follows.

## Operation

### Configure the settings

1. Connect the turbidimeter power cable to an electrical outlet.
2. Push the power button to set the turbidimeter to on.
3. Push **Sipper**.

**4.** Select **On** to enable the sipper operation.

**Note:** When **On** is selected, the sipper unit completes a sipper cycle before each measurement and a purge cycle (if enabled) after each measurement.

**5.** Select and configure each option.

Option	Description
<b>Sip Time</b>	Sets the period of time that the sipper unit pulls the sample through the flow cell before each measurement. Options: 1 to 99 seconds (default: 8 seconds). The pump flow speed is 1 mL/sec.
<b>Settle Time</b>	Sets the period of time that the sample sits in the flow cell before it is measured. Options: 1 to 999 seconds (default: 5 seconds). Use the Settle Time setting to let the air bubbles and turbulence in the sample decrease before each measurement.
<b>Purge Time</b>	Sets the period of time that the sipper unit pulls a user-supplied rinse agent through the flow cell when a purge cycle is done. Options: Off (default) or 1 to 99 seconds (default: 8 seconds). <b>Note:</b> A minimum of 120 mL of rinse agent is necessary to remove sample from the tubing.
<b>Purge Start</b>	Sets the purge cycle to start automatically after a measurement or when the user pushes Purge. Options: Auto (default) or Manual. <b>Note:</b> The Purge Start button is only enabled when Purge Time is set to on.

**6.** Push **OK**.

To exit and not save the changes, push **Cancel**. To change the settings to the factory default settings, push **Default**.

## Flush the flow cell assembly and tubing

Before initial use of a new flow cell assembly or tubing, flush the flow cell assembly and tubing.

**1.** Put the inlet tubing in approximately 400 mL of deionized water.

Make sure that the tube opening is at the bottom of the container.

**2.** Push **Sipper**, then push **Purge** to start a purge cycle.

**3.** Do step **2** again until approximately 360 mL of deionized water is pulled through the flow cell assembly and tubing.

**4.** Push **OK**.

## Measure a sample

**1.** Do a purge as follows:

a. Put the inlet tubing in deionized water.

b. Push **Sipper**, then push **Purge** to start a purge cycle.

c. Push **OK**.

**2.** Put the inlet tubing in the sample.

**3.** Push **Read**. The sipper unit completes a sipper cycle, then the turbidimeter measures the sample.

If the Purge Start setting is set to Auto, a purge cycle is completed after the measurement.

If the Purge Start setting is set to Manual, the Read button changes to a Purge button.

**4.** Push **Purge** to start a purge cycle if applicable.

**5.** Do steps **3** and **4** again until the sample measurements are complete.

**6.** Do a purge. Refer to step **1**.

## Prepare for storage

Before short-term storage (three hours or less), flush the flow cell assembly and tubing. Refer to [Flush the flow cell assembly and tubing](#) on page 11. Keep the flow cell assembly full of distilled or deionized water to prevent air locks and the collection of particles on the parts.

Before long-term storage, remove and clean the flow cell assembly. Do steps [1 to 4](#) of [Clean the flow cell assembly](#) on page 12. Keep the flow cell assembly full of distilled or deionized water to prevent air locks and the collection of particles on the parts.

## Maintenance

### Clean the instrument

Clean the exterior of the instrument with a moist cloth and a mild soap solution and then wipe the instrument dry as necessary.

### Clean the flow cell assembly

#### ▲ CAUTION

Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.

#### NOTICE

Do not touch or scratch the flow cell. Contamination or scratches on the glass can cause measurement errors.

Clean the flow cell assembly before initial use and as necessary to remove contamination.

#### Items supplied by the user:

Laboratory glass-cleaning detergent

Laboratory detergent

Distilled or deionized water

Soft, lint-free cloth

Silicone oil

Oiling cloth

1:1 Hydrochloric acid<sup>1</sup>

1. Remove and disassemble the flow cell assembly. Refer to [Figure 3](#).

2. Clean the glass parts as follows:

- a. Clean the internal and external surfaces of the glass parts with a laboratory glass-cleaning detergent.
- b. Fully rinse the glass parts many times with distilled or deionized water.  
*Note: Steam sterilize the tubing, flow cells, and caps in the flow cell assembly if necessary.*
- c. For low-range turbidity samples or dilution water measurements, clean the internal and external surfaces of the glass parts with 1:1 hydrochloric acid.


<sup>1</sup> Only necessary for low-range turbidity measurements

- Fully rinse the glass parts many times with dilution water (not distilled or deionized water).
3. Dry the exterior surfaces of the glass parts with a soft, lint-free cloth. Keep the interior surfaces wet.
  4. Clean the internal and external surfaces of the plastic parts and tubing with a laboratory detergent and warm water.
  5. Apply silicone oil to the exterior of the flow cell as follows:
    - a. Hold the flow cell by the top only to minimize dirt and fingerprints.
    - b. Clean the exterior surface of the flow cell with a soft, lint-free cloth to remove water spots and fingerprints.
    - c. Apply a small drop of silicone oil from the top to the bottom of the flow cell.
- Note: Use only the provided silicone oil. This silicone oil has the same refractive index as the flow cell glass and masks minor glass scratches.*
- d. Use the oiling cloth to apply the oil equally to the surface of the flow cell. Remove most of the oil. Make sure that the vial is almost dry.

*Note: Put the oiling cloth in a plastic storage bag to keep the cloth clean.*

6. Assemble and install the flow cell assembly. Do the steps in [Figure 3](#) in the opposite direction.

**Figure 3 Disassemble the flow cell assembly**

## Replace the O-rings

If the flow cell assembly has a leak, replace the O-ring in the top cap and bottom cap with the replacement O-rings supplied. Refer to [Figure 2](#) on page 7.

## Replace the sample tubing

At intervals, replace the sample tubing. Contamination (e.g., microbiological growths) is not easy to remove from the interior surface of the tubing.



After the tubing is replaced, flush the tubing. Refer to [Flush the flow cell assembly and tubing](#) on page 11.



## Replace the pump tubing

When damage or wear is seen on the pump tubing, replace the pump tubing as shown in the illustrated steps that follow.

**Item to collect:** Pump tubing (Lagoprene<sup>®</sup>) pre-assembled with peristaltic pump cover and connections

After the tubing is replaced, flush the tubing. Refer to [Flush the flow cell assembly and tubing](#) on page 11.

## Troubleshooting

Error	Description	Solution
Sipper module not connected. Please check cable.	There is no USB connection between the sipper and the turbidimeter.	Examine the USB cable. Make sure that the cable length is not more than 1 m (3.3 ft). Make sure that the USB cable is connected to the turbidimeter and the sipper.
Please check sipper and tube.	The pump tubing is not correctly installed.	Loosen the peristaltic pump cover, then install it again. Refer to the illustrated steps in <a href="#">Replace the pump tubing</a> on page 14. Operate the pump for 15 seconds. Push down the lever and make sure that the tube is correctly installed around the rollers.

## Replacement parts and accessories

### ⚠ WARNING

Personal injury hazard. Use of non-approved parts may cause personal injury, damage to the instrument or equipment malfunction. The replacement parts in this section are approved by the manufacturer.

**Note:** Product and Article numbers may vary for some selling regions. Contact the appropriate distributor or refer to the company website for contact information.

## Replacement parts

Description	Item no.
Flow cells (2x)	4709500
Top or bottom cap, flow cell assembly	4744800
Flow cell cover and frame	9650400
Oiling cloth	4707600
Pump tubing, Lagoprene, pre-assembled with peristaltic pump cover and connections	LZV877
Silicone oil	126936
SIP 10 sipper unit with tubing kit	LQV157.99.60002
Tygon tubing <sup>2</sup> , 1/8-in. ID, 1/4-in. OD, 1/16-in. thickness	11012
Tubing kit, includes Tygon tubing, 1/4-in. OD, 1.70 m (5.6 ft) and LZV877 <i>Note: Y fitting is not included.</i>	LZV940
USB cable, type AB, 1 m (3.3 ft)	LZQ104

## Accessories

Description	Item no.
Micro fiber cloth (for vial cleaning)	LZY945

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P.O. Box 389, Loveland, CO 80539-0389 U.S.A.
(800) 227-4224 (U.S.A. only)