ELECTROSURGICAL PENCIL EN

• Store in a Cool Dry Place • Device Rated for 2.5k Vpeak
WARNING: This medical device cannot be effectively deaned and/or sterilized by the user and therefore cannot be safely reused. It is intended for single use only. Any attempt by the user to clean and resterilize this device may result in bioincompatibility, infection, or other risks of device failure to the patient.

WARNING: Improper electrode installation may result in injury to the patient or operating room personnel by arcing at the electrode pencil connection.

WARNING: Keep active electrodes away from patient when not in use. Accessories should be placed in a clean, dry, nonconductive and highly visible area when not in use.

WARNING: This product is sterilized utilizing Ethylene Oxide. These products can expose you to chemicals including Ethylene Oxide (EtO), which is known to the state of California to cause cancer and/or birth defects or other reproductive harm. For more information go to **www.P65Warnings.ca.gov**.

Intended Use:

This device is intended to be used as a foot switch controlled monopolar accessory in conjunction with a compatible electrosurgical generator for tissue cutting and coagulation.

Indications: The distal end of the pencil has a 3/32” (2.36mm) collet that will accept a variety of electrosurgical electrodes with 3/32” (2.36mm) diameter shafts. The foot controlled pencil plug end is used with a foot control adapter, Bovie Part Number A1205A or equivalent. The foot control pencil must be used with a monopolar footswitch.

Contraindications:

Do not use in the presence of flammable materials. Fire could result.

This device should never be used when:

- There is visible damage to the pencil body or cable,
- It has not been inspected and tested prior to use,
- The device fails the test procedure, or
- The accessory electrode is not properly seated in the pencil electrode receptade.

Testing:

- Plug pencil connector into the appropriate adapter in the electrosurgical unit.
- Turn on ESU power and set power to zero in cut and coagulate functions.
- Hold pencil away from all persons and depress the CUT foot control button. Verify that the ESU signifies the activation of the CUT function.
- Repeat step 3 to test the COAG function.

Important:

- Pencil, cable, connector and electrode should be carefully inspected for defects prior to use.
- Do not allow cable to contact patient skin during activation.
- Do not permit cable to be in parallel and in close proximity to leads of other devices.
- Do not wrap cable around handles of metallic objects such as hemostats or clamp cable to drapes.

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