<div> <div><span>0001311 Rev. 1</span></div> <div><span></span></div> </div>	<div> <div>For product information, contact: <p>Yderligere produktoplysninger fås ved at kontakte: <p>Pour obtenir des informations sur le produit, veuillez contacter<span> </span>: <p>Für Produktinformation kontaktieren Sie: <p>Για πληροφορίες σχετικά με το προϊόν, επικοινωνήστε με: <p>Per informazioni sul prodotto, contattate: <p>製品情報、お問い合わせ先: <p>Para obtener información del producto, póngase en contacto con: <p><b><span><span></span></span> symmetry   olsen®</b> <p>3230 Commerce Center Place <p>Louisville, Kentucky, 40211 <p>(502) 813-5022 → (800) 297-6344 <p>See our products at <p>Se vores produkter på <p>Découvrez nos produits sur <p>Siehe unsere Produkte auf <p>Δείτε τα προϊόντα μας στην ιστοσελίδα <p>È possibile vedere i nostri prodotti all'indirizzo web <p>当社の製品は以下でご覧いただけます <p>Consulte nuestros productos en <p><b>http://www.symmetrysurgical.com</b></p></p></p></p></p></p></p></p></p></p></p></p></p></p></p></p></p></p></p></p></div> </div>
<div> <p><b>℞ ONLY</b></p></div>	
<div> <p><b>OLSEN MEDICAL® Single Use Bayonet Bipolar Irrigating Forceps</b> <p><b>OLSEN MEDICAL® Bipolær bajonetskylletang til engangsbrug</b> <p><b>Pinces bipolaires d'irrigation à baïonnette à usage unique d'OLSEN MEDICAL®</b> <p><b>OLSEN MEDICAL® bipolares Bajonett-Spülungsfassinstrument zur Einmalverwendung</b> <p><b>OLSEN MEDICAL® Μονής Χρήσης Αιχμηρή Διπολική Λαβίδα Καταιόνησις</b> <p><b>Pinze a baionetta bipolari ad irrigazione OLSEN MEDICAL® monouso</b> <p><b>OLSEN MEDICAL® 使い捨てバイオネット洗浄用バイポーラ鉗子</b> <p><b>Pinzas de irrigación bipolares de bayoneta de un solo uso de OLSEN MEDICAL®</b></p></p></p></p></p></p></p></p></div>	
<div> <p><b>ENGLISH</b></p></div>	
<div> <p><b>Description:</b></p> <p>The OLSEN MEDICAL® Single Use Bayonet Bipolar Irrigating Forceps has a nonpyrogenic fluid pathway and is packaged sterile, unless otherwise noted. It consists of an inlet pathway running into a transparent rigid PVC tube attached to the left hand tine of the bayonet forceps. This gives the physician a visual of the fluid pathway which facilitates the local fluid irrigation through the forceps. The device is for Single Use Only and is designed for use with bipolar electro<span>surgical</span> units.</p> <p><b>Intended Use (per EU definition) / Indications (per FDA definition):</b> <p>The OLSEN MEDICAL® Single Use Bayonet Bipolar Irrigating Forceps (Catalog No. 20-0xxx) is a single use product sold sterile and is intended for use in electro<span>surgery</span> for coagulation and irrigation of tissue. This device is intended for use with the OLSEN MEDICAL® Integrated Irrigation Tubing and Bipolar Cord Set or similar design of Bipolar Cord and Irrigation Tubing. The Single Use Bayonet Bipolar Irrigating Forceps are intended for use with a maximum voltage of 500 Volts.</p> <p><b>Handling:</b></p> <ol style="list-style-type: none"><li>Inspect package for damage before use. Discard if damaged or seal is broken.</li> <li>Inspect device after opening in a sterile field. Discard if damaged.</li> <li>Remove and discard the tip protector on all forceps before use.</li> <li>Device use and application to be determined by physician.</li></ol> <p><b>Storage:</b> <p>These devices should be stored in their original labeled boxes or in protected storage to prevent damage during handling. Take precautions to ensure that the sterile devices are stored in areas where it is completely dry, and protected from extremes of temperature and humidity.</p> <p><b>Sterility:</b> <p>Sterility is guaranteed up to the expiration date unless package is opened or seal is broken.</p> <p><b>Sterilization:</b> <p>Single Use Only (intended to be used on one patient during one procedure). Dispose after use. Do not resterilize. Reusing or reprocessing single use devices may cause damage to the device which in turn may cause unnecessary harm to the user and/or patient. Reuse or reprocessing is not recommended.</p> <p><b>Notices:</b></p> <ul style="list-style-type: none"><li>Electrosurgical equipment and accessories are intended to be used by health professional educated in their use.</li> <li>After Use or upon determining product can no longer be used due to wear or damage, safely discard in accordance with established procedures for biohazardous waste.</li> <li>Any serious incident that has occurred in relation to the device should be reported to Symmetry Surgical Inc., the governing body, and the Competent authority of the member state in which the incident occurred.</li> <li>Please refer to the manufacturer of the generator for warnings, precautions, contraindications, undesirable side-effect, measures to be taken, and limitations of use for the electro<span>surgical</span> system and accessories.</li> <li>For applicable IFU/additional information: Visit <a href="http://www.symmetrysurgical.com">www.symmetrysurgical.com</a>.</li></ul> <p><b>Caution:</b></p> <ul style="list-style-type: none"><li>Inspect every device before use. Damaged electro<span>surgical</span> devices may cause electrical burns to the patient or physician.</li> <li>Use a wet sponge to remove eschar build up from the forceps tips.</li> <li>Active electro<span>surgical</span> devices must be placed away from the patient when not in use. A holster is recommended.</li> <li>Consult RF generator Owner's Manual. Use only in the active RF generator output connections.</li> <li>Do not perform electro<span>surgery</span> in the presence of flammable anesthetics or other flammable gases, fluids or objects, or in the presence of oxidizing agents, as fire could result.</li> <li>Avoid unnecessary or prolonged activation. Activate only when touching or immediately next to the intended point of contact.</li> <li>Keep the voltage/power settings at the lowest possible level to achieve the desired effect.</li> <li>After use, safely discard in accordance with established procedures for biohazardous waste.</li></ul> <p><b>Warning:</b> <p>Improper use of any element in electro<span>surgical</span> systems may cause serious burn, shock, and fire or explosion hazard. The use of the device with safely functioning and compatible handles, cords, electro<span>surgical</span> generators and accessories is the user's responsibility. Read and understand all precautions and instructions before attempting to use any electro<span>surgical</span> devices. This product is for SINGLE USE ONLY; DO NOT RESTERILIZE. Use aseptic technique in all phases of handling. Olsen Medical Single Use devices have not been designed to undergo or withstand any form of alteration, such as disassembly, cleaning or resterilization, after a single patient use. Reuse can potentially compromise device performance and any usage beyond the design intent of this single use device may result in unpredictable loss of functionality.</p> <p><b>Warning:</b> <p>Avoid HF output mode/power settings where maximum output voltage may exceed the stated maximum voltage.</p> <p><b>Warning:</b> <p>Do not make any modification to these devices. Modification may cause burn or injury.</p> <p><b>Warning:</b> <p>Due to concerns about the carcinogenic and infectious potential of electro<span>surgical</span> byproducts (such as tissue smoke plume and aerosols which may include toxic gasses and vapors, live and dead cellular material, and viruses), protective eyewear, filtration masks, and effective smoke evacuation equipment should be used in both open and laparoscopic procedures. Contact Symmetry Surgical at <a href="mailto:customerservice@symmetrysurgical.com">customerservice@symmetrysurgical.com</a> for additional information or to inquire about our smoke evacuation solutions.</p></p></p></p></p></p></p></p></p></div>	
<div> <p><b><span><span></span></span> Warning:</b> This product is sterilized utilizing Ethylene Oxide. These products can expose you to chemicals including Ethylene Oxide (EtO), which is known to the state of California to cause cancer and/or birth defects or other reproductive harm. For more information go to <a href="http://www.P65Warnings.ca.gov">www.P65Warnings.ca.gov</a>.</p></div>	
<div> <p><b>Surgical Assembly Instructions:</b></p> <ol style="list-style-type: none"><li>Ensure that the Irrigation Tubing and Bipolar Cord set is prepared for connection to the instrument. If the cord set is included in a kit make sure that the cord is prepared prior to connecting the forceps.</li> <li>Remove the Single Use Bayonet Bipolar Irrigating Forceps from the peel pouch using aseptic technique.</li> <li>Connect the twin pin connector of the OLSEN MEDICAL® Single Use Bayonet Bipolar Irrigating Forceps into the female socket connector of the OLSEN MEDICAL® Integrated Irrigation Tubing and Bipolar Cord Set. Connect the set's male LUER-LOK® connector to the female irrigation port on the forceps.</li></ol></div>	
<div> <p><b>Warranty:</b></p> <p>All products are warranted to be free from defects in workmanship and materials when used for its intended purpose. Any product that proves defective in workmanship or materials will be repaired or replaced at our discretion free of charge. Repair, modification, or alteration to any product that is performed by anyone not authorized by Olsen Medical® will result in the immediate loss of warranty.</p> <p>LUER-LOK® is a registered trademark of Becton, Dickinson, and Company™. OLSEN® is a registered trademark of Olsen Medical.</p></div>	
<div> <p><b>DANSK - DANISH</b></p></div>	
<div> <p><b>Beskrivelse:</b></p></div>	

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