Sterile Distraction Screws

Standard Screws: DS-0012, DS-0014, DS-0016

Quick Start Screws: DS-0022, DS-0024, DS-0026

USER GUIDE

Company Information:

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Product Description:

The TeDan Surgical Innovations, Inc. (TSI) Distraction Screws are single use, STERILE devices available in three screw lengths (12mm, 14mm, and 16mm) and are supplied sterile, 1 screw per pouch and 5 pouches per box.

REF#'s:

	Type	Screw Length (mm)	Qty of Screws/Pouch	Qty of Pouches/Box
DS-0012	Standard	12	1	5
DS-0014	Standard	14	1	5
DS-0016	Standard	16	1	5
DS-0022	Quick Start	12	1	5
DS-0024	Quick Start	14	1	5
DS-0026	Quick Start	16	1	5

Surgical Setting:

Always ensure use of TSI instruments in a strictly sterile clinical environment.

Indication for Use:

Anterior cervical discectomy fusion surgery.

Intended Patients Groups:

Patients for which cervical spinal surgery using an anterior approach is required.

Intended User Groups:

Licensed practicing physicians who have been informed on the use of distraction screws used with TSI retractor systems. It is also expected that clinical assistants are knowledgeable in setting up the retractor and distractor system.

Contraindications:

No known contraindications.

Intended Purpose:

The TSI Distraction Screws are used in combination with distractor frames to distract two adjacent vertebral bodies. Products have been designed to be used with TSI Distractors and Screwdrivers.

Intended Clinical Benefits:

Provides a temporary anchoring mechanism for the distractors during cervical spine surgery.

Instructions for use:

1. Inspect package prior to use and do not use if the package is damaged or sterility has been compromised. If packaging has been damaged or unopened before use, DO NOT USE. Report any damaged packaging to TSI.
2. Select the desired length of distraction screw. The length should be chosen based on patient anatomy and surgeon experience.
3. Start with a pilot hole using the Phantom Twist Drill (CS-0112) or the Phantom Twist Drill with handle (CS-0109).
4. Drill through the Phantom Drill Guide (ie. CS-0110) or Phantom Universal Drill Guide (ie. CS-0119).
5. Using the Phantom Distraction Screw Driver (CS-0114) screw the distraction screw thru the drill guide until two fingers tighten or until screw is fully seated. Phantom Quick Start Screws are only compatible with TSI Phantom Springless Distraction Screw Drivers (CS-0114).
6. Repeat steps 1 thru 4 for the placement of the 2nd screw.
7. Remove the screwdriver and drill guide.
8. The above procedure ensures parallel placement of the distraction screws Do not bend the pins after insertion as this may cause failure of the screw!
9. Place appropriate distractor over the distraction screws and open as required.

ENGLISH

10. Remove distraction screws after use with the Phantom Distraction Screw Driver (CS-0114).
11. TSI distraction screws are single use only products. After the procedure, dispose of the distraction screw in an approved sharps container.
12. If needed, the screw or any fragment thereof can be located by means of an X-Ray.

Reporting Serious Incidents:

Any serious incidents that has occurred in relation to the device should be reported to TSI and the Competent Authority of the Member State in which the user and/or patient is established. Adverse events and incidents can be reported to TSI per contact information on the device packaging.

Warning/Precautions:

1. CAUTION: US Federal law restricts this device to sale by or on the order of a physician.
2. This product is supplied STERILE and Single Use Only. Do not reuse. Reuse may compromise the structural integrity of the device and/or lead to device failure. Reuse also presents biological hazards associated with disease transmission and immune/allergy issues, some of which could cause severe illness or be fatal.
3. Do not bend the pins after insertion as this may cause failure of the screw.
4. The decision to use Distraction Screws is based upon suitability of the patient's bone condition.
5. Phantom Quick Start Screws are only compatible with TSI Phantom Springless Distraction Screw Driver (CS-0114).
6. Do not use this device for any purpose other than the stated intended use.
7. Potential Adverse Events / Residual Risks: Infection / Adverse Reaction, Delay in Surgery, Injury.

Sterilization and Shelf Life:

TSI Distraction Screws are STERILE single use only and are sterilized by ethylene oxide. They may not be re-sterilized or re-used. Sterile packages have a five (5) year shelf life and should be stored at room temperature.

See www.tedansurgical.com/symbolsglossary for the Symbol Glossary.

Company Information:

For more information on this product or any products mentioned in this IFU please contact:


TeDan Surgical Innovations, Inc.
12615 West Airport Blvd, Suite 200
Sugar Land, TX 77478 USA

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TeDan Surgical Innovations B.V.