The GREENBERG™ Retractor and Handrest, A Universal System

Instruction Manual
Featuring the
“3-Level Staircase Concept”

Symmetry Surgical
3034 Owen Drive
Antioch, TN 37013 USA
☎ 1-800-251-3000
FAX: 1-615-964-5566
www.symmetrysurgical.com

EC REP

Symmetry Surgical GmbH
Maybachstraße 10
78532 Tuttlingen, Germany
☎ +49 7461 96490
FAX: +49 7461 77921

CE 2797

LCN 500-1500-85/M

© 2002–2020 Symmetry Surgical
Revised 08/2020

ENGLISH

IMPORTANT INFORMATION

Please Read Before Use

The GREENBERG Retractor and Handrest, A Universal System

Instruction Manual Featuring the “3-Level Staircase Concept”

Table of Contents

	Page
Introduction	3
Description, Indications, Contraindications	3
Warning, Cautions	3
Components of the GREENBERG Neurosurgical Retractor System	3
Suggested Assembly and Use Instructions	
1. Craniotomy Position: Patient Prone	4
2. Craniotomy Position: Patient Lateral (Park Bench)	6
3. Craniotomy Position: Patient Sitting	6
4. Laminectomy Position: Patient Prone	7
An Overview of the Staircase Concept	7
Care and Maintenance	8
References	8
Ordering Information	9
Service and Repair	9
Warranty	9

Introduction

by I.M. Greenberg, M.D.

Since its introduction in 1975, the GREENBERG Retractor and Handrest System has been most frequently used in conjunction with the Gardner Headrest as a support base. Many surgeons, however, prefer to use the Symmetry Surgical™ Skull Clamp because of its mechanical advantages.

The GREENBERG System is a set of instruments which surrounds the operative site. First are a pair of 6 inch (15.2 cm) long, steel primary bars, 0.5 inch (1.3 cm) in diameter, with clamps for attachment to the headrest. Longer secondary bars, 12 inches (30.5 cm), are then clamped to either the primary bars or to each other. FLEXBAR™ Retractor Arms are attached to these secondary bars and serve to support and position various instruments, including handrests, retractor blades, suctioning devices, dissectors, scissors, or drills. Most recently the increasingly sophisticated manner in which this system is used has led to the practice of conceptualizing three progressively-higher levels of instruments, based upon function.

Level I (Figure 1)

This is the level at which the system is attached to either the headrest or the railing of the operating room table and is the physical foundation of the system. The function at this level is to create a stable perimeter about the operative field for mechanical support of normally hand-held instruments. The surgical tasks at this level consist primarily of retraction, dissection and suction.

Level II (Figure 2)

This level has one function—handrests for surgeon support.

Level III (Figure 3)

Structurally, Level III consists of a network of flexible arms on which other instruments will be supported above the handrests and kept in or near the operative site during surgery. The functions at this level are dictated by the specific instruments chosen. These instruments are referred to as “instruments of action.”

In creating new instruments, I have kept in mind the need to perform a surgical procedure with greater ease and control, while using a lesser amount of physical effort, and being more independent from operating room personnel. These instruments and instructions for their use are presented in this manual. Proper utilization of these instruments will permit both physical and conceptual order: physically, by compartmentalizing instrument structure into levels; and conceptually, by assigning specific tasks to those levels.

Description

The GREENBERG Retractor and Handrest System provides a means for self-retaining retraction, support of the surgeon's hands, and support of instrumentation during neurosurgery. It must be used with either a Symmetry Surgical Skull Clamp, catalog no. 19-1042 or 19-1043*, or a Gardner Skull Clamp, catalog no. 19-1020, for cranial procedures.

Indications

The GREENBERG Retractor and Handrest System is indicated for self-retaining retraction, support of the surgeon's hands, and support of surgical instruments.

Contraindications

The GREENBERG Retractor and Handrest System is not designed, sold, or intended for use except as indicated.

WARNING

Failure to read and follow instructions in this manual may lead to device malfunction, which could result in patient injury or death.

CAUTIONS

RESET THE DRAWBAR TO ZERO AFTER EACH USE of the FLEXBAR Arms. Do not store or sterilize FLEXBAR Arms while under tension. See the “Care and Maintenance” section of this instruction manual for preventive maintenance.

Forcing FLEXBAR Arms to move against preset tension can cause the cable to wear and possibly break, and the ball joints to become scored causing FLEXBAR Arms to drift. See the “Care and Maintenance” section of this instruction manual for preventive maintenance.

*Symmetry Surgical no longer sells Symmetry Surgical Skull Clamps.

The GREENBERG Neurosurgical Retractor System

GREENBERG Primary Bar, Figure 4, attaches to skull clamp.
Catalog no. 50-1507

The GREENBERG Primary Bar is the base of the GREENBERG Retractor System. The primary bar is attached to the skull clamp by turning the vise (A) counterclockwise until the opening (B) is sufficient to surround the base or arms of the skull clamp.

GREENBERG Secondary Bar, Figure 5, attaches to primary bar.
Catalog no. 50-1508

The GREENBERG Secondary Bar is used to construct the framework on which retraction, handrests, and instruments of action are attached. The secondary bars attach to the primary bars and to each other by means of a vise grip (A).

GREENBERG Floating Secondary Bar, Figure 6, attaches to primary or secondary bar, or short coupler.
Catalog no. 50-1516

The GREENBERG Floating Secondary Bar is used in the same way as the GREENBERG Secondary Bar; however, the floating secondary bar incorporates a universal joint (A) enabling conical adjustment of the angle of the bar. The vise clamp of the floating secondary bar (B) mounts onto the primary or secondary bar, or short coupler. It can also be mounted onto a flat metal edge of up to 0.50 in. (1.3 cm) thickness.

Turning the universal joint adjustment lever (C) counterclockwise permits adjustment of the rod. To stabilize the rod in the desired position, turn the joint lever clockwise until tightened. The slot in the collar of the universal joint permits 90° positioning of the floating secondary bar while maintaining stability. The joint lever is usually positioned upward for the convenience of the surgeon.

See assembly instructions for using the GREENBERG Floating Secondary Bar in conjunction with the Symmetry Surgical Skull Clamp.

NOTE: Floating secondary bars are not supplied with the GREENBERG Universal Retractor and Handrest Kit, and must be ordered separately; see “Ordering Information” section of this instruction manual.

GREENBERG Long Retractor Arm, 9 inches (22.9 cm) long, attaches to primary or secondary bar.
Catalog no. 50-1509

GREENBERG Short Retractor Arm, 6 inches (15.2 cm) long, attaches to primary or secondary bar.
Catalog no. 50-1511

The GREENBERG Retractor Arm, Figure 7, is a FLEXBAR Arm which can be locked in various positions. A mini-vise (A) at one end of the retractor arm holds all sizes of GREENBERG retractor blades. A clamp located at the other end attaches the retractor arm to either a primary or secondary bar. A tension-setting knob positions the retractor arm. To position the retractor arm, turn the tension-setting knob counterclockwise until the retractor arm becomes limp. Position the retractor arm as desired, then turn the tension-setting knob clockwise to tighten. For minor adjustments to the retractor arm, loosen the tension-setting knob slightly. Reposition the retractor arm, then tighten the knob. **NOTE:** When retractor blades are set in the field and subtle positioning is needed, no change in tension or readjustment of the knob is necessary.

CAUTION: Forcing the FLEXBAR Arm to move against preset tension can cause the cable to fray and eventually break. Also, the ball joints can become scored causing the FLEXBAR Arm to drift. See “Care and Maintenance” section of this instruction manual for preventive maintenance.

GREENBERG Microinstrument Holder, Figure 8, attaches to secondary or floating secondary bar.
Catalog no. 50-1529

The GREENBERG Microinstrument Holder is a FLEXBAR Arm used to support continuous suction and those microinstruments used repeatedly within the retracted field. Secure the microinstrument at a point along its shaft in the maxi-vise (A) or, preferably, secure the microinstrument at its handle so that the microinstrument holder is out of the field, preventing obstruction of the surgical field.

Refer to “GREENBERG Long/Short Retractor Arm” for positioning instructions.

CAUTION: Forcing the FLEXBAR Arm to move against preset tension can cause the cable to wear and eventually break. Also, the ball joints can become scored causing the FLEXBAR Arm to drift. See “Care and Maintenance” section of this instruction manual for preventive maintenance.

**GREENBERG Small Instrument Holder, Figure 9, attaches to secondary or floating secondary bar.
Catalog no. 50-1523**

The GREENBERG Small Instrument Holder is a FLEXBAR Arm used when greater holding force than that offered by the microinstrument holder is needed. The small instrument holder uses a maxi-vise (A) to support instruments with shafts of up to 0.3 inch (7.9 mm) in diameter.

Refer to “GREENBERG Long/Short Retractor Arm” for positioning instructions.

CAUTION: Forcing the FLEXBAR Arm to move against preset tension can cause the cable to wear and eventually break. Also, the ball joints can become scored causing the FLEXBAR Arm to drift. See “Care and Maintenance” section of this instruction manual for preventive maintenance.

**GREENBERG Extra Large and Large Instrument Holders, Figure 10, attach to secondary or floating secondary bar.
Catalog no. 50-1524, Large, 50-1525; Extra Large**

The GREENBERG Extra Large and Large Instrument Holders are FLEXBAR Arms used to support heavy instruments, such as a small pneumatic drill, ultrasonic aspirator, ultrasonic scanning probe, or hand-held laser unit. The extra large and large instrument holders use a cantilever clamp (A) to support tools with shafts of up to 1 inch (2.5 cm) in diameter without drifting.

Refer to “GREENBERG Long/Short Retractor Arm” for positioning instructions.

CAUTION: Forcing the FLEXBAR Arm to move against preset tension can cause the cable to wear and possibly break. Also, the ball joints can become scored causing the FLEXBAR Arm to drift. See “Care and Maintenance” section of this instruction manual for preventive maintenance.

**GREENBERG Short Coupler, Figure 11, attaches to the secondary bar.
Catalog no. 50-1522**

The GREENBERG Short Coupler attaches to the GREENBERG Secondary Bar by means of a locking screw (A). The GREENBERG Handrest, Large, or Small Instrument Holders can be mounted on the GREENBERG Short Coupler post (B).

**GREENBERG Handrest, Figure 12, attaches to short coupler or secondary bar.
Catalog no. 50-1512**

The GREENBERG Handrest provides support needed to help minimize tremor. The head (A) can rotate almost 360°. The handrest mounts on the post of the short coupler or the secondary bar by means of a vise grip (B).

**GREENBERG Pattie Tray, Figure 13
Catalog no. 50-1513**

The GREENBERG Pattie Tray holds Symmetry Surgical Surgical Patties during the surgical procedure. The pattie tray attaches to the mini-vise of the GREENBERG Retractor Arm (A).

**GREENBERG Tapered Retractor Blades, Figure 14
Catalog nos. 50-1530, 50-1531, 50-1532, 50-1533, 50-1534**

The tapered retractor blades are used for more discrete retraction in the distal part of the field, while providing protection to more proximal areas of the brain along the length of the blade.

Craniotomy Position: Patient Prone
Suggested assembly for:

- Level I Retraction, Dissection, Suction**
- Level II Handresting (one handrest for each hand)**
- Level III Instruments of Action**

using the Gardner Skull Clamp (catalog no. 19-1020) or the Symmetry Surgical Skull Clamp, catalog no. 19-1042 or 19-1043* for craniotomy with the patient prone.

***Symmetry Surgical no longer sells Symmetry Surgical Skull Clamps.**

Level I—360° Retraction

Refer to Figures 15-1 through 15-8 when using the Gardner Skull Clamp. Refer to Figures 16-1 through 16-7 when using the Symmetry Surgical Skull Clamp.

NOTE: Position the mounting screws of the skull clamp towards the patient's shoulders, away from the area in which the surgeon is standing or sitting. The base of the skull clamp should be parallel to the floor (see Figure 15-1 or Figure 16-1). FROM THIS POINT, THE PATIENT IS DRAPED AND THE FIELD IS STERILE. The entire assembly process to be shown is carried out by sterile technique, with the GREENBERG Primary Bar clamps being applied on top of the sterile drapes, covering the headrest.

Gardner

Step 1—Attaching Primary Bars

The GREENBERG Primary Bars are mounted perpendicular to the base of the Gardner Skull Clamp and point toward the surgeon (see Figure 15-2).

Step 2—Mounting Secondary Bars

One GREENBERG Secondary Bar is placed on each GREENBERG Primary Bar (see Figure 15-3). These secondary bars are positioned so that they rest on top of the opposite primary bar.

A secondary bar is then mounted onto each of the secondary bars previously mounted (see Figure 15-4).

To complete a frame around the surgical field, a final secondary bar is attached (see Figure 15-5).

Step 3—Attaching Pattie Tray

Attach pattie tray to the secondary bar using a GREENBERG Retractor Arm. The tray is held by the mini-vise.

NOTE: Scalp retraction can be achieved by attaching small towel clips to the flap edge. Rubber bands can be attached to the clips and mounted on the frame. See Figure 15-6.

Step 4—Attaching Retractor Arms (Figure 15-7)

Retractor arms (four are shown) are attached to the frame to enable 360° retraction and dissection of the brain. (Scalp retractors omitted from illustration for simplicity.)

Step 5—Attaching Microinstrument Holders (Figure 15-8)

Attach microinstrument holders to the primary and/or secondary bars. Microsuction and microdissection instruments can now be placed in the field.

Symmetry Surgical

Step 1—Attaching Primary Bars

The GREENBERG Primary Bars are mounted on the vertical arms of the Symmetry Surgical Skull Clamp so that they are parallel with the base of the skull clamp and perpendicular to the patient (see Figure 16-2).

NOTE: The primary bar applied to the ratchet arm must be applied low and must completely encircle that arm. When the primary bar clamp is fully tightened, it must make contact with the skull clamp arm at three points.

Step 2—Mounting Secondary or Floating Secondary Bars

One secondary bar is then mounted onto each primary bar (see Figure 16-3).

To complete a frame around the surgical field, two secondary bars are attached (see Figure 16-4). **NOTE:** Attaching floating secondary bars will allow symmetry. Floating secondary bars are available separately. Using secondary bars will result in an asymmetrical box.

Step 3—Attaching Pattie Tray

Attach pattie tray to the secondary or floating secondary bar using a GREENBERG Retractor Arm. The tray is held by the mini-vise.

NOTE: Scalp retraction can be achieved by attaching small towel clips to the flap edge. Rubber bands can be attached to the clips and mounted on the frame. See Figure 16-5.

Step 4—Attaching Retractor Arms (Figure 16-6)

Retractor arms (four are shown) are attached to the frame to enable 360° retraction and dissection of the brain. (Scalp retractors omitted from illustration for simplicity.)

Step 5—Attaching Microinstrument Holders (Figure 16-7)

Attach microinstrument holders to the primary, secondary, and/or floating secondary bars. Microsuction and microdissection instruments can now be placed in the field.

Level II—Dual Handrests

Step 1—Attaching Short Couplers

(Refer to Figure 17-1 for use with the Gardner Skull Clamp or to Figure 18-1 for use with the Symmetry Surgical Skull Clamp.)

Attach two short couplers (A and B) to the frame to support handrests. Short couplers can also be used to support extra large, large or small instrument holders. (Retraction/Dissection omitted from illustrations for simplicity.)

NOTE: Instead of the short couplers, two secondary bars can be mounted on the frame to support handrests; however, this assembly is less desirable due to limitations in positioning the handrests.

Step 2—Attaching Handrests

(Refer to Figure 17-2 for use with the Gardner Skull Clamp or to Figure 18-2 for use with the Symmetry Surgical Skull Clamp.)

Mount a handrest onto each of the short couplers. It is suggested that the surgeon place the handrests directly on top of the surgical field.

Coarse adjustments are made by moving the handrest back and forth in the operative field. Fine adjustments of the handheld instrument are made by the supported hand holding the instrument between thumb, index, and third finger, utilizing extension, flexion, and rotation (see Figure 19).

To provide physical support and complete independence of movement for each hand at Level II, the use of two handrests is recommended for most procedures; with surgical activity taking place in the trough formed between the handrests (see Figure 20).

Level III—Instruments of Action

Level III mechanically supports and stores additional instruments in or near the operative site during surgery. For example, secondary and floating secondary bars attach to the structure at Level I and are positioned above Level II. Instrument holders of various sizes attach to the secondary and/or floating secondary bars, with instruments such as pneumatic drills, ultrasonic aspirators, ultrasound probes, handheld lasers, and scissors attached to the instrument holders. These instruments are referred to as “instruments of action.” The surgical functions at Level III are determined by the specific instruments of action chosen.

Because the instruments of action are supported independently in or near the surgical field, the surgeon can return to these instruments and use them as needed. This lessens the surgeon’s reliance on other surgical team members. By adjusting the tension of the instrument holders and secondary or floating secondary bars, the surgeon is able to work within approximately a two-inch cube (see Figure 21). When planned and executed correctly, this cube encompasses the operative site so that many activities in the procedure are accomplished without readjustments.

Step 1—Attaching Floating Secondary Bar to Frame

(Refer to Figure 22-1 for use with the Gardner Skull Clamp or to Figure 23-1 for use with the Symmetry Surgical Skull Clamp.)

Attach a floating secondary bar to the frame to hold instruments of action.

Step 2—Attaching Extra Large, Large and Small Instrument Holders

(Gardner Skull Clamp: refer to Figures 22-2 through 22-6. Symmetry Surgical Skull Clamp: refer to Figures 23-2 through 23-6.)

Extra large, large and small instrument holders, mounted on the floating secondary bar, support pneumatic drills, ultrasonic aspirators, ultrasound probes, and handheld lasers. Additional floating secondary bars and instrument holders can be added as needed to support instruments of action at Level III.

NOTE: To help support the weight of tubing, use a sterile rubber band as shown in Detail A.

Craniotomy Position: Patient Lateral (Park Bench)

Suggested assembly for:

- Level I Retraction, Dissection, Suction
- Level II Handresting (one handrest for each hand)
- Level III Instruments of Action

using the Symmetry Surgical Skull Clamp, catalog no. 19-1042 or 19-1043* for craniotomy with patient in lateral, or “Park Bench” position.

*Symmetry Surgical no longer sells Symmetry Surgical Skull Clamps.

Level I—Retraction

THE PATIENT IS DRAPED AND THE FIELD IS STERILE. The entire assembly process shown is carried out by sterile technique, with the GREENBERG Primary Bar clamp being applied on top of the sterile drapes covering the headrest.

Step 1—Attaching Primary and Floating Secondary Bars

The primary bar is mounted to the Symmetry Surgical Skull Clamp on top of the sterile drapes. Attach a floating secondary bar to the primary bar as shown in Figure 24. **NOTE:** The floating secondary bar is not included in the GREENBERG Universal Retractor and Handrest Kit and must be purchased separately. See “Ordering Information” section of this instruction manual.

Step 2—Attaching Retractor Arms, Microinstrument Holders (Figure 25)

Attach retractor arms and microinstrument holders to the floating secondary bar for retraction/dissection.

Step 3—Attaching Additional Retractor Arms and Microinstrument Holders (Figure 26)

If additional supported instrumentation from the opposite side of the craniotomy opening is required, an additional primary and floating secondary bar can be added as shown. Continuous microsuction is frequently mounted at this level with the advantage of being out of the traffic pattern of handheld instruments entering and leaving the field. Additional instruments for retraction/dissection can be mounted if needed.

Level II—Handrests

Step 1—Attaching Primary Bar (Figure 27)

Attach a primary bar to the Symmetry Surgical Skull Clamp, and one to the operating room table.

Step 2—Attaching Floating Secondary Bar (Figure 28)

Attach a floating secondary bar to the primary bar that was mounted onto the operating room table.

Step 3—Attaching Short Couplers (Figure 29)

Attach one short coupler to the floating secondary, and one to the primary bar.

Step 4—Attaching Handrests

One handrest is mounted on each side of the operative field, as shown in Figure 30.

Step 5—Counterretraction (Figure 31)

At this point, counterretraction and suction can be added.

Level III—Instruments of Action

Step 1—Attaching Primary Bar (Figure 32)

Mount a primary bar to the Symmetry Surgical Skull Clamp as shown.

Step 2—Attaching Floating Secondary Bar (Figure 33)

Attach a floating secondary bar to the primary bar.

Step 3—Attaching Instruments of Action

Secure desired instruments of action to floating secondary bar using extra large or large and small instrument holders. A pneumatic drill is shown in Figure 34; however, ultrasonic aspirators, ultrasound probes, and handheld lasers can all be supported at Level III.

Craniotomy Position: Patient Sitting

Suggested assembly for:

- Level I Retraction, Dissection, Suction
- Level II Handresting (one handrest for each hand)
- Level III Instruments of Action

using either the Gardner Skull Clamp or the Symmetry Surgical Skull Clamps, when patient is sitting:

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Level I—Retraction

Bring patient to a sitting position and apply Gardner or Symmetry Surgical Skull Clamp. FROM THIS POINT, THE PATIENT IS DRAPED AND THE FIELD IS STERILE. The entire process shown is carried out using sterile

technique with the GREENBERG Primary Bar clamp being applied on top of the sterile drapes covering the skull clamp.

Step 1—Attaching Primary Bar (Figure 35)

When using the Symmetry Surgical Skull Clamp, the primary bar must be positioned on the thickest vertical support, as shown. When using the Gardner Skull Clamp, the primary bar can be applied to either corner.

Step 2—Attaching Secondary or Floating Secondary Bars

Assembly continues with two secondary or floating secondary bars as shown in Figure 36. From this point on, the Gardner and Symmetry Surgical Skull Clamp applications are the same.

NOTE: The floating secondary bar is not included in the GREENBERG Universal Retractor and Handrest Kit and must be purchased separately. See “Ordering Information” section of this instruction manual.

Step 3—Attaching Retractor Arms (Figure 37)

Apply retractor arms relative to the direction needed for retraction/dissection.

Step 4—Attaching Pattie Tray and Additional Retractors

If bilateral retraction is necessary, requiring 360° of opportunity, additional secondary bars are applied as shown in Figure 38. The pattie tray is mounted in a convenient place.

Step 5—Attaching Microinstrument Holders and Microinstruments (Figure 39)

Microsuction and microdissector device(s) are mounted using the microinstrument holders. (Retraction has been omitted from the illustration for simplicity.)

Level II—Handrests

Possible variations for positioning the handrests follow.

Variation 1 (Figure 40)

The primary bar (I) is mounted onto the railing of the operating room table, on top of the sterile draping. The handrest (J) is attached to the primary bar. Coarse adjustments are made by moving the handrest back and forth from the operative field. Fine adjustments of the handheld instrument are made by the supported hand holding the instrument between thumb, index, and third finger, utilizing extension, flexion, and rotation.

Variation 2

Attach another primary bar to the railing of the operating room table, on top of the sterile draping. Mount a second handrest next to the first. The handrests may be positioned perpendicular to, and slightly oblique from, the operative field (see Figure 41), or in a straight line approximately parallel to the floor (see Figure 42). Coarse adjustment of the position of the hand-held instrument is accomplished by moving the hand across the horizontal support of the handrest. For dealing with superficial structures, the hand should slide laterally, while for deep structures, it slides medially. Fine adjustment remains the same.

Level III—Instruments of Action

Step 1—Attaching Secondary Bars

Mount two additional secondary bars as shown in Figure 43.

Step 2—Attaching Single Instrument of Action

A retractor arm with instrument of action (microscissors) is attached. For steadiness during use, instruments of action should pass over the top of the handrest so that the fingers are within reach of the instruments' handles while the hand is supported by the handrest.

Step 3—Attaching Additional Instruments of Action

A second instrument of action (pneumatic drill) can be mounted using a large instrument holder (see Figure 44). Steadiness is achieved by the hand resting on the handrest and the drill supported by the instrument holder. Thus, two advantages are presented: steadiness and accuracy similar to a drill press, and the maneuverability of a handheld instrument.

Another microscissors may be mounted on the left side and can function alone (see Figure 45) or in conjunction with a microscissors mounted at Level III on the right side. A microscissors would then be available from either side of the operative field.

For some procedures, it may be more comfortable for the surgeon to mount the pneumatic drill on the left side (see Figure 46).

Laminectomy Position: Patient Prone

Suggested assembly for:

Level I Retraction, Dissection, Suction

Level II Handresting (one handrest for each hand)

Level III Instruments of Action

when patient is in prone position for laminectomy surgery.

Level I—Retraction

THE PATIENT IS DRAPED AND THE FIELD IS STERILE. The primary bar and all other attachments are applied using sterile technique.

Step 1—Attaching Primary Bar (Figure 47)

The primary bar is mounted to the operating room table on top of the sterile drapes.

Step 2—Attaching Floating Secondary and Secondary Bars

Attach a floating secondary to the primary bar, positioning it approximately perpendicular to the floor. Attach a secondary bar to the floating secondary bar. Position it perpendicular to the floating secondary bar, across the patient, as shown in Figure 48. **NOTE:** The floating secondary bar is not included in the GREENBERG Universal Retractor and Handrest Kit and must be purchased separately. See “Ordering Information” section of this instruction manual.

Step 3—Attaching Retractor Arms, Microinstrument Holders, and Microinstruments (Figure 49)

Attach retractor arms. Mount microsuction and microdissector devices using the microinstrument holders.

Level II—Handrests

Step 1—Attaching Primary and Floating Secondary Bars (Figure 50)

Mount a second primary bar to the operating room table on top of the sterile drapes. Attach a floating secondary bar to the primary bar.

Step 2—Attaching Handrest

Mount a handrest to the floating secondary as shown in Figure 51. **NOTE:** To increase the options for positioning the handrest, attach a short coupler to the floating secondary; attach the handrest to the short coupler.

Coarse adjustments are made by moving the handrest back and forth from the operative field. Fine adjustments of the handheld instrument are made by the rested hand holding the instrument between thumb, index, and third finger, utilizing extension, flexion, and rotation.

Level III—Instruments of Action

Step 1—Attaching Floating Secondary Bar

Attach a floating secondary bar to the secondary bar that was added during Level I, see Figure 52. (Retractor arms and microinstruments have been omitted from illustrations for simplicity.)

Step 2—Attaching Extra Large, Large and Small Instrument Holders

Mounted on the floating secondary bar, extra large, large and small instrument holders support instruments of action. In addition to the pneumatic drill shown in Figure 53, Level III can support ultrasonic aspirators, ultrasound probes, and handheld lasers.

Figure 54 shows the system completely assembled for a complex laminectomy.

An Overview of the Staircase Concept (Figure 55)

The goal of the operating surgeon should be to compartmentalize work performance into three levels.

Level I includes retraction and dissection and has a low profile directly next to the operative field.

Level II is above Level I, separate and distinct from Level I and consists of handresting, with one handrest for each hand.

Level III is above Level II and consists of instruments mounted above the handrests. These instruments are within easy reach so the surgeon can move from Level II to Level III with minimal hand travel distance.

By the same token, the surgeon, while resting the hand at Level II, will be able to move down and make subtle adjustments during the retraction/dissection process at Level I.

The key is the concept of the surgeon's hands supported at Level II with movements from the handrest down to Level I, up to Level III.

The stability for most of the surgeon's activities will come from the handheld instruments at Level II resting on the handrest during surgery. *This platform is the hub, or focal point, of all the surgeon's work.*

With the increased use of microsurgical and advanced surgical technique, there has developed the necessity of conceptual models of surgery. The idea of separate levels of instrument structure and surgical performance is more than a way to describe the assembly of instruments.^(1,2,3,6) It compartmentalizes the operating surgeon's performance and brings order and neatness into the surgical process. By adhering to the Staircase Concept,⁽³⁾ tools seem to fit more easily where they are being used, performance becomes more systematic and there is minimal interference from "things getting in the way."^(6,8,9)

The increasing use of Level III for the mechanical support of instruments of action greatly increases the surgeon's efficiency. The key to Level III, with its many options, is to use only the part or parts of the system which are helpful for a given procedure.^(6,8,9) Instrument positions should be planned in advance so as to maximize efficiency in the use of space.^(5,7)

These instruments, their suggested uses, and the Staircase Concept of structure and function have developed to concepts of surgery, here condensed into seven ideas:

1. Work analysis of each step of a given procedure.
2. Decreasing repetitive movements (particularly instrument exchanges) during surgery.^(8,9)
3. Use of mechanical support for each hand to reduce tremor and fatigue while increasing precision.^(1,2,3)
4. Reliance upon mechanical support for instruments in the field.^(4,5,6)
5. Decreasing reliance upon other operating room personnel so as to reduce miscommunication and enable the surgeon to work at his own tempo.^(4,5,6,7,8,9)
6. Maintain versatility of the instruments so that new instruments can be added easily.
7. Reduce the number and required training of assisting personnel and thus reduce costs.^(7,8,9)

These goals are attainable. Surgeons can be more ergonomically efficient, and the instruments described in the preceding chapters are versatile, simple, and reliable. The techniques, of course, take time to learn, but as experience is gained, the surgeon will find that the time required to set up the system will diminish and the use of mechanically supported instruments will increase.

Care and Maintenance

For purposes of care and maintenance, it may be helpful to separate the instruments into bars (primary, secondary, and floating secondary bars), and FLEXBAR Arms (long and short retractor arms, handrests, microinstrument holders, extra large, large and small instrument holders).

Bars

Cleaning

After each surgery, use a soft brush and mild detergent to remove all residual debris from each instrument. It is recommended that the instruments and parts be ultrasonically cleaned.

Lubrication

It is extremely important that moveable parts be properly lubricated to keep them functional. It is recommended that all components be immersed in a water-soluble lubricant. PRESERVE[®] Concentrate (catalog no. 43-1033) is recommended for instruments to be steam sterilized.

FLEXBAR Arms

Cleaning

After each surgery, use a soft brush and mild detergent to remove all residual debris from each instrument. It is recommended that the instruments and parts be ultrasonically cleaned.

FLEXBAR Arms consist of a series of hollow cylinders through which a multi-strand steel cable runs. Attached to the cable is a swaged screw which is threaded into a drawbar. See Figures 56 and 57.

Setting the Drawbar to Zero

WARNING: RESET THE DRAWBAR TO ZERO AFTER EACH USE to help prevent "drifting" of the FLEXBAR Arms and damage to the drawbar. Do not store or sterilize FLEXBAR Arms while the cable is under tension.

NOTE: The FLEXBAR Arm is received by the customer with the drawbar at zero setting; the jointed arm has little or no slack and the drawbar is totally recessed in the well of the tension-setting knob. The FLEXBAR arm is stiffened during surgery by rotating the tension setting knob clockwise. The hub of the drawbar will protrude slightly as the cable is pulled. After much use, the drawbar may begin to protrude excessively (see Figures 58–59).

To relieve this condition, it is necessary to reset the drawbar to zero as follows:

1. Hold the FLEXBAR Arm of the instrument in the left hand and, with the forefinger of the right hand, spin the tension-setting knob counterclockwise until it meets resistance. See Figure 60.
2. Hold the FLEXBAR Arm in the left hand and, with the right hand, spin the entire clamping assembly clockwise. This will advance the threaded cable into the drawbar. See Figure 61–62.

Inspecting the Cable and Drawbar

1. With the left hand, hold the FLEXBAR Arm at the most proximal ball joint, and with the right hand, spin the entire clamping assembly counterclockwise for eight complete turns to loosen the ball joints. Inspect the cable for wear. If frayed (Figure 63–64), return for repair or replacement. **NOTE:** Fraying is most likely to occur where the cable meets the swaged screw.

2. At this point, the drawbar is at zero setting (fully recessed in the tension knob). Pull the tension-setting knob back to expose the threaded section of the drawbar and inspect the drawbar for wear. If threads are worn and/or metal filings can be seen (Figures 65–66), return the instrument for repair or replacement.

3. Push the drawbar into the slot. Hold the FLEXBAR Arm of the instrument in the left hand and, with the forefinger of the right hand, spin the entire clamping assembly clockwise. This will advance the threaded cable into the drawbar. When resistance is met, turn back one-half turn. See Figures 67–68.

Lubrication

It is extremely important that moveable parts be properly lubricated to keep these parts functional. It is recommended that all components be immersed in a water-soluble lubricant. PRESERVE Concentrate (catalog no. 43-1033) is recommended for instruments to be steam sterilized.

Sterilization

All components of the GREENBERG Retractor System require cleaning and sterilization before each use, including the first use.

The instruments can be steam or ethylene oxide sterilized. Refer to the sterilizer manufacturer's instructions for correct time, temperature, and pressure settings.

References

1. Greenberg, I.M.: "Self-Retaining Retractor and Handrest System for Neurosurgery," NEUROSURGERY, 8:205–208, 1981
2. Greenberg, I.M.: "New Options for the Neurosurgeon, Multiple Instrumentation in a Single Microsurgical Field," NEUROSURGERY, 9:566–583, 1981
3. Greenberg, I.M.: "Staircase Concept of Instrument Placement in Microsurgery," NEUROSURGERY, 9:696–702, 1981
4. Greenberg, I.M.: "Cerebral Aneurysm Rupture During Neurosurgery," NEUROSURGERY, 15:243–245, 1984
5. Greenberg, I.M.: "Rotation of Scissors While Fixed in the Microsurgical Field: A Technique for Increasing the Surgeon's Efficiency," MICROSURGERY, 5:218–227, 1984
6. Greenberg, I.M.: "Self-Retaining Retractors," in Wilkins, R.H., Rengachary, S.S., (eds): NEUROSURGERY, New York, McGraw-Hill Book Co., 1985

7. Greenberg, I.M.: "Intraoperative Ultrasonography with a Cystoscope for the Biopsy of a Deep-Seated Brain Lesion: Case Study," *NEUROSURGERY*, 19:49–58, 1986

8. Greenberg, I.M.: "Proximity Storage of Surgical Instruments in the Operative Field," *NEUROSURGERY*, 20:950–953, 1987

9. Greenberg, I.M.: "Neurosurgical Instrumentation," *Video Journal of Neurosurgery*, Vol. 2, Issue 3.

Ordering Information For Levels I And II

Catalog Number	Description
50-1500	GREENBERG Universal Retractor and Handrest Kit, which includes:
50-1507	GREENBERG Primary Bar, 2 each
50-1508	GREENBERG Secondary Bar, 4 each
50-1509	GREENBERG Long Retractor Arm, 4 each
50-1512	GREENBERG Handrest, 1 each
50-1513	GREENBERG Pattie Tray, 1 each
50-1514	GREENBERG Sterilizing Case, 1 each
50-1530	GREENBERG Tapered Retractor Blade, malleable, 0.75 in. (1.9 cm), 1 each
50-1531	GREENBERG Tapered Retractor Blade, malleable, 0.63 in. (1.6 cm), 1 each
50-1532	GREENBERG Tapered Retractor Blade, malleable, 0.50 in. (1.3 cm), 1 each
50-1533	GREENBERG Tapered Retractor Blade, malleable, 0.38 in. (9.5 mm), 1 each
50-1534	GREENBERG Tapered Retractor Blade, malleable, 0.25 in. (6.4 mm), 1 each
50-1507	GREENBERG Primary Bar, 1 each
50-1511	GREENBERG Short Retractor Arm, 2 each
50-1512	GREENBERG Handrest, 1 each
50-1516	GREENBERG Floating Secondary Bar, 3 each
50-1522	GREENBERG Short Coupler, 2 each
50-1529	GREENBERG Microinstrument Holder, 2 each
50-1530	GREENBERG Tapered Retractor Blade, malleable, 0.75 in. (1.9 cm), 1 each
50-1531	GREENBERG Tapered Retractor Blade, malleable, 0.63 in. (1.6 cm), 1 each
50-1532	GREENBERG Tapered Retractor Blade, malleable, 0.50 in. (1.3 cm), 3 each
50-1533	GREENBERG Tapered Retractor Blade, malleable, 0.38 in. (9.5 mm), 3 each
50-1534	GREENBERG Tapered Retractor Blade, malleable, 0.25 in. (6.4 mm), 3 each

Level III

To achieve Level III, add to the above components the following

50-1507	GREENBERG Primary Bar, 1 each
50-1514	GREENBERG Sterilizing Case, 1 each
50-1516	GREENBERG Floating Secondary Bar, 1 each
50-1523	GREENBERG Small Instrument Holder, 1 each
50-1524	GREENBERG Large Instrument Holder, 1 each
50-1525	GREENBERG Extra Large Instrument Holder, 1 each

Also available separately:

50-1515	GREENBERG Maxi-Vise Adapter
---------	-----------------------------

Service and Repair

For service or repair outside the United States, contact your local Symmetry Surgical representative.

Inside the United States, send instruments for service or repair to:

Symmetry Surgical
3034 Owen Drive,
Antioch, TN 37013 USA

Always include the purchase order number and a written description of the problem.

Warranty

FLEXBAR Arms are warranted for one (1) year from the date of purchase. GREENBERG Neurosurgical Retractor System components are guaranteed to be free from defects in both materials and workmanship. **Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed. Suitability for use of this medical device for any particular surgical procedure should be determined by the user in conformance with the manufacturer's instructions for use. There are no warranties that extend beyond the description on the face hereof.**

[™] GREENBERG is a trademark of I.M. Greenberg, M.D.

® PRESERVE is a registered trademark of Symmetry Surgical

[™] FLEXBAR is a trademark of Flexbar Machine Corp.

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Referenties

1. Greenberg, I.M.: "Self-Retaining Retractor and Handrest System for Neurosurgery," NEUROSURGERY, 8:205-208, 1981

2. Greenberg, I.M.: "New Options for the Neurosurgeon, Multiple Instru- mentation in a Single Microsurgical Field," NEUROSURGERY, 9:566-583, 1981

3. Greenberg, I.M.: "Staircase Concept of Instrument Placement in Micro- surgery," NEUROSURGERY, 9:696-702, 1981

4. Greenberg, I.M.: "Cerebral Aneurysm Rupture During Neurosurgery," NEUROSURGERY, 15:243-245, 1984

5. Greenberg, I.M.: "Rotation of Scissors While Fixed in the Microsurgical Field: A Technique for Increasing the Surgeon's Efficiency," MICROSUR- GERY, 5:218-227, 1984

6. Greenberg, I.M.: "Self-Retaining Retractors," in Wilkins, R.H., Renga- chary, S.S., (eds): NEUROSURGERY, New York, McGraw-Hill Book Co., 1985

7. Greenberg, I.M.: "Intraoperative Ultrasonography with a Cystoscope for the Biopsy of a Deep-Seated Brain Lesion: Case Study," NEUROSUR- GERY, 19:49-58, 1986

8. Greenberg, I.M.: "Proximity Storage of Symmetry Surgical Instruments in the Operative Field," NEUROSURGERY, 20:950-953, 1987

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