




<b>EN / English</b>	Instructions for Use   Non-Sterile Reusable Endocavity Needle Guides
<b>DE / Deutsch</b>	Gebrauchsanweisungen   Unsterile wiederverwendbare Nadelführungen für Endokavitäten
<b>FR / Français</b>	Mode d'emploi   Guides-aiguilles endocavitaires réutilisables non stériles
<b>BG / Български</b>	Инструкции за употреба   Нестерилни водачи за ендокавитарна игла за многократна употреба
<b>CS / Čeština</b>	Návod k použití   Nesterilní endokavitární vodiče jehel pro opakované použití
<b>DA / Dansk</b>	Brugsanvisning   Ikke-sterile, genanvendelige nålestyr til indre kropshulrum
<b>EL / Ελληνικά</b>	Οδηγίες χρήσης   Μη αποστειρωμένοι επαναχρησιμοποιήσιμοι οδηγοί βελόνας ενδοκοιλότητας
<b>ES / Español</b>	Instrucciones de uso   Guías de agujas endocavitarias no estériles reutilizables
<b>ET / Eesti Keel</b>	Kasutusjuhend   Mittesteriilsed korduskasutatavad õonesisesed nõelajuhihukud
<b>FI / Suomi</b>	Käyttöohje   Ei-steriilit, uudelleenkäytettävät ontelonsisäiset neulanohjaimet
<b>HR / Hrvatski</b>	Upute za upotrebu   Nesterilne višekratne vodilice igle za endokavitete
<b>HU / Magyar</b>	Használati utasítás   Nem steril, újrafelhasználható endokavitális tűvezető
<b>IT / Italiano</b>	Istruzioni per l'uso   Guida aghi endocavitari riutilizzabili non sterili
<b>KO / 한국어</b>	사용 설명서   비멸균 재사용 가능 내부 체강 바늘 가이드
<b>LT / Lietuvių k.</b>	Naudojimo instrukcijos   Nesterilūs daugkartiniai „Endocavity“ adatų kreipikliai
<b>LV / Latviešu</b>	Lietošanas instrukcija   Nesterilas atkārtoti lietojamas endokavitālo adatu vadotnes
<b>NL / Nederland</b>	Gebruiksaanwijzing   Niet-steriele, herbruikbare endocavitaire naaldgeleiders
<b>NO / Norsk</b>	Bruksanvisning   Ikke-sterile gjenbrukbare endokavitets nåleføring
<b>PL / Polski</b>	Instrukcja użycia   Niejałowe wielorazowe prowadnice igieł wewnątrzjamowych
<b>PT / Português</b>	Instruções de utilização   Guias reutilizáveis não esterilizadas para transdutor endocavitário
<b>RO / Română</b>	Instrucțiuni de utilizare   Ghidaje pentru ace de endocavitate reutilizabile nesterile
<b>RU / Русский</b>	Инструкции по применению   Нестерильные многоцветные направляющие для внутрисполостных игл
<b>SK / Slovensky</b>	Návod na na použitie   Nesterilné endokavitálne vodidlá ihel na opakované použitie
<b>SL / Slovenščina</b>	Navodila za uporabo   Nesterilna vodila za endokavitetne igle za večkratno uporabo
<b>SR / Srpski</b>	Uputstvo za upotrebu   Nesterilni vodiči igle za endo-šupljine, za višekratnu upotrebu
<b>SV / Svenska</b>	Brugsanvisning   Icke-sterila återanvändbara endokavitära nålledare
<b>TR / Türkçe</b>	Kullanım İçin Talimatlar   Steril Olmayan Yeniden Kullanılabilir Endokavite İğne Kılavuzları
<b>ZH / 中國人</b>	使用说明   非无菌可重复使用的腔内导针器。
<b>JA / 日本語</b>	使用手順   非滅菌の再利用可能なエンドキャビティニードルガイド

#### INTENDED PURPOSE/INTENDED USE

This is a medical device used for guiding a needle or catheter during a diagnostic ultrasound/imaging procedure to perform a precise needle placement to an anatomical target of human body, at specific depths.

#### INTENDED USER

Needle guides are used by healthcare professionals trained and qualified on using these types of devices.

#### ⚠ BEFORE STARTING PROCEDURE

- Prepare probe per system's operational manual.
- Verify needle guide fits properly and is securely attached to probe. Do not use if damaged or does not fit properly.
- Always verify needle follows guidelines prior to use.
- Use appropriate needle length to reach target area.
- Always use a probe cover to minimize any contaminants.

#### INSTRUCTIONS FOR USE

1. Place gel inside cover and/or on probe scanning surface. Place probe cover over probe. Secure cover with bands.
  2. Place needle guide (styles vary) onto probe as shown.
  3. Make sure needle guide is securely attached to probe (a).
  4. Place second probe cover over probe and needle guide.
  5. Secure cover with bands.
- After procedure, to remove the needle guide and/or probe cover follow the steps shown for mounting on probe, in reverse order.

#### DISPOSAL

After use, dispose of all single-use components as infectious waste.

#### MALFUNCTION

In case of malfunction or changes in performance of the needle guide that may affect safety, stop procedure, remove defective needle guide, and replace with an intact one.

*Notice: Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user is established.*

#### CLEANING, DISINFECTING & STERILIZING

##### REUSABLE NEEDLE GUIDES

Reusable needle guide is provided non-sterile. Clean/sterilize needle guide before each patient use, following recommended instructions.

##### CLEANING

Reusable devices must be thoroughly cleaned before they are subjected to the sterilization process.

- Remember, you cannot achieve sterilization or high-level disinfection unless the assembly is cleaned first.
1. Rinse excess soil from device.
  2. Prepare the enzymatic solution (like Enzol), following manufacturer instructions for proper dilution.
  3. Fully immerse the device in the detergent and allow soaking for 10 minutes.
  4. After the soak, brush device thoroughly using a soft bristle brush. Clean lumens and holes using an appropriate size brush. A syringe might be used to flush the hard-to-reach areas of the device.
  5. Thoroughly rinse and dry the device. Always use R/O (reverse osmosis) or DI (de-ionized) water for rinsing the device, tap water may recontaminate the reprocessed device. Use a clean soft cloth to dry the device.
  6. Thoroughly examine all surfaces to make certain that no visible bioburden remains. If any visible parts observed on the device, repeat cleaning and drying, until the ENTIRE device is clean.

##### RECOMMENDATIONS FOR HIGH LEVEL DISINFECTION/STERILIZATION OF PLASTIC OR ANODIZED ALUMINUM REUSABLES

After cleaning by using the above recommendations, soak for a minimum of 5 minutes in a neutral pH Enzol Enzymatic detergent such as made by Johnson & Johnson®.

Sterilize with Cidex Activated Dialdehyde solution such as made by Johnson & Johnson®, or any equivalent validated sterilization method adopted by the sterilization department for like devices.

##### STEAM STERILIZATION OF STAINLESS STEEL REUSABLES

After cleaning by using the above recommendations, steam sterilize (sterilizer type - gravity) in an autoclave at minimum 132 degrees Celsius (270 degrees Fahrenheit) for a minimum of 15 minutes. Dry the device for a minimum of 25 minutes. When sterilizing multiple instruments in one autoclave cycle, ensure that the sterilizer's maximum load is not exceeded.

#### WARNING

- Long narrow cannulations and blind holes require particular attention during cleaning.
  - Do not reuse or reprocess devices labeled for single-use.
  - After reprocessing the reusable device follow above instructions, wrap in CSR or equivalent, and store in designated location to ensure sterility is not compromised prior to reuse.
- #### INSPECTION AND FUNCTION TESTING
- Visually inspect for wear and damage such as cracks or breakage. If damage is evident, discontinue use of device.
  - For hinged instruments, check for smooth movement of hinge without excessive "play."
  - Locking mechanisms should be checked for action.

#### DEFECTIVE REUSABLE DEVICES

A defective device can be reused only if it is reconditioned under the responsibility of the manufacturer to comply with the general safety and performance requirements.

#### LIMITATIONS ON REPROCESSING

Repeated processing has minimal effect on these devices. End of life is normally determined by wear and damage due to use.

#### DISCLAIMER

The instructions provided above have been validated by the medical device manufacturer as being capable of preparing a medical device for re-use. It is the user's responsibility to qualify any deviations from recommended method of processing, and properly evaluate for effectiveness and potential adverse consequences.

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Po čiščenju z uporabo zgornjih priporočil namakajte najmanj 5 minut v pH nevtralnem encimskem detergentu Inzol, kakršnega proizvaja Johnson & Johnson ©. Sterilizirajte z raztopino dialdehida, aktivirane s Cidexom, kot jo na primer proizvaja Johnson & Johnson ©, ali s katerokoli enakovredno preverjeno metodo sterilizacije, ki jo je za podobne pripomočke sprejel oddelek za sterilizacijo.
Aspen Surgical Products, Inc.