Symmetry Surgical™ Skull Tong

Catalog No. 66-1025

Symmetry Surgical
3034 Owen Drive,
Antioch, TN 37013 USA
☎ 1-800-251-3000
Fax: 1-615-964-5566

EC REP

Symmetry Medical Sheffield & Cheltenham
Unit 1, Stella Way Malvern View Park
Bishops Cleeve,
Cheltenham, GL52 7DQ UK

LCN 173019-001/M
© 2002–2012 Symmetry Surgical
♻ Revised 07/12

ENGLISH

IMPORTANT INFORMATION

Please Read Before Use

Symmetry Surgical™ Skull Tong Catalog No. 66-1025

Description

The Symmetry Surgical Skull Tong, sold nonsterile, is used for application of skeletal traction to the cervical spine in children and adults. The tong can be applied under antiseptic rather than sterile conditions, eliminating the need to shave, incise, or drill. A threaded hole at each end accommodates the skull pin housing for advancement of the cone-shaped pins through the scalp to the outer table. The points are tilted in the direction of pull to prevent disengagement when traction is applied to the skull. One of the points has an enclosed spring which serves as a pressure indicator. A test bar furnished with the tong allows the spring-loaded point to be checked for proper action, and provides protection for the points when the tong is not in use.

Inspection

Inspect each instrument before and after use to detect imperfections. Do not use damaged instruments; send them for repair or replacement.

Indications

The Symmetry Surgical Skull Tong provides skeletal traction immobilization of the neck in the presence of a fracture or dislocation of the cervical spine. This device can also be used to apply traction in cervical spondylosis or to widen the interspace for insertion of the bone plug in anterior interbody fusion.

Contraindications

Any fracture of the skull that might be aggravated by the tong must be ruled out before application.

The tong should not be applied to extremely thin or diseased skull.

Infection of the scalp contraindicates use of the skull tong.

Compatibility of the Symmetry Surgical Skull Tong with magnetic resonance imaging (MRI) has not been determined.

WARNINGS

To prevent neurological injury, the patient's head must not be moved unnecessarily during or after attachment of the tong.

The tong must have a firm grip in the skull before the patient is placed in traction. Traction must be applied in the plane of articulating facets (in the long axis of the cervical spine).

Once traction has been achieved, the scalp must be watched carefully for signs of infection. If infection occurs at the application site, skull traction must be immediately placed in a new location or an alternate traction technique instituted.

Loosening of the pins may occur. Recheck the pins after 24 hours and tighten if necessary.

Take precautions to prevent the patient from unscrewing the points.

Complications

Disregarding instructions could result in the pins penetrating the inner table of the skull.

Osteomyelitis of the skull, cerebral infarctions, and intracerebral hemorrhage have all been reported following cranial penetration by other devices used in caliper skull traction.

Scalp, epidural, subdural, and intracerebral infections have been reported with the use of cranial tongs.

Epidural hematomas, cervical spinal injury, and fusiform intracranial middle cerebral aneurysm have occurred with skull traction.

Instructions for Use

Remove the pins from the tong and place in antiseptic solution.

Clean the scalp and spray the application sites with liquid antiseptic. Rub the antiseptic in at the temporal ridges and inject a local anesthetic. Once the anesthetic has taken effect, reapply the antiseptic.

Reassemble the pins in the tong and apply the pins just below the temporal ridges, which mark the greatest lateral diameter of the skull.

Advance the points alternately through the intact skin by turning the knurled knobs. The spring-loaded pin will compress upon entering bone and the distal end of the point will begin to protrude. Advance the pins until this end protrudes, becoming flush, 0.00 in. \pm 0.02 in. (0.00 mm \pm 0.5 mm) with the end of the knob, indicating that the spring is exerting approximately 30 pounds \pm 5 pounds (13.59 kg \pm 2.27 kg) of pressure between the pins.

Tilt the tong back and forth several times to ensure proper seating of the pins. Check protrusion of the spring-loaded pin and retighten if necessary. Apply traction weight and take precautions to prevent rotation of the head.

Recheck pins after 24 hours and retighten if necessary.

Reprocessing

Refer to Symmetry Surgical Instructions for Use 204233 for updated information regarding reprocessing the instrument. **SEE NOTE BELOW.**

Copies of this instruction sheet are available by requesting IFU 204233 from Symmetry Surgical Customer Service at:

Telephone: 800 251 3000
Fax: +1 615 964 5566
E-mail:
ssi-customerservice@symmetrysurgical.com

Note: The tong may be applied under antiseptic conditions if necessary by removing the pins and soaking them in liquid germicide until chemical sterilization is achieved.

Service and Repair

For service or repair outside the United States, contact your local Symmetry Surgical representative.

In the United States, send instruments for service or repair to:

Symmetry Surgical
3034 Owen Drive,
Antioch, TN 37013 USA

Always include a repair purchase order number and a written description of the problem.

Warranty

Symmetry Surgical warrants that this medical device is free from defects in both materials and workmanship. **Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed. Suitability for use of this medical device for any particular surgical procedure should be determined by the user in conformance with the manufacturer's instructions for use. There are no warranties that extend beyond the description on the face hereof.**

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