Bovie®

Bovie Medical Corporation
5115 Ulmerton Road
Clearwater, Florida 33760 USA

Emergo Europe
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands
SRN #: NL-AR-000000116
+31.70.345.8570
EmergoEurope@ul.com


Bovie®

Neuro-Pulse™ Surgical Nerve Locator

REF

CAT. 0003Y

R_x ONLY

STORAGE
10 to 75 % non-condensing
OPERATE
30 to 75 % non-condensing
TRANSPORT
10 to 100 % including condensation

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EN NEURO-PULSE™
SURGICAL NERVE LOCATOR

Store At Room Temperature

WARNINGS:

Do not use this device in the presence of flammable gases/materials or oxygen-enriched environments. Fire could result.

This medical device cannot be effectively cleaned and/or sterilized by the user and therefore cannot be safely reused. It is intended for single use only. Any attempt by the user to clean and resterilize this device may result in bioincompatibility, infection, or other risks of device failure to the patient.

Do not use this device at or around the heart or the central circulation system.

WARNING: This product is sterilized utilizing Ethylene Oxide. These products can expose you to chemicals including Ethylene Oxide (EtO), which is known to the state of California to cause cancer and/or birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

CAUTIONS:

Do not use if package is damaged or opened.

This is a single use device. Do not resterilize.

Use 1/2 mA initial setting. Limit nerve contact to one second maximum.

Instructions for the use of NEURO-PULSE™ Surgical Nerve Locator

Illustration

1. Remove the Neuro-Pulse™ (see Illustration) from the sterile package. Next, remove both the band which

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secures the needle wire and the piece of white foam from the probe.

2. **IMPORTANT:** Test the unit by touching the probe and the needle together. This will activate the LED at the distal end of the unit. If the LED is not activated, the unit is not functioning properly and should not be used.

Repeat this test after any change in the amperage (mA).

3. Insert the needle (the positive ground) into subcutaneous tissue near the proposed operative site. When the probe contacts tissue and the ground needle is in place, the circuit is completed. The LED indicator at the distal end of the instrument then lights indicating the instrument is operative.

IMPORTANT: The lighted LED indicator DOES NOT indicate a nerve has been located, only that the circuit has been completed. The muscle response is the indicator that a nerve has been located.

4. Select the amperage by rotating the probe index until the marker is near aligned with the desired output setting, and the engagement is felt. Use the 1/2 mA setting initially. If the muscle does not respond, increase the amperage output as required.

IMPORTANT: Muscle response to this stimulation is instantaneous. Nerve contact should not exceed one second. This instrument is not designed to determine the degree of viability of a nerve. Its intended use is for the identification of exposed motor nerves.

Medical Purpose / Indication

The intended use of the Nerve Locator is to identify the general location of exposed motor nerves while at the same time reducing the possibility of accidental nerve damage or severance.

Operating Principle

Nerve Locator is surgically invasive and is intended for transient use. The units are compact and portable and have three power settings (0.5 mA, 1.0 mA, 2.0 mA at 3.0 V). The yellow LED illuminates when the needle and probe comes in contact with tissue

to indicate the device is active. The LED indicator DOES NOT let the surgeon know when a nerve has been located, only that the circuit has been completed. The Muscle Response is the indicator that an exposed nerve has been located. These are single-use devices and are contraindicated for use at or around the heart or the central circulation system. They are further contraindicated for use in the presence of flammable gases and materials.

Patient Population

- No restriction

Site of Use and Site Conditions

- Site of use – No restrictions
- Site Condition – Aseptic

Intended User Profile

- Education
 - Trained medical professional
- Knowledge
 - Read and understand the supplied IFU
 - Understands hygiene
- Experience
 - Some training on techniques or training under surveillance/supervision
- Permissible Impairments
 - Mild reading/vision impairment or corrected vision to 20/20

Intended Conditions for Use

- Environment
 - Physician Procedure Room, Surgery center or Hospital, intended for professional use only
 - Sterile medical device
 - Single patient use
- Conditions of visibility
 - Clear and unobstructed view of point of use
- Frequency of Use
 - Single Use - Single patient multiple activation
- Mobility
 - Hand-Held

IMPORTANT: After use, use proper disposal method.

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