
Indications for Use: This product is a disposable electrosurgical pencil which is adaptable to common brands of electrosurgical generators for cutting and coagulation to remove tissue and control bleeding by use of high frequency current. This model can be utilized to remove smoke from the surgical field.

Figure A: Connection to Smoke Evacuator	Figure B: Connection to Electrosurgical Generator	Figure C: Removal of Blade Cover	Figure D: Removal of Electrode Blade	Figure E: Modes of Operation	Figure F: Smart Shark Remote Activator Device
<p>Remove Bovie® Orca™ pencil from packaging. When removing the device from the packaging, ensure that the device does not come into contact with any part of the non-sterile packaging. Connect the provided tubing to a surgical smoke evacuator or other compatible vacuum source. Be sure to avoid placing the tubing in areas through which personnel will travel.</p> 	<p>Plug the 3-prong power cord into the monopolar port on an electrosurgical generator. Prior to beginning a procedure, ensure that all power settings on the generator are at the appropriate levels. The intensity setting on the electrosurgical generator should be set to the lowest value that can achieve the desired result.</p> 	<p>The Orca pencil has a translucent blade cover that provides additional protection to both the patient and user during use. This cover can be removed by firmly grasping the shaft of the pencil and pulling the cover off. Avoid contact with the electrode and buttons when maneuvering the electrode blade cover. Build-up of eschar may hinder the device's ability to perform. To avoid injury, do not activate the device while cleaning the electrode.</p> 	<ol style="list-style-type: none"> Unplug the pencil from the electrosurgical generator. Use a surgical clamp to pull electrode from the socket. Insert replacement electrode into socket. Ensure electrode is entirely seated into socket prior to activating pencil. 	<p>The Orca pencil has 2 modes of operation. The modes of operation are CUT and COAGULATION. The YELLOW button on top of the pencil shaft activates the CUT mode and the BLUE button on top of the pencil shaft activates the COAGULATION mode. When finished with the device, turn off the smoke evacuation unit as well as the electrosurgical generator. Discard the pencil assembly (pencil, blade, tubing and power cord) per protocol for contaminated items, as it should now be considered a biological hazard.</p> 	<p>Bovie® Smart Shark Remote Activator Device can be used with the Orca pencil in order to maximize the smoke evacuator's filter life. The Remote Activator Device ensures that the smoke evacuation unit is activated only when the cut and coagulation buttons on the Orca pencil have been pressed and for a few seconds after their release.</p> 

IMPORTANT INFORMATION, Please Read Before Use: CONTRAINDICATIONS

DO NOT USE device in patients with electronic implants without first consulting a qualified professional in advance

⚠ CAUTIONS

CAUTION: Intended for use with a maximum voltage of 5kVp

CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a physician.

Inspect every device before use. Damaged device may cause electrical burns to the patient or physician. Active electrosurgical devices must be placed away from the patient when not in use. A holster is recommended.

Use a wet sponge to remove eschar build up from active electrode tips.

Do not perform electrosurgery in the presence of flammable anesthetics or other flammable gases, fluids or objects, or in the presence of oxidizing agents, as fire could result.

Avoid unnecessary or prolonged activation. Activate only when touching or immediately next to intended point of contact.

Electrodes must be compatible with the pencil and ensured that the insulated sleeve fits securely inside the handle. Improperly seated electrodes may cause burn or injury.

Do not bend, cut, or modify electrode before or during surgery.

Keep the voltage/power settings at the lowest possible level to achieve the desired effect.

After use, safely discard in accordance with established procedures for biological hazards

⚠ WARNINGS

Improper use of any element in electrosurgical systems may cause serious burn, shock, and fire or explosion hazard. The use of the device with safely functioning and compatible handles, cords, electrosurgical generators and accessories is the user's responsibility. Read and understand all precautions and instructions before attempting to use any electrosurgical devices. This product is not intended to be used for fluid removal. This product is for SINGLE USE ONLY; DO NOT RESTERILIZE. Use aseptic technique in all phases of handling. Symmetry Surgical single use devices have not been designed to undergo or withstand any form of alteration, such as disassembly, cleaning or resterilization after a single patient use. Reuse can potentially compromise device performance and any usage beyond the design intent of this single use device may result in unpredictable loss of functionality.

This product is sterilized utilizing Ethylene Oxide. These products can expose you to chemicals including Ethylene Oxide (EtO), which is known to the state of California to cause cancer and/or birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

Symmetry Surgical Inc.
3034 Owen Drive, Antioch, TN 37013, USA
USA & Canada 800.251.3000 Fax: 800.342.3272
customerservice@symmetriesurgical.com
International: +1615.964.5532
internationalorders@symmetriesurgical.com

CE 0197

REF ISEP1000

Zhejiang Shuyou Surgical Instrument Co., Ltd.
3689#, East Yangguang Avenue, Dipu street Anji
313300 Huzhou city, zhejiang province China

EC REP

Luxus Lebenswelt GmbH
Kochstr.1,47877,Willich,Germany
TEL: 0049-1715605732

MC- 23437 Rev.1

Do not resterilize

Do not reuse

EC REP

Authorized Representative in the European Community

Rx ONLY

Rx Only

STERILE **EO**

Sterilized using Ethylene Oxide

Manufacturer

Date of Manufacture: (YYYY-MM-DD)

CE Marking

QTY

Packaging Unit

LOT


Batch Code

REF


Catalogue Number

Consult Instructions for Use

Humidity Limitations

Temperature Limit

Do not use if package is damaged

Separate Collection

Use-By-Date: (YYYY-MM-DD)

This product contains no detectable Latex

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