

INSTRUCTIONS FOR USE



aspen surgical.

Sterile Needle Guide

Needle Guide Starter Kit

Director XT™ & RetroFit XT™

English

Intended Use

These are medical devices used for guiding a needle or catheter during a diagnostic ultrasound/imaging procedure to perform a precise needle placement to an anatomical target of human body, at specific depths.

Intended User/Patient Target Groups


Intended to be used by Healthcare Professionals on patients requiring medical procedures.

Instructions for Use

Prepare probe per system's operational manual.

Reusable Bracket

- Attach bracket (styles vary) to the probe, aligning needle guide receiver / locator indicator on the bracket (a) with the scanning side of the probe (b).
- Make sure bracket is securely attached to probe (c).

Multi-Angle Bracket

Changing Bracket Angle

Select bracket angle by pressing locator pad and moving to desired angle. Release pad to lock into position.

- Single-Use Needle Guide**
- Place gel inside cover and/or on probe scanning surface. Place probe cover (d) over probe and bracket assembly. Secure cover with bands.
 - Place needle guide (e) onto mounting bracket over probe cover as shown with locking pad in elevated position.
 - Lock guide into position by pushing down locking pad and up on lock button (f) until completely closed.
 - For a quick release of the needle, press needle guide pad (g) and remove needle from guide.
 - After use, to remove needle guide, push up locking pad and down on lock button (f).
 - The device is intended for continuous use for less than 60 minutes.
 - Do not use if package is damaged.

Caution

Always verify needle follows guidelines prior to use. Use Verification Kit when provided to verify the alignment of the needle with the system software line (refer to system's manual instructions). After use, dispose of probe cover and needle guide. Do NOT use on patients! If needle guide does not fit properly, do not use. Use appropriate needle length to reach target area.

Reusable bracket is provided non-sterile. Clean/sterilize bracket before use, following recommended instructions.

Warning

Single Use Devices Only: Reusing or reprocessing single use devices may cause damage to the device which in turn may cause unnecessary harm to the user and/or patient. Reuse or reprocessing is not recommended. Reusable devices are provided non-sterile and must be processed prior to use to include first use. Processing is defined on specific Instructions provided with each device.

Disposal

Dispose or recycle general surgical devices per facility protocol and all applicable federal, state, regional, and/or local laws and regulations. Notice: After use or upon determining product can no longer be used due to wear or damage, safely discard in accordance with established procedures for biohazardous waste.

Cleaning, Disinfecting & Sterilizing Reusable Brackets

Cleaning

Reusable devices must be thoroughly cleaned before they are subjected to the sterilization process. Remember, you cannot achieve sterilization or high-level disinfection unless the assembly is cleaned first.

- Rinse excess soil from device.
- Prepare the enzymatic solution (like Enzol), following manufacturer instructions for proper dilution.
- Fully immerse the device in the detergent and allow soaking for 10 minutes.
- After the soak, brush device thoroughly using a soft bristle brush. Clean lumens and holes using an appropriate size brush. A syringe might be used to flush the hard-to-reach areas of the device.
- Thoroughly rinse and dry the device. Always use R/O (reverse osmosis) or DI (de-ionized) water for rinsing the device, tap water may recontaminate the reprocessed device. Use a clean soft cloth to dry the device.
- Thoroughly examine all surfaces to make certain that no visible bioburden remains. If any visible parts observed on the device, repeat cleaning and drying, until the ENTIRE device is clean.

Recommendations for High Level Disinfection/Sterilization of Plastic or Anodized Aluminum Reusables

After cleaning by using the above recommendations, soak for a minimum of 5 minutes in a neutral pH Enzol Enzymatic detergent such as made by Johnson & Johnson®.

Sterilize with Cidex Activated Dialdehyde solution such as made by Johnson & Johnson®, or any equivalent validated sterilization method adopted by the sterilization department for like devices.

Steam Sterilization of Stainless Steel Reusables

After cleaning by using the above recommendations, steam sterilize (sterilizer type - gravity) in an autoclave at minimum 132 degrees Celsius (270 degrees Fahrenheit) for a minimum of 15 minutes. Dry the device for a minimum of 25 minutes. When sterilizing multiple instruments in one autoclave cycle, ensure that the sterilizer's maximum load is not exceeded.

Warning

- Long narrow cannulations and blind holes require particular attention during cleaning.
- After sterilization, appropriately package and store the device to ensure that sterility is not compromised prior to reuse.
- Always use a probe cover to minimize any contaminants.

Inspection and Function Testing

- Visually inspect for damage and wear. For hinged instruments, check for smooth movement of hinge without excessive play.
- Locking mechanisms should be checked for action.

Limitations on Reprocessing

Repeated processing has minimal effect on these devices. End of life is normally determined by wear and damage due to use.

Disclaimer

The instructions provided above have been validated by the medical device manufacturer as being capable of preparing a medical device for re-use. It is the user's responsibility to qualify any deviations from recommended method of processing, and properly evaluate for effectiveness and potential adverse consequences.

NOTICE: Any serious incident that has occurred in relation to the device should be reported to Aspen Surgical Products, Inc., the governing body, and the Competent authority of the member state in which the incident occurred.

Notice: For applicable IFU/additional information: Visit aspensurgical.com

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