Dr. Don's Arm Canoe and Suspension

Sterile

Dr. Don's® Arm Canoe 711500

Dr. Don's® Arm Canoe - V 711501

Dr. Don's® Arm Suspension Kit 711502

Dr. Don's® Arm Suspension Kit - Montana 711503

Dr. Don's® Arm Suspension Kit - Northern 711504

Aspen Surgical Products, Inc.
6945 Southbelt Drive SE
Caledonia, MI 49316 USA
Phone +1 616-698-7100
Toll-Free 888-364-7004
aspensurgical.com

Emergo Europe
Westervoortsejdijk 60
6827 AT Arnhem
The Netherlands

For Country of Origin See Product Label

English

Intended Use

Intended to support and maintain the position of an anesthetized patient during surgery.

Intended User/Patient Target Groups

Intended to be used by Healthcare Professionals on patients requiring medical procedures.

Instruction for Use:

- Present sterile wrapped goods to sterile field following appropriate sterile technique.
- Unwrap sterile CSR wrap.
- Place arm in foam holder after it has been aseptically prepped for surgery.
- Wrap with self-adhering wrap tightly.
- Hook to the distraction device and apply the manufacturer's recommended weight to the device for optimum shoulder distraction.
- Device is intended for limited contact duration, less than or equal to 8 hours, less if damaged.
- Discard after use.

Precautions:

- For use with the ARTHREX® Lateral Decubitus Shoulder Traction Tower and ARTHREX® 3-Point Shoulder Distraction System to achieve lateral position distraction allowing joint access.
- When hooking to the positioning device ensure weight is maintained on the device to avoid the arm from dropping.
- Ensure wrap is properly secured to the arm to avoid possible slippage and contamination of the sterile field.
- Single Use Only (intended to be used on one patient during one procedure). Dispose after use. Do not resterilize.
- Do not use if product is damaged.

Warning

Reuse of device could result in infection/contamination and/or device failure which could lead to patient harm.

Disposal

Dispose or recycle general surgical devices per facility protocol and all applicable federal, state, regional and/or local laws and regulations.

Notice to Users and/or Patients in EU

Any serious incident that has occurred in relation to the device should be reported to Aspen Surgical Inc., the governing body, and the Competent authority of the member state in which the incident occurred.

Arthrex® is a registered trademark of Arthrex®, Inc. ("Arthrex®"). Aspen Surgical Products, Inc. is not affiliated with, endorsed by, or sponsored by Arthrex®.

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Emergo Consulting (UK) Limited