




Transpak®



Symmetry Surgical Inc.
3034 Owen Drive
Antioch, TN 37013 USA
☎ 1-800-251-3000
Fax: 1-615-964-5566
www.symmetrysurgical.com



 **symmetry surgical**
Symmetry Surgical GmbH
78532 Tuttlingen, Germany
Maybachstraße 10
☎ +49 7461 96490
Fax: +49 7461 77921



These instructions pertain to Class I devices

IFU-LCN SYM8600-001 Rev A
© 2017 Symmetry Surgical



INSTRUCTIONS FOR USE

Intended Use

The TransPak® case is intended for use as a soiled surgical instrument point of use and transport vessel.

Indications for Use

None.

Product Description

The TransPak® case is designed to reduce leakage, minimize potential for injury and exposure to potentially infectious materials or fluids and help prevent damage to soiled surgical instruments during transport.

 **WARNINGS**

NOT INTENDED FOR USE IN AUTOCLAVES. KEEP AWAY FROM HOT SURFACES.

DIRECTIONS FOR USE

1. Wear appropriate Personal Protective Equipment (PPE), as required by facility policies, when handling soiled instruments.
2. Place soiled instruments inside TransPak®.
3. Keep instruments moist until they are cleaned. A towel moistened with water placed over the instruments may be used. Instruments that cannot be cleaned immediately may be treated with an instrument cleaner according to the device and the instrument cleaner manufacturers' written Instructions For Use (IFU). Liquids used to soak contaminated items at the point of use should be discarded before transport.
4. Close and latch the lid to prevent spillage, contamination or exposure to potential hazards.
5. Place TransPak® on the designated cart for transport to the instrument decontamination area.
6. Use caution during transport to minimize spilling and other potential hazards.
7. Transpak may be cleaned manually or in washer/decontaminator cycle suitable for processing plastics.

CLEANING

TransPak® cases are compatible with enzymatic and nonenzymatic detergents and the below list of FDA cleared liquid chemical germicides (LCG) for manual reprocessing medical devices. Always read and follow both the medical device manufacturer and LCG manufacturer's reprocessing instructions.

| | |
|-------------------------------------|--------------------------|
| Cidex Plus | ASP/J&J |
| MedSci 3% Glutaraldehyde | MedSci, Inc. |
| Metricide | Metrex Research, Inc. |
| Metricide Plus 30 Long-Life | Metrex Research, Inc. |
| Metricide 28 Long-Life | Metrex Research, Inc. |
| Omnicide Long Life | Cottrell Limited |
| Omicide Plus | Cottrell Limited |
| Procide 14 | Cottrell Limited |
| Sporicidin | Sporicidin International |
| Wavicide | Wave Energy Systems |
| <i>Ortho-phthaldehyde solutions</i> | |
| Cidex OPA Solution | ASP/J&J |
| Opaciden Solution | Ciden Technologies, LLC |
| <i>Hydrogen Peroxide Solutions</i> | |
| Acecide | Minntech Corp |
| EndoSpox Plus | Cottrell Limited |
| Peract 20 | Minntech Corp |
| Resert XL | Steris Corporation |
| Sporox | Reckitt & Colman |

WARRANTY

Symmetry Surgical warrants that this medical device is free from defects in both materials and workmanship for one (1) year from the date of purchase. **Any other express or implied warranties, including warranties of merchantability or fitness for a particular purpose, are hereby disclaimed. Suitability for use of this medical device for any particular surgical procedure should be determined by the user in conformance with the manufacturer's instructions for use. There are no warranties that extend beyond the description on the face hereof. Abuse or misuse of the product or failure to comply with the instructions for use shall void this warranty.**

Copies of the IFU are available by requesting them from Symmetry Surgical Customer Service at:

Phone: 1-800-251-3000

Email: customerservice@symmetrysurgical.com

Website: www.symmetrysurgical.com

| Product Name | Manufacturer |
|---------------------------------|---------------------------|
| <i>Glutaraldehyde solutions</i> | |
| Aldahol III | Healthpoint LTD |
| Banicide Advanced | Pascal Company, Inc. |
| Cetylclde-G | Cetylite Industries, Inc. |
| Cidex Activated Dialdehyde | ASP/J&J |
| Cidex Formula 7 | ASP/J&J |

LABELING SYMBOLS



Symbol for Manufacturer; the name and address of the Manufacturer is next to this symbol



Caution-Consult accompanying documents




Lot number or batch number



Catalog number



Consult instructions for use



European Representative



In USA contact

Symmetry Surgical Inc.
3034 Owen Drive
Antioch, TN 37013 USA
Phone: 1-800-251-3000
Fax: 1-615-964-5566
www.symmetrysurgical.com



In Europe contact

Symmetry Surgical GmbH
Maybachstraße 10
78532 Tuttlingen, Germany
Phone: +49 7461 96490
Fax: +49 7461 77921