<p>For product information, contact:</p> <p>Yderligere informationsopsyninger fås ved at kontakte: For obtenir des informations sur le produit, veuillez contacter :</p> <p> Für Produktinformationen kontaktieren Sie:</p> <p>Για πληροφορίες σχετικά με το προϊόν, επικοινωνήστε με:</p> <p>Per informazioni sul prodotto, contattate:</p> <p>製品情報、お問い合わせ先:</p> <p>Para obter informações sobre o produto, contacte:</p> <p>Para obtener información del producto, póngase en contacto con:</p>	<p>symmetry olsen[®]</p>
<p>3230 Commerce Center Place Louisville, Kentucky, 40211 (502) 813-5022 ⇔ (800) 297-6344</p> <p>See our products at Se vores produkter på Découvrez nos produits sur Siehe unsere Produkte auf Αιτή το προϊόν μας στην ιστοσελίδα</p> <p>È possibile vedere i nostri prodotti all'indirizzo web 当社の製品は以下でご覧いただけます Consulte os nossos produtos em Consulte nuestros productos en</p>	
<p>http://www.symmetryurgical.com</p>	

R ONLY

Reusable Irrigating Bipolar Forceps Genanvendelig bipolar skylletang Pinces bipolaires d'irrigation réutilisables Wiederverwendbare bipolare Spülungsfassinstrumente Επαναχρησιμοποίησημη Διπολική Λαβίδα Καταιόνησης Pinze bipolari ad irrigazione riutilizzabili 再使用可能洗浄用バイポーラ鉗子 Fórceps de irrigação bipolar reutilizáveis Pinzas bipolares de irrigación reutilizables

ENGLISH

Description:

Reusable devices are provided non-sterile and must be processed prior to use to include first use. Processing is defined on specific Instructions provided with each device. Each reusable device is designed to withstand up to 20 sterilization cycles when properly cared for and sterilized according to the instructions below. The irrigating tube facilitates local fluid irrigation through the forceps.

Intended Use (per EU definition) / Indications (per FDA definition):
Olsen Medical® Irrigating Bipolar Forceps are intended for use in coagulating and irrigating tissue during surgical procedures and are intended for use with a maximum voltage of 500 Volts.

Handling:

- Remove and discard the tip protector on forceps before use.
- Inspect device before each use. Discard or send for repair if damaged.
- Device use and application to be determined by physician.
- Clean and decontaminate the device according to the instructions below:
- Begin decontamination immediately after a surgical procedure. Wipe off all debris with a moist sponge or gauze pad and flush lumen with water. Do not use saline.
- Clean devices with mild detergent (Tergazyme or a similar type solution) that is safe for plastics.
- Completely immerse the Irrigating Bipolar Forceps in Ultrasonic Cleaner following the recommendations of the ultrasonic manufacturer regarding cycle times, detergents, proper placement of the instruments, etc.
- Rinse thoroughly with running tap water to remove any debris or cleaning solution.
- Use a syringe or pressurized water to rinse the irrigation tube and verify that fluid flow is unobstructed. If obstruction exists use a stylet to clean out the lumen or repeat the ultrasonic cycle. Blow compressed air through the lumen to dry.
- Wipe dry and sterilize according to the instructions below.

Storage:

The Reusable Irrigating Bipolar Forceps should be stored individually in a protective container in a clean, dry, moisture free area with suitable protection from extreme temperature and humidity. Keep the stylet inserted in the lumen for non-sterile storage.

Sterility:

This device is sold as non-sterile and must be steam sterilized prior to each use. Other methods of sterilization are not recommended and their use will void the warranty.

Steam Sterilization:

Sterilization Cycle	Temperature	Time	Minimum Dry Time	Device Configuration
Prevacuum	270°F (132°C)	4 Minutes	20 Minutes	Wrapped/ Individually Pouched
Gravity Displacement	270°F (132°C)	10 Minutes	10 Minutes	Wrapped

Notices:

- Electrosurgical equipment and accessories are intended to be used by health professional educated in their use.
- After Use or upon determining product can no longer be used due to wear or damage, safely discard in accordance with established procedures for biohazardous waste.
- Any serious incident that has occurred in relation to the device should be reported to Symmetry Surgical Inc., the governing body, and the Competent authority of the member state in which the incident occurred.
- Please refer to the manufacturer of the generator for warnings, precautions, contraindications, undesirable side-effects, measures to be taken, and limitations of use for the electrosurgical system and accessories
- For applicable IFU/additional information: Visit www.symmetryurgical.com.

Caution:

- For effective sterilization of the irrigating bipolar forceps remove the stylet from the irrigation tube.
- Ensure that the Irrigation Tubing and Bipolar Cord set is properly connected to the irrigating bipolar forceps before turning on the power to the generator.
- Inspect every device before each use. Damaged electrosurgical devices may cause electrical burns to the patient or physician.
- Active electrosurgical devices must be placed away from the patient when not in use. A holster is recommended.
- Use a wet sponge to remove eschar build up from the forceps tips. Use a stylet as needed to keep the lumen free of debris.
- Consult RF generator Owner’s Manual. Use only in the active RF generator output connections.
- Do not perform electrosurgery in the presence of flammable anesthetics or other flammable gases, fluids or objects, or in the presence of oxidizing agents, as fire could result.
- Avoid unnecessary or prolonged activation. Activate only when touching or immediately next to the intended point of contact.
- Do not bend, cut, or modify any electrosurgical device before or during surgery.
- Keep the voltage/power settings at the lowest possible level to achieve the desired effect.
- Establish a protocol for device examination and maintenance.

Warning:

Improper use of any element in electrosurgical systems may cause serious burn. The use of the device with safely functioning and compatible handles, cords, electrosurgical generators and accessories is the user’s responsibility. Read and understand all precautions and instructions before using any electrosurgical device.

Warning:

Avoid HF output mode/power settings where maximum output voltage may exceed the stated maximum voltage.

Warning:

Do not make any modification to these devices. Modification may cause burn or injury.

Warning:

Due to concerns about the carcinogenic and infectious potential of electrosurgical byproducts (such as tissue smoke plume and aerosols which may include toxic gasses and vapors, live and dead cellular material, and viruses), protective eyewear, filtration masks, and effective smoke evacuation equipment should be used in both open and laparoscopic procedures. Contact Symmetry Surgical at customerservice@symmetryurgical.com for additional information or to inquire about our smoke evacuation solutions.

Warranty:

All products are warranted to be free from defects in workmanship and materials when used for its intended purpose. Any product that proves defective in workmanship or materials will be repaired or replaced at our discretion free of charge. Repair, modification, or alteration to any product that is performed by anyone not authorized by Olsen Medical® will result in the immediate loss of warranty.

OLSEN® is a registered trademark of Olsen Medical.

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