INSTRUCTIONS FOR USE: AntiFogOR™ AntiFog Solution for Endoscopic Optics

EN

INDICATIONS FOR USE / CLINICAL BENEFITS

The AntiFogOR is a sterile, single-use, disposable medical device intended to prevent condensation on the distal lens of endoscopes and laparoscopes. The AntiFogOR consists of a 6ml vial of solution and a sponge with adhesive on the back.

CONTRAINDICATIONS

This device is not intended for use except as indicated.

CAUTIONS / WARNINGS

- AntiFogOR was designed, tested and manufactured for single patient use only.
- Reuse or reprocessing of the device may lead to its failure and subsequent patient injury.
- Reprocessing/re-sterilization of this device may create the risk of contamination and patient infection.
- Upon completion of the procedure, AntiFogOR should be discarded in biohazardous waste as per standard procedure.

INSTRUCTIONS FOR USE

1. Open the sterile pouch, remove the AntiFogOR bottle and place it on the sterile field.
2. Attach the sponge to the patient drape or instrument table.
3. Open the AntiFogOR bottle and wet the sponge with the solution.
4. Press the distal tip of the endoscope onto the sponge and carefully turn it to apply a membrane of the AntiFogOR solution.
5. Remove the endoscope from the sponge.

Any serious incident that has occurred in relation to this device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Figure 1

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Şekil 1

	Product conforms with requirements of European Medical Device Regulation 2017/745.		Caution: US Federal Law restricts this device to sale by or on the order of a physician.
	Caution, consult accompanying documents.		Single Sterile Barrier System.
	Do not re-sterilize.		Do not use if package is damaged.
	Keep away from sunlight.		Consult Instructions for Use.
	Keep dry.		Not made with natural rubber latex.
	Sterilized using irradiation.		Do not reuse.
	Single sterile barrier system with protective packaging outside.		Medical Device.
	Use-by date.		Store between 15°C to 25°C temperature.
Trademark protected.			

Distributed exclusively in USA by: Symmetry Surgical, 3034 Owen Drive, Antioch, TN 37013, USA
Email: customerservice@symmetrysurgical.com, Telephone: USA & CANADA +1 800-251-3000 (Monday-Friday, 7:00AM - 5:30PM CST)

The O R Company Pty Ltd, 4/47 Wangara Road, Cheltenham, VIC, 3192, AUSTRALIA.
www.theorcompany.com

EC REP OBELIS SA, Bd Général Wahis 53, 1030 Brussels, BELGIUM,
Tel: +(32) 2 7325954

Issued date: 22-07-2022

PM_PS003_101_V7