# PretzelFlex™

Part of the EndoFlex™ Range

Surgical Retraction System

EN Instructions for use

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







16.

17.

Description	Part no.
5mm PretzelFlex Retractor	SHR7001
5mm Sterilisation Sleeve	SHR7101
3mm PretzelFlex Retractor	SHR7301
3mm Sterilisation Sleeve	SHR7331

### INDICATIONS FOR USE

This device is designed as an organ and tissue retractor for use in minimally invasive surgical procedures to elevate or retract organs and tissues to provide better access as well as visualization of surgical sites.

### CONTRAINDICATIONS

None Known.

### WARNINGS

Surgical Innovations Ltd devices are to be used in accordance with these instructions for use. Read all sections of this insert prior to use. Improper use of the instrument may cause a serious injury. In addition, improper care and maintenance of the instrument may render the instrument non-sterile and/or damaged prior to patient use and cause a serious injury to the patient or health care provider.

Do not clean device using an automated cleaning system. Clean device manually as specified within this Instruction for Use.

### CAUTIONS

If there are any variations between these instructions and either your facility's policies and/or your cleaning/sterilisation equipment manufacturer's instructions, those variations should be brought to the attention of the appropriate responsible hospital personnel for resolution before proceeding with cleaning and sterilizing your devices.

Use of the device for a task other than that for which it is intended will usually result in a damaged or broken device.

Prior to use, inspect device to ensure proper function and condition. Do not use devices if they do not satisfactorily perform their intended function or have physical damage.

Avoid mechanical shock or oversteering the devices.

It is important when using the retractor to move organs or tissue that the shape and load is supported evenly by the shaft and predefined platform shape of the device. Place the load in the orientation and shape in the diagram. Refer to Figure 16. Incorrect orientation could result in instrument damage and patient injury.

When clamping the 3mm device to a bench clamping system ensure the device is clamped in the clamp area shown in Figure 1 and Figure 17.

Incorrect clamping in any other area could result in instrument damage and patient injury.

Always use caution when inserting or removing devices through cannula. Lateral pressure on the device during removal can damage the working tip, shaft of the device.

Use a neutral pH (6-8) solution when cleaning the instrument. Use of highly acidic or alkaline detergents will cause pitting or breakage of the instruments.

The retractor is not entirely gas tight. To minimize gas leakage through the instrument, always ensure the luer fitting is covered and that the instrument is not left un-actuated in the patient.

The warranty may be void if the instrument is continually exposed to high acid or high alkaline solutions, detergents, stain removers etc.

Do not use instruments if they do not perform satisfactorily in an operational test. Avoid mechanical shock or oversteering the instruments which will cause damage.

Do not dismantle the device.

When sterilizing or storing this device, always use the protective/sterilizing sleeve provided. Failure to use the sleeve may result in premature device failure. This device should never be folded or bent to fit into a small sterilisation tray.

Only the cleaning and sterilisation processes which are defined within these instructions for use have been validated.

### Examples

1. Sterilisation using any other method than stated may damage the instrument.
2. Do not place the instrument in direct contact with live electrical instruments.
3. During use it is possible to generate considerable load on the instrument, which can under extreme circumstances break. **Do not over tighten the actuation knob, once the retractor has assumed its actuated shape only rotate the knob a maximum of 1/4 turn (90°) in clockwise direction.**
4. In the event of a breakage it is important to carefully remove the instrument and carefully examine to ensure that all components, specifically the actuating segments are accounted for.

### COMPATIBILITY

5mm retractors are designed for use with 5mm laparoscopic port access systems with 5.5mm minimum internal diameter.

3mm retractors are designed for use with 3mm laparoscopic port access systems with 3.5mm minimum internal diameter.

Both 3mm & 5mm devices are compatible with the Fastclamp Endoscopic Clamping System.

### HOW SUPPLIED

The device is packaged non-sterile. Cleaning and sterilisation must occur prior to use.

### LIMITATIONS ON REPROCESSING

The life of the device is validated by the manufacturer to be 100 reprocessing cycles, subject to cleaning process and inspection/maintenance criteria in re-processing section below.

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## RE-PROCESSING

### PRE-PROCESSING INSTRUCTIONS

1. Initiate cleaning of device within 2 hours of use.
2. Transport devices via the institutions established transport procedure.
3. Remove excess gross soil as soon as possible after use by rinsing the device in 6 gallons of tap for a minimum of 3 minutes. All devices must be processed in the relaxed straight configuration. All flush ports shall remain in the fully open position. Refer to Figure 2

### MANUAL CLEANING

1. Ensure all pre-processing instructions are followed prior to cleaning.
2. Prepare the enzymatic / neutral pH detergent solution, per manufacturer's instructions, utilizing tap water with a minimum temperature of 36°C (97°F). **Note:** This detergent preparation shall be used to complete cleaning stages 3, 4 and 5. It is recommended that the detergent solution is changed when it becomes grossly contaminated (bloody and/or turbid).
3. Place device in the open/relaxed position, with flush port open, and completely immerse in the detergent solution and allow device to soak for a minimum of 20 minutes. Actuate all movable parts during the initiation of the soak time. Refer to Figure 3
4. Using a small firm bristled brush, remove all visible soil from the device. Actuate device while brushing, paying particular attention to hinges, crevices and other difficult to clean areas. Refer to Figure 4
5. Suck or flush the internal lumen of the instrument with neutral pH detergent by using the flushing port located on the handle. A minimum of 3 complete flushes with 50ml of neutral pH detergent or a sufficient volume of detergent and number of flushes should be used until all fluid exiting the instrument is clear of particulate matter. Refer to Figure 5 (same as Figure 2)
6. If visible soil is detected during the final lumen flush, re-perform brushing and flushing of the lumen.
7. Rinse the device by completely immersing in 4 gallons of tap water with a temperature range of 27°C to 44°C (81°F to 111°F), for a minimum of 2 minutes to remove any residual detergent or debris. Refer to Figure 6
8. Place the device in an ultrasonic wash with enzymatic / neutral pH detergent solution (prepared as per the manufacturer's instruction) for a minimum of 10 minutes, at a minimum temperature of 39°C (102°F). Ensure sufficient volume of detergent solution is prepared to fully submerge the device. Refer to Figure 7
9. Following the ultrasonic wash, the device shall be thoroughly rinsed internally and externally with a minimum of 4 gallons of osmosis/deionized (RODI) water at room or ambient temperature for a minimum of 2 minutes. The internal flushing of the lumen should be through the flushing port located on the handle/shaft. In both instances internal and external flushing of the device should be repeated a minimum of 3 times ensuring all fluid exiting the device runs clear. Refer to Figure 8
10. Dry the device with a clean, lint-free towel.
11. For lumen devices, manipulate the device to allow rinse water to drain from the lumen.
12. Visually examine each device for cleanliness.
13. If visible soil remains, repeat cleaning procedure.
14. Before sterilising, lubricate the device with instrument milk or a steam permeable/ water soluble lubricant, following the lubricant manufacturer's instructions. Refer to Figure 9 Before sterilising the 3mm device pay particular attention to ensuring lubricant is directed into the spring clutch inside the actuation knob as shown in Figure 9.
15. Let devices drip dry for 3 minutes before packaging for sterilisation.

### INSPECTION/MAINTENANCE

Proper care and handling is essential for satisfactory performance of any surgical device. The previous cautions should be taken to ensure long and trouble-free service from all your surgical devices. Inspect devices before each use for broken, cracked, tarnished surfaces, movement of hinges, and chipped or worn parts. If any of these conditions appear, do not use the device. Return devices to an authorized repair service centre for repair or replacement.

### PACKAGING

Devices can be loaded into dedicated packaging systems. Sterilisation wrap material must be cleared for the applicable sterilisation modality by your country's regulatory body. Use in accordance with packaging manufacturer's sterilisation instructions.

Wrap devices in accordance with local protocols and procedures (See sample preparation, section Pre-vacuum Steam Sterilisation Parameters below for details of wrap specification).

### STERILISATION (REFER TO FIGURE 10)

All devices must be processed in the completely open position (i.e. flushports open) to allow sterilitant contact of all surfaces.

#### Parameters for Steam Sterilisation (Autoclave)

Sterilisation of instruments may be accomplished by steam autoclave. Time and temperature parameters required to steam sterilisation vary according to type of sterilizer, cycle design and packing material. The following are the manufacturer's guidelines for product sterilisation.

The use of "flash" sterilisation is not recommended, as it will shorten the life of endoscopic instruments.

Do not sterilize instruments at temperatures over 141°C (285°F)

The Pretzel retractors must only be sterilized in the fully relaxed position. The Flush ports should be un-obstructed by luer cap. Before sterilisation the "protective/sterilizing sleeve" must be located on the flexible, segmented portion of the instrument. Failure to follow these instructions may significantly impair the life of the retractor.

Note: Care must be taken to protect the flexible portion from damage, do not manually bend the segment string.

#### Pre-vacuum (Dynamic Air Removal) Steam Sterilisation Parameters:

Preconditioning Pulses: 3

Temperature: Minimum 132°C (270°F)

Exposure Time: Minimum 3 minutes

Dry Time: 30 minutes

Sample Preparation: Individually wrapped in two layers of 1-ply polypropylene wrap (Kimguard KC600 or equivalent using sequential square wrapping technique).

### STORAGE

After sterilisation, instruments should remain in sterilisation wrap and be stored in a clean, dry cabinet or storage case. Always store the retractor in its protective/sterilisation sleeve.

Following sterilisation, the shelf life of the sterile wrapped device must be determined by referring to the wrap manufacturer's guidance on shelf life with respect to sterility maintenance.

## OPERATING INSTRUCTIONS

### INSTRUMENT INSERTION

1. Ensure instrument is cleaned and sterilized in accordance with cleaning and sterilisation guidelines prior to use.
2. Insufflate the peritoneal cavity in the normal manner. Carefully guide the fully relaxed instrument through a laparoscopic port into the operating site. Refer to Figure 11
- 3.a. For 5mm Device rotate the actuation knob clockwise and form to the pre-defined shape. Ensure that the tip of the retractor and the segments are always in the field of endoscopic vision during formation. **DO NOT over tighten the actuation knob. Once the retractor has assumed its actuated shape only rotate knob maximum of ¼ turn (90°) in clockwise direction.** In the unlikely event of tissue becoming "snagged" between the segments relax the retractor, rotate slightly and re-activate. Refer to Figure 12
- 3.b. For 3mm Device rotate the actuation knob clockwise and form to the pre-defined shape. Ensure that the tip of the retractor and the segments are always in the field of endoscopic vision during formation. A slipping clutch is included within the actuation knob to prevent over tightening. The Actuation knob will slip and make an audible clicking when the correct unrestricted tension is achieved. Refer to Figure 12.
4. When actuating in limited space, insert section A of the flexible portion and actuate to approximately 90°. Refer to Figure 13
5. Insert the section B and further actuate the instrument ensure that the tip remains approximately 2cm away from the port. Refer to Figure 14
6. Insert section C and actuate the instrument until fully formed, ensure that all segments are inserted prior to forming. Take care to ensure the port does not bow during formation, if this occurs slacken the instrument and insert more segments. Refer to Figure 15
7. When using the retractor to move organs and tissue ensure that the shape and load is supported by the shaft of the device. Place the load on the orientation and shape as shown in the diagram. Refer to Figure 16
8. When instrument is in the position required, clamp to the operating table using an appropriate clamping instrument, such as Surgical Innovations Ltd "Fast Clamp". Ensure the clamp is fixed to the shaft of the device. When clamping the 3mm device ensure the device is clamped onto the larger diameter area next to the handle. Refer to Figure 1 and Figure 17.
9. DO NOT throw away the sterilisation tube.

### INSTRUMENT REMOVAL

1. Unclamp the instrument from the table and turn the actuation knob anti-clockwise to the fully relaxed position. Ensure the direction of removal is in line with the port.
2. Carefully withdraw the instrument from the port. If there is resistance during this process gently rotate the instrument during withdrawal. Do not repeatedly insert and withdraw the instrument as this may cause damage to the port and valve system
3. Rinse the device from visible soil immediately after use or within a maximum of 2 hours.

## WARRANTY

The Surgical Innovations Ltd device carries a warranty to be free from defects in material and workmanship. Any product which proves to be defective in material or workmanship will be repaired or replaced at our discretion. This warranty does not cover damage to the instrument by overstress or mechanical shock.

Repair, alteration or modification of any product by persons other than Surgical Innovations Ltd, or products subjected to misuse or abuse will result in the immediate loss of warranty. If Surgical Innovations Ltd instruments are damaged by accident or when used for a purpose other than that originally intended, a repair charge will apply.

Surgical Innovations Ltd offers a 2 year or 100 cycles guarantee, whichever is reached first.

## REPAIR SERVICE

If the device requires a service, return it to your authorised Surgical Innovations Ltd representative or the address stated. If the repair is covered under warranty, it will be repaired or replaced at no charge when requested in writing. A nominal service charge will be made for repaired instruments outside the warranty.

**Note:** The warranty is void if the instrument is repaired by any other repair service, subjected to abuse, or improperly reprocessed.

**Note:** All devices being returned for maintenance, repair, etc. must be cleaned and sterilised per these instructions prior to shipment

### CONTACT INFORMATION:

For repair services and/or advice for use concerning this instrument, please contact your local Surgical Innovations representative, or use the contact address on the Instruction for Use.

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