INSTRUCTIONS FOR USE

Disposable Light
Handle Covers
Non-Sterile

aspensurgical.

MD R_X ONLY

EC REP

Aspen Surgical Products Inc.
6945 Southbelt Drive SE
Caledonia, MI 49316
Phone +1 616-698-7100
Toll-Free 888-364-7004
aspensurgical.com

Emergo Europe
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

For Country of Origin See Product Label

English

Intended Use

Intended to provide a sterile interface between the operating room surgical lamps (non-sterile) and the Surgeon or Nurse (sterile).

Intended User/Patient Target Groups

Device is non-patient contacting and is to be used by Healthcare Professionals.

Provides a means of manipulating O.R. Lights.

Instructions for Use:

- Verify securement of the light handle cover to the adaptor.
- Do not use if product appears visually damaged or degraded.
- Device is intended for use less than or equal to 16 hours.

Sterilization Instructions:

- Place device into a kit or sterilization tray.
- Using Ethylene Oxide sterilization method, validate according to ISO 11135, sterilize according to the following parameters:

Cycle Parameter	Specification	Duration
Preconditioning		
Temperature	110°F +/- 10°F	
Humidity	65% +/- 15% RH	10.35 Hrs - 72 Hrs
Chamber		
Vacuum	21" +/- 1.0 inHgA	
Humidification/Dwell	22.5" +/- 1.0 HgA	60 +30/-0 minutes
Gas Inject/Dwell	10.3" +/- 0.5 Hg rise	360 +5/-0 minutes
Jacket Temperature	135° +/- 5° F	
Aeration*		
Temperature	120° +/- 10° F	24 Hrs - 72 Hrs Max

*Proper aeration is essential to eliminating the risk of residuals to meet ISO 10993-7 requirements.

Warning

Reuse of device could result in infection/contamination and/or device failure which could lead to patient harm.

Disposal

Dispose or recycle general surgical devices per facility protocol and all applicable federal, state, regional, and/or local laws and regulations.

Notice to Users and/or Patients in EU

Any serious incident that has occurred in relation to the device, should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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Emergo Consulting (UK) Limited