For product information, contact:

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### ENGLISH

**Description:**
Reusable devices are provided non-sterile and must be processed prior to use to include first use. Processing is defined on specific Instructions provided with each device. Each reusable device is designed to withstand up to 20 sterilization cycles when properly cared for and sterilized according to the instructions below. These devices are designed to function as active electrodes or as accessories to active electrodes. This family of devices includes, but is not limited to: monopolar and bipolar forces, pens/pencils and electrode handles, electrodes, cords and adapters.

**Intended Use (per EU definition) / Indications (per FDA definition):**
Monopolar Cord and Adapters - Olsen Medical® electrosurgical cords and adapters are designed to bring high frequency electrical power from the electrosurgical generator to the active electrode and are intended for use with a maximum voltage of 3500 Volts (Vpeak).
Bipolar Forceps and Cords - Olsen Medical® Bipolar Forceps and cords are intended for use in coagulating tissue during surgical procedures and are intended for use with a maximum voltage of 500 Volts (Vpeak)
Monopolar Forces - Olsen Medical® Monopolar Forceps are intended for use in cutting and coagulating tissue during surgical procedures and are intended for use with a maximum voltage of 2000 Volts (Vpeak).
Lietz Electrodes - Olsen Medical® Lietz Electrodes are intended for use in coagulating and cutting tissue during surgical procedures and are intended for use with a maximum voltage of 2000 Volts (Vpeak).
Monopolar Electrodes - Olsen Medical® Electrodes are intended for use in coagulating and cutting tissue during surgical procedures and are intended for use with a maximum voltage of 3500 Volts (Vpeak).
Pens/Pencils and Handles - Olsen Medical® pens/pencils and handles are designed to hold active electrodes and are intended for use with a maximum voltage of 3000 Volts (Vpeak). The distal end of the pencil has a 3/32" (2.36mm) collet that will accept a variety of electrosurgical electrodes with 3/32" (2.36mm) diameter shaft.

**Handling:**
1. Remove and discard the tip protector on all devices before use.
2. Inspect device before each use. Discard or send for repair if damaged.
3. Device use and application to be determined by physician.
4. Cleaning and re-processing steps: reference Symmetry Surgical Instructions for Use (IFU) #204233 for updated information regarding reprocessing the instrument.
• IFU # 204233 is available at: www.Symmetrysurgical.com/IFU

**Storage:**
Reusable devices must be stored individually in a protective container in a clean, dry, moisture free area with suitable protection from extreme temperature and humidity.

**Sterility:**
1. These devices are sold as non-sterile and must be steam sterilized prior to each use. Other methods of sterilization are not recommended and their use will void the warranty.
2. Instruments designed to come apart must be disassembled prior to cleaning.
3. Arrange all devices to allow access of steam to all surfaces. Open hinged devices and ensure devices are disassembled if it is recommended.

| Sterilization Cycle  | Temperature   | Time       | Minimum Dry Time | Device Configuration          |
|----------------------|---------------|------------|------------------|-------------------------------|
| Prevacuum Flash      | 270°F (132°C) | 3 Minutes  | –                | Sterilization Tray            |
| Prevacuum            | 270°F (132°C) | 4 Minutes  | 20 Minutes       | Wrapped/ Individually Pouched |
| Gravity Displacement | 270°F (132°C) | 10 Minutes | 10 Minutes       | Wrapped                       |

**Notices:**

- Electrosurgical equipment and accessories are intended to be used by health professional educated in their use.
- After Use or upon determining product can no longer be used due to wear or damage, safely discard in accordance with established procedures for biohazardous waste.
- Any serious incident that has occurred in relation to the device should be reported to Symmetry Surgical Inc., the governing body, and the Competent authority of the member state in which the incident occurred.
- Please refer to the manufacturer of the generator for warnings, precautions, contraindications, undesirable side-effect, measures to be taken, and limitations of use for the electrosurgical system and accessories
- For applicable IFU/additional information: Visit www.symmetrysurgical.com.

**Caution:**

- Inspect every device before each use. Damaged electrosurgical devices may cause electrical burns to the patient or physician.
- Active electrosurgical devices must be placed away from the patient when not in use. A holster is recommended.
- Use a wet sponge to remove eschar build up from the forceps tips.
- Consult RF generator Owner's Manual. Use only in the active RF generator output connections.
- Electromagnetic compatibility (EMC) – Olsen Medical devices are considered medical electrical equipment. Medical electrical equipment may require precautions regarding EMC and may need to be installed or put into service according to the EMC information provided in the operator’s manual of the generator.
- Do not use electrosurgical equipment unless properly trained to use it in the specific procedure being undertaken. Use of electrosurgical by physicians without such training may result in serious, unintended patient injury.
- Do not perform electrocautery in the presence of flammable anesthetics or other flammable gases, fluids or objects, or in the presence of oxidizing agents, as fire could result.
- Avoid unnecessary or prolonged activation. Activate only when touching or immediately next to the intended point of contact.
- Non-insulated switching contacts of hand controlled forceps must not come in close proximity to or touch the skin or bodily tissue to prevent a burn to the patient or physician.
- Electrode must be compatible with the pen/pencil or handle and ensure that the insulated sleeve fits securely inside the handle. Improperly seated electrode may cause burn or injury.
- Do not bend, cut, or modify any electrosurgical device before or during surgery.
- Keep the voltage/power settings at the lowest possible level to achieve the desired effect.
- Establish a protocol for device examination and maintenance.

**Warning:**
Improper use of any element in electrosurgical systems may cause serious burn. The use of the device with safely functioning and compatible handles, cords, electrosurgical generators and accessories is the user’s responsibility. Read and understand all precautions and instructions before using any electrosurgical device.

Warning: Avoid HF output mode/power settings where maximum output voltage may exceed the stated maximum voltage.

Warning: Do not make any modification to these devices. Modification may cause burn or injury.

Warning: Due to concerns about the carcinogenic and infectious potential of electrosurgical vapors (such as tissue smoke plume and aerosols which may include toxic gasses and vapors, live and dead cellular material, and viruses), protective eyewear, filtration masks, and effective smoke evacuation equipment should be used in both open and laparoscopic procedures. Contact Symmetry Surgical at customerservice@symmetrysurgical.com for additional information or to inquire about our smoke evacuation solutions.

**Warranty:**
All products are warranted to be free from defects in workmanship and materials when used for its intended purpose. Any product that proves defective in workmanship or materials will be repaired or replaced at our discretion free of charge. Repair, modification, or alteration in any product that is performed by anyone not authorized by Olsen Medical® will result in the immediate loss of warranty.

OLSEN® is a registered trademark of Olsen Medical.

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