Symmetry Surgical® QUAD-LOCK® Sterilization Container System

Symmetry Surgical Inc.

3034 Owen Drive
Antioch, TN 37013 USA
☎ 1-800-251-3000
Fax: 1-615-964-5566
www.symmetrysurgical.com

EC REP

Symmetry Surgical GmbH
78532 Tuttlingen, Germany
Maybachstraße 10
☎ +49 7461 96490
Fax: +49 7461 77921

LCN 206296-001/L

© 2010–2017 Symmetry Surgical
♻️ Revised 06/17

QUAD-LOCK® Container System Instructions for Use

IMPORTANT INFORMATION

Please Read Before Use

SYMMETRY SURGICAL QUAD-LOCK® Container System

Indications

The Symmetry Surgical Quad-Lock® Container System is indicated for use by hospitals and by health care facilities to:

- organize and protect stainless steel, aluminum, titanium, plastic* and silicone surgical instruments during sterilization
- allow sterilization of the contained instruments by prevacuum steam, ethylene oxide, and low temperature hydrogen peroxide sterilization cycles (validated modalities and parameters are shown in Table 1)
- maintain sterility of the instruments and accessories after sterilization by prevacuum, ethylene oxide or low temperature hydrogen peroxide for up to 180 days during storage and transport within the healthcare facility as long as the container has not been opened and the integrity of the filters and container have not been compromised.

Note: The system has not been tested for maintenance of sterility after transportation outside the health care facility.

The Quad-Lock system is intended for use with surgical instruments and accessories made with stainless steel, aluminum, titanium, plastic* and silicone. However not all medical device/instrument materials are compatible with each sterilization modality and some restrictions are noted below. Refer to the instrument/sterilizer manufacturer's instructions for sterilization recommendations and limitations.

*The sterilization of plastic instruments was confirmed by Symmetry Surgical for steam sterilization due to the high temperature exposure. Please reference the sterilization equipment manufacturer's information regarding material compatibility for each sterilizer/cycle type.

For effective sterilization and drying of any size Symmetry Surgical Quad-Lock® Sterilization Container, the recommended maximum combined weight of the single container and its contents is defined in Table 1.

For each sterilization cycle type the qualified lumen quantity and size (minimum diameter and maximum length) that has been validated for that cycle is shown in Table 1.

Container system component compatibility by sterilization cycle type is detailed in Table 2.

Product Description

The Symmetry Surgical Quad-Lock® Container System provides a means to organize and protect stainless steel, aluminum, titanium, plastic*, and silicone general surgical instruments and accessories during sterilization and storage. The filter lids, security lids, and container bottoms are made from anodized aluminum (see Figure 1). The security lids help ensure sterility during storage and protect the filters when the containers are stacked. Testing of assembled, loaded containers up to the weights specified in Table 1 has shown that the contents can be effectively sterilized when the validated sterilization parameters are followed. Hospitals should refer to AAMI/AORN guidelines on weights and weight limit. Please refer to the Indications section for specific information regarding the sterilization of lumen instruments.

Precautions

1. Use only the filters available from Symmetry Surgical that are designed for use with the Symmetry Surgical Quad-Lock® Container System. The use of other filters might compromise the sterilization process. The disposable filters, paper indicator cards, and security seals are for single use only.
2. Personnel should wear all appropriate personal protective equipment (PPE) as defined in hospital operating procedures or as recommended by the sterilizer manufacturer.
3. Colored filter lids may be used with steam sterilization ONLY. The grey filter lids can be used with all sterilization modalities. See Table 1 for complete compatibility information.

4. Container stacking can only be done when using steam sterilization. Do not stack more than three containers within the steam sterilizer.
5. Process Indicators and Tamper Evident Seals: Per AAMI ST79, process indicators are used to indicate that the container has been exposed to the sterilization process and to distinguish between processed and unprocessed containers. Use of indicators should be in accordance with the facility's policies and procedures. Steam, Ethylene Oxide, and Hydrogen Peroxide sterilization chemical indicators should be used for each sterilization load.

DO NOT REUSE. These products include the following:

Catalog Number Description

50-9008 Security Seal
50-9011 Filter, Disposable, Round
Indicator cards

CAUTION: Symmetry Surgical® Single Use devices have not been designed to undergo or withstand any form of alteration, such as disassembly, cleaning or re-sterilization, after a single patient use. Reuse can potentially compromise device performance and patient safety.

Inspection

CAUTION: Only appropriately trained personnel must perform the inspection of the sterilization system components.

Inspect the container parts during cleaning and assembly.

Do not use if there are visible signs of excessive wear (e.g. cracking, peeling, or flaking) of the internal anodized layer, damaged, or defective parts; replace these products. Worn or damaged components could compromise the sterile barrier or the use of the sterile contents. The gaskets cannot be repaired

A. Inspecting the Filter Lid

1. Visually check the lid; warpage is not acceptable. Ensure the rim of the lid and the gasket area has not been damaged in a manner that will compromise the ability of the lid to maintain a seal when lid is placed upon container bottom.
2. Visually inspect all gaskets inside the lid. Cuts or cracks in a gasket are not acceptable. The gasket in the perimeter of the lid must be seated in the channel along its entire length.
3. Test the gasket for elasticity by compressing it with a finger. The gasket must return to its original shape without visible indentation.
4. Do not use if there are visible signs of excessive wear (e.g. cracking, peeling, or flaking).

B. Inspecting the Security Lid

Visually check the security lid; warpage is not acceptable.

C. Inspecting the Container Bottom

1. Inspect the rolled upper edge of the container bottom. Sharp edges or other defects that could damage the gasket in the lid when the parts are assembled are not acceptable.
2. Visually inspect all gaskets inside the perforated bottom. Cuts or cracks in a gasket are not acceptable.
3. Do not use if there are visible signs of excessive wear (e.g. cracking, peeling, or flaking).

D. Inspecting the Filter Retainer

These instructions pertain to filter retainers located in both the lid and the perforated container bottom.

1. Visually check the retainer; warpage is not acceptable.
2. Visually inspect the gasket on the bottom side of the retainer (see Figure 2). Cuts or cracks in a gasket are not acceptable.
3. Test the action of the spring latch of the retainer by pressing the latch.
4. Do not use if there are visible signs of excessive wear (e.g. cracking, peeling, or flaking).

Observe the latch for smooth movement (see Figure 3).

Replacement filter retainers are available:

Catalog Number Description

50-9044 Filter Retainer

QUAD-LOCK® Container System Instructions for Use

Cleaning

The containers, lids, and accessories must be cleaned and thoroughly dried before use, including the first use.

CAUTION: Do not use abrasive cleaning agents, metal brushes, or abrasive cleaning pads on the container components. The use of abrasive materials can damage the anodized aluminum components.

A. Disassembly

1. Remove and discard all paper indicator cards.
2. Remove the filter retainers and discard the disposable filters.
3. Remove all contents from the container. Remove the instruments from the baskets.
4. Remove the security lid from the container lid by pressing the release latch inside the security lid, as shown in Figure 4. Lift the container lid out by pulling up on one of the locks.
5. Perform cleaning, as outlined in one of the next two sections.

B. Manual Cleaning

1. Prepare a neutral pH enzymatic cleaning solution in accordance with the manufacturer's instructions.
2. Soak soiled components of the container system, including accessories, for 5 minutes.
3. With the products submerged under water, use a soft bristle brush to remove all traces of blood and debris; pay close attention to the container locks, the security lid latches, and any hard to reach areas, textured surfaces, or crevices. Open and close the container locks and press and release the security lid latches during cleaning to free trapped blood and debris.
4. Rinse the products thoroughly with warm tap water. Pay close attention to the container locks and the security lid latches. Open and close the container locks and press and release the security lid latches while rinsing.
5. Prepare a neutral pH detergent in accordance with the manufacturer's instructions.
6. With products submerged under water, use a soft sponge to thoroughly clean the products.
7. Rinse the products thoroughly with warm tap water. Failure to remove the detergent residue can affect the anodized aluminum finish of the container.
8. Dry all components thoroughly with a soft dry cloth.

C. Mechanical Cleaning

1. Prepare a neutral pH enzymatic cleaning solution in accordance with the manufacturer's instructions.
2. Soak soiled components of the container system, including accessories, for 5 minutes.
3. With the products submerged under water, use a soft bristle brush to remove all traces of blood and debris; pay close attention to the container locks, the security lid latches, and any hard to reach areas, textured surfaces, or crevices. Open and close the container locks and press and release the security lid latches during cleaning to free trapped blood and debris.
4. Rinse the products thoroughly with warm tap water. Pay close attention to the container locks and the security lid latches. Open and close the container locks and press and release the security lid latches while rinsing.
5. Place components of the container system, including accessories in a suitable validated washer/disinfector. Follow the washer/disinfector manufacturer's instructions for loading the components for maximum cleaning exposure.
6. Process the container components and accessories using a standard washer/disinfector instrument cycle according to the manufacturer's instructions. The following minimum wash cycle parameters are recommended:

Cycle	Description
1	Pre-wash --- Cold Softened Tap Water --- 2 minutes
2	Enzyme Spray & Soak --- Hot Softened Tap Water --- 1 minute
3	Rinse --- Cold Softened Tap Water
4	Detergent Wash --- Hot Tap Water (64-66°C/146-150°F) --- 2 min.
5	Rinse --- Hot Purified Water (64-66°C/146-150°F) --- 1 minute
6	Hot Air Dry (116°C/240°F) --- 7 to 30 minutes

Notes:

- The washer/disinfector manufacturer's instructions should be followed.
- A washer/disinfector with demonstrated efficacy (e.g. FDA approval, validated to ISO 15883) should be used.
- Dry time is shown as a range because it is dependent upon the load size placed into the washer/disinfector.
- Many manufacturers pre-program their washer/disinfectors with standard cycles and they may include a thermal low-level disinfection cycle after the detergent wash. The thermal disinfection cycle should be performed to achieve a minimum value $A_0 = 600$ (e.g. 90°C/194°F for 1 minute according to ISO 15883-1) and is compatible with Symmetry instruments.
- If a lubrication cycle is available that applies a water-soluble lubricant such as Preserve®, Instrument Milk or equivalent it is acceptable to use on Symmetry Instruments unless otherwise indicated.

Assembly for Use

A. Assembling the Wire Baskets

Assemble the accessories in the baskets to produce the configurations desired. Use only Symmetry Surgical® QUAD-LOCK® Sterilization Container accessories, screws, and nuts with the wire baskets to assure compatibility. Installing the Instrument Holding Studs and Studs with Slots

- Install either type of stud by positioning the shaft within the basket, and screwing the nut onto the shaft from below the wire basket (Figure 5).
- Both the instrument holding studs and the studs with slots can be used to hold ringed instruments (Figure 6).
- Studs with slots are designed to hold divider plates (Figure 7). In general, use two studs for each divider.

Installing Instrument Support Sliding Bars

- First slide the plastic bar(s) onto the rail of the instrument support (Figure 8).
- Then use the stainless steel screws and nuts to install the instrument support in the wire basket, as desired (Figure 9).

Installing the Spring Holding Plates and Limiting Bars

Use the stainless steel screws and nuts to install the spring holding plates and limiting bars in the wire basket, as desired (Figures 10 and 11). Position instruments in the spring holding plates so they are securely held by the spring.

Installing the Fixation Clamps

Install the fixation clamps by sliding the spring-loaded retainer through the mesh of the basket, with the small tab going through the mesh (Figure 12). Position instruments in the clamps (Figure 13).

Installing the Drape Clamps

Slide the clamps over the drape-covered edges of the basket (Figure 14).

QUAD-LOCK® Container System Instructions for Use

B. Inserting Filters

Insert new filter(s) in each lid and perforated container bottom before each use. The polypropylene filters (catalog no. 50-9011) are for **single use only; DO NOT REUSE.**

1. Remove the filter retainer by squeezing the spring latch in the handle (see Figure 15), and lifting the retainer off.
2. Visually inspect the new filter by holding it up to the light. If there are pinholes or tears, discard the filter.
3. Place a new filter on the spindle.
4. Inspect the gasket on the bottom side of the retainer and test the action of the spring latch as outlined in *Inspection*.
5. Test the action of the spring latch of the retainer by pressing the latch. Observe the latch for smooth movement. If the latch does not operate smoothly, replace the product.
6. Reinstall the retainer by centering it over the spindle (ensure that the printed side of the filter faces up). Squeeze the spring latch and partially seat the retainer so it engages the spindle. Release the spring latch. With the flat of the hand, push the retainer fully onto the spindle until a "click" is heard.
7. Inspect the entire edge of the retainer to ensure that the gasket is firmly compressed against the lid (see Figure 16).

If the gasket is not compressed all around, set that retainer aside and use another retainer. Repeat Steps 4 and 15 through 17.

C. Loading the Instruments

CAUTION: When using a perforated container bottom

1. Position the instruments in the wire baskets in accordance with hospital routines. Assemble the instruments to allow uniform exposure to the sterilization agent. Ensure that the instruments are placed so that water will not pool inside the instrument, but will drip downward.

CAUTION: Do not wrap the containers or the instrument baskets in any type of hospital wrap prior to sterilization. No absorbent sheets or cloths are required in the container bottom for sterilization.

2. When using a perforated container bottom, place a wire base in the container before adding instruments or baskets.
3. Place the instruments and wire basket(s) into the container bottom.

CAUTION: Adhere to the following instructions for usage and placement of the Pin Mat within the containers. Inappropriate use of the pin mat could cause moisture to remain in the container after sterilization.

- Do not use a pin mat in containers that will be stacked in the sterilizer.
- Do not place the pin mat into the container bottom before adding instruments, baskets, or the wire base.

4. Make sure that there is a minimum of 2 inches (51 mm) of free space between the instruments and the top edge of the container bottom.
5. Place the container filter lid on top of the container bottom. Snap the locks down to lock the lid in place. Make sure that the lid securely compresses the gasket against the container bottom. Make sure that the locks are secure by testing that they require sufficient force to open.

6. Insert a security seal into the channel of the container lid lock and secure the seal (see Figure 17). It is recommended to use two seals per container, located on opposite sides.

7. Use the paper indicator cards to record important information, and insert the indicator cards into the appropriate holders on the container bottom. Insert the metal identification plate(s) into the appropriate holders, if desired.

Sterilization

CAUTION: Do not sterilize the container with the security lid on. The security lid will prevent steam from entering the perforations in the container lid and could prevent effective sterilization of the contents.

CAUTION: When sterilizing with Ethylene Oxide or Hydrogen Peroxide do not stack the containers in the sterilizer. The container has not been validated for stacking with these cycle types.

CAUTION: When sterilizing with STERRAD Low Temperature Hydrogen Peroxide cycles do not use the silicone mats. These mats are not compatible with these cycle types and may cause cycle termination:

Protective Silicone Mat		Not Compatible with Sterrad Cycles
508550	Half	
508551	Three-Quarter	
508552	Full	

CAUTION: After sterilization, allow the containers to cool before handling. If containers are handled while hot, use heat proof gloves. If the optional silicone handles are available, hold the containers by the silicone handle covers to prevent burns.

Use a validated, properly maintained and calibrated sterilization equipment. Effective sterilization can be achieved using the cycle parameters and load configurations described in Table 1.

General Precautions

Storage

CAUTION: Do not place hot containers in a cold draft, or on a cold floor. This can cause the container to deform or may cause condensation.

CAUTION: Place the heaviest containers on the bottom of the stack. Over stacking may damage the container.

CAUTION: Do not stack the containers higher than 60cm for storage.

After sterilization cycle and when the containers are completely cooled, an optional security lid may be placed atop the container lid to provide additional protection for the filters. Align the arrow of the security lid with one of the notches on the lid (See Figure 18). Press the security lid down until it clicks. Store the sterilized containers in a clean, dry, secure environment. Under these conditions, the contents of the unopened containers will maintain sterile for 180 days.

Storage of the Disposable Filters and Indicator Cards

Store the disposable filters and the paper indicator cards in a cool, dry place.

Opening the Container

1. Verify the identification on the outside of the container. Verify that the sterilization indicator has changed color, when exposed to each cycle type..
2. Cut the security seals and remove them from the container lid locks.
3. Lift all four locks on the container lid and lift the lid, if using a security lid, lift the lid with the security lid in place away from the container bottom.
4. Remove the contents of the container under aseptic conditions.

Service

For service outside the United States, contact your local Symmetry Surgical representative.

Inside the United States, call the Symmetry Surgical at 1 800 251 3000.

Summary of Testing

Sterilization performance studies were conducted for the Symmetry Surgical Quad-Lock® Container Systems for the cycles described in Table 1 with a total loaded container weight as listed in Table 2 including baskets, instruments, and accessories. All met the acceptance criteria for effective sterilization. Shelf life testing has been conducted to support 180 days for all cycle types. Reuse testing up to 100 cycles showed no effect on the function of the container after processing. Further details on container performance testing can be obtained from Symmetry Surgical. Qualified container contents and lumen size and quantity noted in Table 1 **were qualified by Symmetry Surgical.** Please reference the sterilization equipment manufacturer's information regarding material compatibility for each sterilizer/cycle type.

QUAD-LOCK® Container System Instructions for Use

Warranty

Symmetry Surgical, Inc., warrants that this medical device is free from defects in both materials and workmanship for one (1) year from the date of purchase. Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed. Suitability for use of this medical device for any particular surgical procedure should be determined by the user in conformance with the manufacturer's instructions for use. There are no warranties that extend beyond the description on the face hereof. Abuse or misuse of the product or failure to comply with the instructions for use shall void this warranty.

Quad-Lock® is a registered trademark of Symmetry Surgical

QUAD-LOCK® Container System Instructions for Use

Table 1
Sterilization Parameters

Sterilization Method	Cycle Type	Cycle Parameters	Instrument Type	Qualified Lumen Size/Maximum Quantity	Load Requirements (Maximum Weight)
Steam	Prevacumm/ Dynamic air removal	Temperature: 132°C Exposure Time: 4 minutes Drying Time: 30 minutes Stackable – Yes, up to three containers	Surgical instruments with conjoined surfaces which meet, touch or unite, silicone brackets, lumen devices, and silicone mats	Full and Three Quarter Sizes: 3 mm x 400mm Half Size: 3mm x 250mm size 4 lumen per container	Full, Three Quarter and Half Sizes: 25 lbs
Ethylene Oxide	100% ETO	Concentration: 725mg/L Temperature: 55°C Exposure Time: 60 minutes Humidity: 50-70% Aeration time 8 hours at 110°F (43°C) Stackable – No	Surgical instruments with conjoined surfaces which meet, touch or unite, silicone brackets, lumen devices, and silicone mats	Full and Three Quarter Sizes: 3 mm x 400mm Half Size: 2mm x 250mm 4 Lumen per container	Full, Three Quarter and Half Sizes: 25 lbs
Low Temperature Hydrogen Peroxide Sterilizers by Model Name/Cycle Type Sterilization Parameters for each cycle type are set by the manufacturer					
STERRAD 100NX	Standard	Fixed Cycle Type* Stackable – No	Surgical instruments with conjoined surfaces which meet, touch or unite, silicone brackets, and lumen devices.	Full Size: 2mm x 400mm Three Quarter and Half Sizes: 2mm x 250mm 5 Lumen per container	Full Size: 21.4 lb Three Quarter and Half Sizes: 13.85 lbs
	DUO	Fixed Cycle Type* Stackable – No	Surgical instruments with conjoined surfaces which meet, touch or unite, silicone brackets, and lumen devices.	Full Size: 2mm x 400mm Three Quarter and Half Sizes: 2mm x 250mm 5 Lumen per container	Full Size: 21.4 lbs Three Quarter and Half Sizes: 13.2 lbs
	FLEX	Fixed Cycle Type* Stackable - No	Porous Lumens	Full, Three Quarter and Half Sizes:	Full Size: 16.65 lbs

QUAD-LOCK® Container System Instructions for Use

				One 1mm x 850mm Flexible Endoscope per container	Three Quarter Size: 12.40 lbs Half Size: 10.75 lbs
STERRAD NX	Standard	Fixed Cycle Type* Stackable - No	Surgical instruments with conjoined surfaces which meet, touch or unite, silicone brackets, and lumen devices.	Full Size: 1mm x 400mm Three Quarter and Half Size: 2mm x 250mm 5 Lumen per container	Full Size: 21.4 lbs Three Quarter and Half Sizes: 13.85 lbs
	Advanced	Fixed Cycle Type* Stackable - No	Surgical instruments with conjoined surfaces which meet, touch or unite, silicone brackets, and lumen devices.	Full, Three Quarter, and Half Sizes: 2mm x 250mm 5 Lumen per container	Full Size: 21.4 lbs Three Quarter and Half Sizes: 13.85 lbs
STERRAD 100S	Standard	Fixed Cycle Type* Stackable - No	Surgical instruments with conjoined surfaces which meet, touch or unite, silicone brackets, and lumen devices.	Full Size: 3mm x 400mm Three Quarter and Half Size: 2mm x 250mm 5 Lumen per container	Full, Three Quarter, and Half Sizes: 18.2 lbs
STERIS V-PRO™ 1	Lumen	Cycle time = 55 minutes Stackable - No	Surgical instruments with conjoined surfaces which meet, touch or unite, silicone brackets, lumen devices, and silicone mats	Full Size: 3mm x 400mm Three Quarter and Half Sizes: 2mm x 250mm 10 lumen per container	Full, Three Quarter, and Half Sizes: 19.65 lbs
STERIS V-PRO™ 1 PLUS	Lumen or Non-Lumen	Cycle time = 55 minutes Stackable - No	Surgical instruments with conjoined surfaces which meet, touch or unite, silicone brackets, lumen devices, and silicone mats	Full Size: 3mm x 400mm Three Quarter and Half Sizes: 2mm x 250mm 10 lumen per container	Full, Three Quarter, and Half Sizes: 19.65 lbs

QUAD-LOCK® Container System Instructions for Use

STERIS V-PRO™ maX	Lumen or Non- Lumen	Cycle time = 55 minutes Stackable - No	Surgical instruments with conjoined surfaces which meet, touch or unite, silicone brackets, lumen devices, and silicone mats	Full Size: 3mm x 400mm Three Quarter and Half Sizes: 2mm x 250mm 10 lumen per container	Full, Three Quarter, and Half Sizes: 19.65 lbs
----------------------	------------------------	--	---	---	--

*For all STERRAD manufactured sterilizers all cycle parameters are preset by the manufacturer and cannot be reset by the user.

Table 2 Container System Compatibility by Cycle Type					
Product Code		Pre-Vacuum Steam	100% Ethylene Oxide	STERRAD 100S NX Standard NX Advanced 100NX Standard 100NX FLEX 100NX DUO	STERIS V-Pro 1 Plus Lumen
Container Component Compatibility					
Filter Bottom Container Base					
Full Length	Height				
50-8726	100mm (4 inch)	YES	YES	YES	YES
50-8727	135mm (5 inch)	YES	YES	YES	YES
50-8728	150mm (6 inch)	YES	YES	YES	YES
50-8729	200mm (8 inch)	YES	YES	YES	YES
50-8730	260mm (10 inch)	YES	YES	YES	YES
Three-Quarter Length					
50-8736	100mm (4 inch)	YES	YES	YES	YES
50-8737	135mm (5 inch)	YES	YES	YES	YES
50-8738	150mm (6 inch)	YES	YES	YES	YES
Half Length					
50-8731	100mm (4 inch)	YES	YES	YES	YES
50-8732	135mm (5 inch)	YES	YES	YES	YES
50-8733	150mm (6 inch)	YES	YES	YES	YES
50-8734	200mm (8 inch)	YES	YES	YES	YES
50-8735	260mm (10 inch)	YES	YES	YES	YES
Solid Bottom Container Base					
Full Length					
50-8739	100mm (4 inch)	YES	NO	NO	NO
50-8740	135mm (5 inch)	YES	NO	NO	NO
50-8741	150mm (6 inch)	YES	NO	NO	NO
50-8742	200mm (8 inch)	YES	NO	NO	NO
50-8743	260mm (10 inch)	YES	NO	NO	NO
Three-Quarter Length					
50-8749	100mm (4 inch)	YES	NO	NO	NO
50-8750	135mm (5 inch)	YES	NO	NO	NO
50-8751	150mm (6 inch)	YES	NO	NO	NO
Half Length					
50-8744	100mm (4 inch)	YES	NO	NO	NO
50-8745	135mm (5 inch)	YES	NO	NO	NO
50-8746	150mm (6 inch)	YES	NO	NO	NO

QUAD-LOCK® Container System Instructions for Use

Table 2 Container System Compatibility by Cycle Type					
Product Code		Pre-Vacuum Steam	100% Ethylene Oxide	STERRAD 100S NX Standard NX Advanced 100NX Standard 100NX FLEX 100NX DUO	STERIS V-Pro 1 Plus Lumen
Container Component Compatibility					
50-8747	200 (8 inch)	YES	NO	NO	NO
50-8748	260 (10 inch)	YES	NO	NO	NO
Filter Lids					
Product Code	Color				
Full Length					
50-8944	Grey	YES	YES	YES	YES
50-8945	Yellow	YES	NO	NO	NO
50-8946	Green	YES	NO	NO	NO
50-8947	Blue	YES	NO	NO	NO
50-8948	Red	YES	NO	NO	NO
50-8949	Black	YES	NO	NO	NO
Three-Quarter Length					
50-8956	Grey	YES	YES	YES	YES
50-8957	Yellow	YES	NO	NO	NO
50-8958	Green	YES	NO	NO	NO
50-8959	Blue	YES	NO	NO	NO
50-8960	Red	YES	NO	NO	NO
50-8961	Black	YES	NO	NO	NO
Half Length					
50-8950	Grey	YES	YES	YES	YES
50-8951	Yellow	YES	NO	NO	NO
50-8952	Green	YES	NO	NO	NO
50-8953	Blue	YES	NO	NO	NO
50-8954	Red	YES	NO	NO	NO
50-8955	Black	YES	NO	NO	NO
Security Lids (used only after container sterilization is complete)					
Product Code	Color				
Full Length					
50-8869	Grey	YES	YES	YES	YES
50-8870	Yellow	YES	YES	YES	YES
50-8871	Green	YES	YES	YES	YES
50-8872	Blue	YES	YES	YES	YES
50-8873	Red	YES	YES	YES	YES
50-8874	Black	YES	YES	YES	YES
Three-Quarter Length					
50-8881	Grey	YES	YES	YES	YES
50-8882	Yellow	YES	YES	YES	YES
50-8883	Green	YES	YES	YES	YES
50-8884	Blue	YES	YES	YES	YES
50-8885	Red	YES	YES	YES	YES
50-8886	Black	YES	YES	YES	YES
Half Length					
50-8875	Grey	YES	YES	YES	YES
50-8876	Yellow	YES	YES	YES	YES
50-8877	Green	YES	YES	YES	YES
50-8878	Blue	YES	YES	YES	YES
50-8879	Red	YES	YES	YES	YES
50-880	Black	YES	YES	YES	YES

QUAD-LOCK® Container System Instructions for Use

Table 2 Container System Compatibility by Cycle Type					
Product Code		Pre-Vacuum Steam	100% Ethylene Oxide	STERRAD 100S NX Standard NX Advanced 100NX Standard 100NX FLEX 100NX DUO	STERIS V-Pro 1 Plus Lumen
Container Component Compatibility					
Basket Configurations					
Product Code	Height				
Full Length					
50-8500	30mm (1 inch)	YES	YES	YES	YES
50-8501	50mm (2 inch)	YES	YES	YES	YES
50-8502	70mm (3 inch)	YES	YES	YES	YES
50-8503	100mm (4 inch)	YES	YES	YES	YES
Three-Quarter Length					
50-8508	30mm (1 inch)	YES	YES	YES	YES
50-8509	50mm (2 inch)	YES	YES	YES	YES
50-8510	70mm (3 inch)	YES	YES	YES	YES
50-8511	100mm (4 inch)	YES	YES	YES	YES
Half Length					
50-8504	30mm (1 inch)	YES	YES	YES	YES
50-8505	50mm (2 inch)	YES	YES	YES	YES
50-8506	70mm (3 inch)	YES	YES	YES	YES
50-8507	100mm (4 inch)	YES	YES	YES	YES
Wire Base Configurations					
Product Code	Length	Height			
50-8512	Half	26mm (1 inch)	YES	YES	YES
50-8513	Three-Quarter	26mm (1 inch)	YES	YES	YES
50-8514	Full	26mm (1 inch)	YES	YES	YES
Protective Mat					
Product Code	Length				
50-8550	Half	YES	YES	NO	YES
50-8551	Three-Quarter	YES	YES	NO	YES
50-8552	Full	YES	YES	NO	YES

TABLE 3 CASE CONFIGURATIONS				
Bottom (Length X Width)	Filter Lid (Length X Width)	Security Lid (Length X Width)	Content Accessories See Table 1 for instruments and lumen specifications	Required Accessories
Small Containers (Half), Perforated Bottoms				
Small, perforated, with filter retainer 270 mm X 270 mm Height varies as shown:	Small, 4 Latch, perforated, with filter retainer 280 mm X 276 mm Grey REF 508950	Small 292 mm X 290 mm Any one of the colors below will be compatible	Wire Basket, Small REF 50-8504 – 50-8507 Pin Mat, Small REF 50-8551* Any other instrument holding	2 Filters REF 50-9010 Security seal (2 recommended) REF 50-9008

QUAD-LOCK® Container System Instructions for Use

Height: 100 mm -- REF 50-8731 Height: 135 mm -- REF 50-8732 Height: 150 mm -- REF 50-8733 Height: 200 mm -- REF 50-8734 Height: 260 mm -- REF 50-8735	Yellow REF 50-8951* Green REF 50-8952* Blue REF 50-8953* Red REF 50-8954* Black REF 50-8955*	Grey REF 50-8875 Yellow REF 50-8876 Green REF 50-8877 Blue REF 50-8878 Red REF 50-8879 Black REF 50-8880	accessory except Divider Plate 460 mm X 20 mm (REF 50-9029)	2 protective sleeve handles REF 50-9012 – 50-9017 Wire Base, Small REF 50-8512
Small Containers (Half), Non-perforated Bottoms				
Small, non-perforated 270 mm X 270 mm Height varies as shown: Height: 100 mm -- REF 50-8744 Height: 135 mm -- REF 50-8745 Height: 150 mm -- REF 50-8746 Height: 200 mm -- REF 50-8747 Height: 260 mm -- REF 50-8748	Small, 4 Latch, perforated, with filter retainer 280 mm X 276 mm Grey REF 508950 Yellow REF 50-8951* Green REF 50-8952* Blue REF 50-8953* Red REF 50-8954* Black REF 50-8955*	Small 292 mm X 290 mm Any one of the colors below will be compatible Grey REF 50-8875 Yellow REF 50-8876 Green REF 50-8877 Blue REF 50-8878 Red REF 50-8879 Black REF 50-8880	Wire Basket, Small REF 50-8504 – 50-8507 Pin Mat, Small REF 50-8551** Any other instrument holding accessory except Divider Plate 460 mm X 20 mm (REF 50-9029)	1 Filter REF 50-9010 Security seal (2 recommended) REF 50-9008 2 protective sleeve handles REF 50-9012 – 50-9017
Medium (Three quarter) Containers, Perforated Bottoms				
Medium, perforated, with filter retainers 440 mm X 270 mm Height varies as shown: Height: 100 mm -- REF 50-8736 Height: 135 mm -- REF 50-8737 Height: 150 mm -- REF 50-8738	Medium, 4 Latch, perforated, with filter retainer 445 mm X 279 mm Grey REF 508956 Yellow REF 50-8957* Green REF 50-8958* Blue REF 50-8959* Red REF 50-8960* Black REF 50-8961*	Medium 456 mm X 295 mm Any one of the colors below will be compatible Grey REF 50-8881 Yellow REF 50-8882 Green REF 50-8883 Blue REF 50-8884 Red REF 50-8885 Black REF 50-8886	Wire Basket, Medium REF 50-8508 – 50-8511 Pin Mat, Medium REF 50-8552** Any other instrument holding accessory except Divider Plate 460 mm X 20 mm (REF 50-9029)	2 Filters REF 50-9010 Security seal (2 recommended) REF 50-9008 2 protective sleeve handles REF 50-9012 – 50-9017 Wire Base, Medium REF 50-8513
Medium Containers (Three Quarter), Non-perforated Bottoms				
Medium, non-perforated 440 mm X 270 mm Height varies as shown: Height: 100 mm -- REF 50-8749 Height: 135 mm -- REF 50-8750 Height: 150 mm -- REF 50-8751	Medium, 4 Latch, perforated, with filter retainer 445 mm X 279 mm Grey REF 508956 Yellow REF 50-8957* Green REF 50-8958* Blue REF 50-8959* Red REF 50-8960* Black REF 50-8961*	Medium 456 mm X 295 mm Any one of the colors below will be compatible Grey REF 50-8881 Yellow REF 50-8882 Green REF 50-8883 Blue REF 50-8884 Red REF 50-8885 Black REF 50-8886	Wire Basket, Medium REF 50-8508 – 50-8511 Pin Mat, Medium REF 50-8552** Any other instrument holding accessory except Divider Plate 460 mm X 20 mm (REF 50-9029)	1 Filter REF 50-9010 Security seal (2 recommended) REF 50-9008 2 protective sleeve handles REF 50-9012 – 50-9017
Large (Full) Containers, Perforated Bottoms				
Large, perforated with filter retainers 560 mm X 275 mm Height varies as shown: Height: 100 mm -- REF 50-8726 Height: 135 mm -- REF 50-8727	Large, 4 Latch, perforated, with filter retainer 570 mm X 280 mm Grey REF 508944 Yellow REF 50-8945* Green REF 50-8946* Blue REF 50-8947* Red REF 50-8948*	Large 582 mm x 288 mm Any one of the colors below will be compatible Grey REF 50-8869 Yellow REF 50-8870 Green REF 50-8871 Blue REF 50-8872	Wire Basket, Large REF 50-8500 – 50-8503 Pin Mat, Large REF 50-8550** Any other instrument holding accessory	4 Filters REF 50-9010 Security seal (2 recommended) REF 50-9008 2 protective sleeve handles REF 50-9012 – 50-9017 Wire Base, Large

QUAD-LOCK® Container System Instructions for Use

Height: 150 mm -- REF 50-8728 Height: 200 mm -- REF 50-8729 Height: 260 mm -- REF 50-8730	Black REF 50-8949*	Red REF 50-8873 Black REF 50-8874		REF 50-8514
Large (Full) Containers, Non-perforated Bottoms				
Large, non-perforated 560 mm X 275 mm Height varies as shown: Height: 100 mm -- REF 50-8739 Height: 135 mm -- REF 50-8740 Height: 150 mm -- REF 50-8741 Height: 200 mm -- REF 50-8742 Height: 260 mm -- REF 50-8743	Large, 4 Latch, perforated, with filter retainer 570 mm X 280 mm Grey REF 508944 Yellow REF 50-8945* Green REF 50-8946* Blue REF 50-8947* Red REF 50-8948* Black REF 50-8949*	Large 582 mm x 288 mm Any one of the colors below will be compatible Grey REF 50-8869 Yellow REF 50-8870 Green REF 50-8871 Blue REF 50-8872 Red REF 50-8873 Black REF 50-8874	Wire Basket, Large REF 50-8500 – 50-8503 Pin Mat, Large REF 50-8550** Any other instrument holding accessory	2 Filters REF 50-9010 Security seal (2 recommended) REF 50-9008 2 protective sleeve handles REF 50-9012 – 50-9017

*Compatibility of colored lids has only been confirmed for steam sterilization. Do **NOT** use in any other sterilization cycle types

Mats are **NOT compatible with STERRAD sterilization cycles


QUAD-LOCK® Container System Instructions for Use

1

- A. Security Lid
- B. Security Lid latch
- C. Filter lid
- D. Container bottom
- E. holders for small indicator cards or metal tags
- F. Protective sleeve, handle
- G. Holders for large indicator cards or metal tags
- H. Container locks

2


Underside of the filter retainer

- A. Outer gasket
- B. Inner gasket

QUAD-LOCK® Container System Instructions for Use

3

- A. Filter retainer
- B. Retainer handle
- C. Spring latch


4

- A. Security lid
- B. security lid latch
- C. Container lock
- D. Underside of container lid


QUAD-LOCK® Container System Instructions for Use

5

- A. Shaft of the instrument holding stud (seen through the wire basket)
- B. Nut of instrument holding stud
- C. Bottom of wire basket


6

- A. Studs with slots


QUAD-LOCK® Container System Instructions for Use

7

- A. Instrument studs with slots
- B. Divider plate
- C. Wire basket


8

- A. Plastic divider
- B. Instrument support sliding bar


QUAD-LOCK® Container System Instructions for Use

9

- A. Screw
- B. Instrument support sliding bar
- C. Plastic divider
- D. Wire basket

10

- A. Spring Instrument Straight Holding Plate
- B. Spring Instrument Angled Holding Plate


QUAD-LOCK® Container System Instructions for Use

11

A. Limiting bar


12

A. Instrument clamp
B. Tab


QUAD-LOCK® Container System Instructions for Use

13

A. Clamps (3 shown)

14

A. Drape clamps
B. Wire basket (beneath the drape)


QUAD-LOCK® Container System Instructions for Use

15

- A. Filter retainer
- B. Retainer handle
- C. Spring latch
- D. Disposable filter
(seen through the retainer)
- E. Inside of the container lid


16

- A. Gasket is firmly compressed against the lid

QUAD-LOCK® Container System Instructions for Use

17

- A. Security Seal
- B. Container lock
- C. Channel in the lock
- D. Label indicator cards (small)
- E. Protective sleeve on handle
- F. Label indicator card (large)

18

- A. Security lid
- B. Arrow outside the security lid latch
- C. filter lid
- D. Notch in the filter lid

QUAD-LOCK® Container System Instructions for Use

-  US Representative
-  Quantity
-  Manufacturer
-  Do not reuse
-  Made in

DANISH SYMMETRY SURGICAL QUAD-LOCK®-Beholdersystem

VIGTIG INFORMATION

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SYMMETRY SURGICAL QUAD-LOCK®-Beholdersystem

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