EN **AUTOCALVABLE 3-BUTTON POWERCONTROL HANDPIECE**

REFA902 – For use with A952 ESU and DERM 101/102/941/942 desiccators

- Store At Room Temperature
- Collet .092” Min - .094” Max
- Device is rated for 6300 Vpeak

WARNING:

Do not use in the presence of flammable materials. Fire could result.

Improper electrode installation may result in injury to the patient or operating room personnel by arcing at the electrode pencil connection.

To avoid burn or shock injury, turn off the desiccator at the mains power.

To avoid burn or shock injury, do not activate pencil while wiping electrode tip during use.

Keep active accessories away from the patient when not in use. Accessories and electrodes should be placed in a clean, dry, nonconductive and highly visible area when not in use.

Inspect prior to use. The desiccator should be replaced if any physical damage such as nicks, cracks or cuts is noted on the hardware or cable.

Due to concerns about the carcinogenic and infectious potential of electrospun byproducts (such as tissue smoke plume and aerosols (which may include toxic gasses and vapors, live and dead cellular material, and viruses), protective eyewear, filtration masks, and effective smoke evacuation equipment should be used in both open and laparoscopic procedures. Contact Aspen Surgical at customer-service@aspensurgical.com for additional information or to inquire about our smoke evacuation solutions. Reusable devices are provided non-sterile and must be processed prior to use to include first use. Processing is defined on specific instructions provided with each device.

CAUTION:

This product is NON-STERILE. The use of a sterile sheath or sterilizing prior to use is recommended. The device is designed for use only with the A952 ESU and DERM 101/102/941/942 desiccators.

Use the lowest possible setting on the A952 ESU and DERM 101/102/941/942 desiccators that will deliver an optimal surgical effect.

NOTICE:

Electroscopic equipment and accessories are intended to be used by health professionals educated in their use.

Please refer to the manufacturer of the generator for warnings, precautions, contra-indications, undesirable side-effects, measures to be taken, and limitations of use for the electrospurgical system and accessories.

For applicable IFU/additional information: Visit <https://www.aspensurgical.com/Resources/Documents/IFUs>. After Use or upon determining product can no longer be used due to wear or damage, safely discard in accordance with established procedures for biohazardous waste.

Any serious incident that has occurred in relation to the device should be reported to Aspen Surgical Products Inc., the governing body, and the Competent authority of the member state in which the incident occurred.

DIRECTIONS FOR USE:

- Place electrode into the handpiece. Ensure the electrode is fully inserted.
- Plug the handpiece into the monopolar output jacks on the lower left of the front of the unit. The plug is designed to fit in only one direction.
- To remotely adjust the power settings on the unit press the up or down button on the handpiece. De-pressing and holding the up or down buttons will cause the power settings to change more rapidly.
- To activate the monopolar output, depress the “ON” button on the handpiece. This is the button located nearest the electrode.

A5. After use, remove electrode and sterilize the pencil following the recommended procedure as follows.

CLEANING INSTRUCTIONS:

- Remove all gross matter (blood, tissue) by wiping down the pencil using a gauze pad or soft cloth dampened in warm tap water until all visible debris is removed.
- Wipe the pencil using a gauze pad or soft cloth dampened in an enzymatic detergent mixture diluted according to the manufacturer’s instructions. The pencil may also be gently scrubbed with a soft cleaning brush.
- IMPORTANT: Abrasive material should not be used in the area of the barcode when cleaning the device. A cleaning removing agents from pencil by wiping with tap water dampened cloth.

STERILIZATION INSTRUCTIONS:

Number of uses when sterilized by steam autoclave: 25
1. Steam autoclave at 250° F (121° C, gravity) for 20 minutes using an autoclavable pouch or wrapped cloth method. Consult the AAMI recommended practice for steam sterilization using the wrapped cloth method. Allow product to cool to 120° C for 30 minutes before use.

2. Flash autoclave method is not recommended.

Intended Use:

The Device is an accessory to be compatible with any standard RF generator that accepts a standard 3-pin monopolar connector and serves as the pathway to deliver RF energy via an assortment of RF electrodes to coagulate different kinds of tissues.

Indications for Use:
Electrospurgical Pencils and Pens are used in conjunction with electrospurgical electrodes and generator to deliver RF Energy used to cut and excise tissue or to coagulate blood vessels/soft tissues.

Patient Population:

Infant to geriatric.

Contraindications:

None.

Clinical Benefits:

The primary patient benefits are to facilitate a surgical procedure through providing a means for cutting, coagulation, and ablation of tissue.

Site of Use and Site Conditions:

- Site of use – Soft tissue
- Site conditions – Aseptic

Intended User:

–Education – Trained physician, physician’s assistant, nurse, nurse practitioner. No maximum

- Knowledge: no maximum
- Minimum:
 - Understands electrospurgery and electrospurgical techniques;
 - Read and understands supplied IFU (Accompanying Document)
 - Understands hygiene
- Maximum:
 - There is no maximum

- Experience:
 - Minimum:
 - Some training on techniques or training under surveillance/supervision
 - No special experience needed
 - Maximum:
 - There is no maximum
- Permissible Impairments:
 - Mild reading / vision impairment or vision correction to 20/20
 - Hearing impaired by 40% resulting in 60% of normal hearing at 500 Hz to 2.0 kHz.

Intended Conditions for Hospital:

- Doctor office, Surgery center or Hospital, intended for professional use only
- Li. clear and unobstructed view at point of use
- Ambient Temperature 10° C to 40° C

Device Life:

The A902 Pencil is reusable and designed to withstand a maximum of 25 gravity autoclave sterilization cycles

25 uses

Disposal:
Biohazard – Dispose of Pencil according to facility procedures

DA **AUTOKLAVERBAR 3-KNAPSAFBRYDERKONTROLLHANDSTYCKE**
REFA902 – TTI brug med A952 ESU og Derm 101/102/941/942 desiccators

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1. Steam autoclave at 250° F (121° C, gravity) for 20 minutes using an autoclavable pouch or wrapped cloth method. Consult the AAMI recommended practice for steam sterilization using the wrapped cloth method. Allow product to cool to 120° C for 30 minutes before use.

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25 uses
1. Steam autoclave at 250° F (121° C, gravity) for 20 minutes using an autoclavable pouch or wrapped cloth method. Consult the AAMI recommended practice for steam sterilization using the wrapped cloth method. Allow product to cool to 120° C for 30 minutes before use.

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Knowledge: no maximum

- Minimum:
 - Understands electrospurgery and electrospurgical techniques;
 - Read and understands supplied IFU (Accompanying Document)
 - Understands hygiene
- Maximum:
 - There is no maximum

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