INSTRUCTIONS FOR USE: InsufflatOR Tubing Set™ Insufflation Tubing Set with 0.1µ Filter

EN

INDICATIONS FOR USE

The InsufflatOR Tubing Set™ is intended to filter and deliver insufflation gas (CO₂) from an insufflation pump to an insufflation needle or suitable laparoscopic port as a means to provide peritoneal insufflation and allow for the manipulation of instruments during laparoscopic procedures.

CONTRAINDICATIONS

The device is not intended for use except as indicated, and when laparoscopic procedures are contraindicated.


CAUTIONS / WARNINGS

- Verify compatibility of the InsufflatOR Tubing Set™ and connectors prior to initiation of the procedure.
- The InsufflatOR Tubing Set™ was designed, tested and manufactured for single patient use only. Reuse or reprocessing of this device may lead to its failure and subsequent patient injury.
- Do not re-sterilize the device. Re-sterilization may compromise the integrity of the device and create the risk of contamination and unintended patient infection.

INSTRUCTIONS FOR USE

Read and become familiar with the instructions and cautions before using the device. You should have a thorough understanding of the principles and techniques of insufflation and surgical laparoscopy before using the product.

- Remove the InsufflatOR Tubing Set™ from its sterile packaging. Inspect the device for damage caused by shipping. Discard if any damage is noted. Otherwise place on sterile field.
- Pass the filter and the short section of the tubing from the sterile field and attach to the insufflator using the appropriate connector (Figure 1).
- Connect long section of tubing to the Veress / Insufflation needle (Figure 1).

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Caution, consult accompanying documents.

Consult Instructions for Use.

Keep away from sunlight.

Keep dry.

Not made with natural rubber latex.

Method of sterilization using ethylene oxide.

Single sterile barrier system

Use by date

Caution: US Federal Law restricts this device to sale by or on the order of a physician.

Do not use if package is damaged

Store between 15°C to 25°C temperature.

Do not re-sterilize.

Product conforms with essential requirements in European Medical Devices Directive 93/42/EEC. Notified body: DNV Product Assurance AS, with CE2460

Single sterile barrier system with protective packaging outside

Medical Device.

Australian Sponsor: The O R Company Pty Ltd, 4/47 Wangara Road, Cheltenham, VIC, 3192, AUSTRALIA.

Distributed exclusively in USA by: Symmetry Surgical, 3034 Owen Drive, Antioch, TN 37013, USA. Email: customerservice@symmetrysurgical.com Telephone: USA & Canada +1 800-251-3000 (Monday - Friday, 7:00AM - 5:30PM CST)

Advanced Medical Design Co. Ltd, 4F-5F, No. 29, Wuquan 5th Road, Wugu District, 24888 New Taipei City, TAIWAN www.amdintl.com.tw

OBELIS S.A, Bd. Général Wahis, 53, 1030 Brussels, Belgium Telephone: +3227325954, Fax: +3227326003, E-mail: mail@obelis.net