## INSTRUCTIONS FOR USE: InsufflatOR Tubing Set™ Insufflation Tubing Set with 0.1µ Filter

EN

### INDICATIONS FOR USE

The InsufflatOR Tubing Set™ is intended to filter and deliver insufflation gas (CO<sub>2</sub>) from an insufflation pump to an insufflation needle or suitable laparoscopic port as a means to provide peritoneal insufflation and allow for the manipulation of instruments during laparoscopic procedures.

### CONTRAINDICATIONS

The device is not intended for use except as indicated, and when laparoscopic procedures are contraindicated.


### CAUTIONS / WARNINGS

- Verify compatibility of the InsufflatOR Tubing Set™ and connectors prior to initiation of the procedure.
- The InsufflatOR Tubing Set™ was designed, tested and manufactured for single patient use only. Reuse or reprocessing of this device may lead to its failure and subsequent patient injury.
- Do not re-sterilize the device. Re-sterilization may compromise the integrity of the device and create the risk of contamination and unintended patient infection.

### INSTRUCTIONS FOR USE

Read and become familiar with the instructions and cautions before using the device. You should have a thorough understanding of the principles and techniques of insufflation and surgical laparoscopy before using the product.

1. Remove the InsufflatOR Tubing Set™ from its sterile packaging. Inspect the device for damage caused by shipping. Discard if any damage is noted. Otherwise place on sterile field.
2. Pass the filter and the short section of the tubing from the sterile field and attach to the insufflator using the appropriate connector (Figure 1).
3. Connect long section of tubing to the Veress / Insufflation needle (Figure 1).

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