For product information, contact:
Yderligere produktoplysninger fås ved at kontakte:
Pour obtenir des informations sur le produit, veuillez contacter :
Für Produktinformation kontaktieren Sie:
Για πληροφορίες σχετικά με το προϊόν, επικοινωνήστε με:
Per informazioni sul prodotto, contattate:
製品情報、お問い合わせ先:
Para obter informacões sobre o produto, contacte:
Para obtener información del producto, póngase en contacto con:

symmetry® olsen®

3230 Commerce Center Place
Louisville, Kentucky, 40211, USA
(502) 813-5022 <> (800) 297-6344
See our products at
Se vores produkter på
Découvrez nos produits sur
Siehe unsere Produkte auf
Δείτε τα προϊόντα μας στην ιστοσελίδα
È possibile vedere i nostri prodotti all'indirizzo web
当社の製品は以下でご覧いただけます
Consulte os nossos produtos em
Consulte nuestros productos en
<http://www.symmetryurgical.com>

Symmetry Surgical GmbH
Maybachstraße 10, 78532
Tuttingen, Germany
+49 7461 96490

R ONLY

Reusable Electrosurgical Devices
Genanvendelige elektrokirurgiske Instrumenter
Appareils d'électrochirurgie réutilisables
Wiederverwendbare elektrochirurgische Geräte
Επαναχρησιμοποιήσιμες Ηλεκτροχειρουργικές Συσκευές
Dispositivi elettrochirurgici riutilizzabili
再使用可能電気外科装置
Dispositivos electrocirúrgicos reutilizáveis
Dispositivos electroquirúrgicos reutilizables

ENGLISH

Description:
Reusable devices are provided non-sterile and must be processed prior to use to include first use. Processing is defined on specific instructions provided with each device. Each reusable device is designed to withstand up to 20 sterilization cycles when properly cared for and sterilized according to the instructions below. These devices are designed to function as active electrodes or as accessories to active electrodes. This family of devices includes, but is not limited to: monopolar and bipolar forceps, pens/pencils and electrode handles, electrodes, cords and adapters.

Intended Use (per EU definition) / Indications (per FDA definition):
Monopolar Cord and Adapters - Olsen Medical® electrosurgical cords and adapters are designed to bring high frequency electrical power from the electrosurgical generator to the active electrode and are intended for use with a maximum voltage of 3500 Volts.
Bipolar Forceps and Cords - Olsen Medical® Bipolar Forceps and cords are intended for use in coagulating tissue during surgical procedures and are intended for use with a maximum voltage of 500 Volts.
Monopolar Forceps - Olsen Medical® Monopolar Forceps are intended for use in cutting and coagulating tissue during surgical procedures and are intended for use with a maximum voltage of 2000 Volts.
Lietz Electrodes - Olsen Medical® Lietz Electrodes are intended for use in coagulating and cutting tissue during surgical procedures and are intended for use with a maximum voltage of 2000 Volts.
Monopolar Electrodes - Olsen Medical® Electrodes are intended for use in coagulating and cutting tissue during surgical procedures and are intended for use with a maximum voltage of 3500 Volts.
Pens/Pencils and Handles - Olsen Medical® pens/pencils and handles are designed to hold active electrodes and are intended for use with a maximum voltage of 3000 Volts.
The distal end of the pencil has a 3/32" (2.36mm) collet that will accept a variety of electrosurgical electrodes with 3/32" (2.36mm) diameter shaft.

Handling:

- Remove and discard the tip protector on all devices before use.
- Inspect device before each use. Discard or send for repair if damaged.
- Device use and application to be determined by physician.
- Cleaning and reprocessing steps: reference Symmetry Surgical Instructions for Use (IFU) #204233 for updated information regarding reprocessing the instrument.

- IFU # 204233 is available at: www.Symmetryurgical.com/IFU

Storage:
Reusable devices must be stored individually in a protective container in a clean, dry, moisture free area with suitable protection from extreme temperature and humidity.

Sterility:

- These devices are sold as non-sterile and must be steam sterilized prior to each use. Other methods of sterilization are not recommended and their use will void the warranty.
- Instruments designed to come apart must be disassembled prior to cleaning.
- Arrange all devices to allow access of steam to all surfaces. Open hinged devices and ensure devices are disassembled if it is recommended.

Steam Sterilization:

Sterilization Cycle	Temperature	Time	Minimum Dry Time	Device Configuration
Prevacuum Flash	270°F (132°C)	3 Minutes	–	Sterilization Tray
Prevacuum	270°F (132°C)	4 Minutes	20 Minutes	Wrapped/ Individually Pouched
Gravity Displacement	270°F (132°C)	10 Minutes	10 Minutes	Wrapped

Notices:

- Electrosurgical equipment and accessories are intended to be used by health professional educated in their use.
- After Use or upon determining product can no longer be used due to wear or damage, safely discard in accordance with established procedures for biohazardous waste.
- Any serious incident that has occurred in relation to the device should be reported to Symmetry Surgical Inc., the governing body, and the Competent authority of the member state in which the incident occurred.
- Please refer to the manufacturer of the generator for warnings, precautions, contraindications, undesirable side-effect, measures to be taken, and limitations of use for the electrosurgical system and accessories
- For applicable IFU/additional information: Visit www.symmetryurgical.com.

Caution:

- Inspect every device before each use. Damaged electrosurgical devices may cause electrical burns to the patient or physician.
- Active electrosurgical devices must be placed away from the patient when not in use. A holster is recommended.
- Use a wet sponge to remove eschar build up from the forceps tips.
- Consult RF generator Owner's Manual. Use only in the active RF generator output connections.
- Do not perform electrosurgery in the presence of flammable anesthetics or other flammable gases, fluids or objects, or in the presence of oxidizing agents, as fire could result.
- Avoid unnecessary or prolonged activation. Activate only when touching or immediately next to the intended point of contact.
- Non-insulated switching contacts of hand controlled forceps must not come in close proximity to or touch the skin or bodily tissue to prevent a burn to the patient or physician.
- Electrode must be compatible with the pen/pencil or handle and ensure that the insulated sleeve fits securely inside the handle. Improperly seated electrode may cause burn or injury.
- Do not bend, cut, or modify any electrosurgical device before or during surgery.
- Keep the voltage/power settings at the lowest possible level to achieve the desired effect.
- Establish a protocol for device examination and maintenance.

Warning:
Improper use of any element in electrosurgical systems may cause serious burn. The use of the device with safely functioning and compatible handles, cords, electrosurgical generators and accessories is the user's responsibility. Read and understand all precautions and instructions before using any electrosurgical device.

Warning:
Avoid HF output mode/power settings where maximum output voltage may exceed the stated maximum voltage.

Warning:

Do not make any modification to these devices. Modification may cause burn or injury.

Warning:

Due to concerns about the carcinogenic and infectious potential of electrosurgical byproducts (such as tissue smoke plume and aerosols which may include toxic gases and vapors, live and dead cellular material, and viruses), protective eyewear, filtration masks, and effective smoke evacuation equipment should be used in both open and laparoscopic procedures. Contact Symmetry Surgical at customerservice@symmetryurgical.com for additional information or to inquire about our smoke evacuation solutions.

Warranty:

All products are warranted to be free from defects in workmanship and materials when used for its intended purpose. Any product that proves defective in workmanship or materials will be repaired or replaced at our discretion free of charge. Repair, modification, or alteration to any product that is performed by anyone not authorized by Olsen Medical® will result in the immediate loss of warranty.

OLSEN® is a registered trademark of Olsen Medical.

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