## INSTRUCTIONS FOR USE

Suture Aid Booties  **Non-Sterile**

Aspen Surgical Products Inc.  
6945 Southbelt Drive SE  
Caledonia, MI 49316  
Phone +1 616-698-7100  
Toll-Free 888-364-7004  
aspensurgical.com

EMERGO  
Westervoortsedijk 60  
6827 AT Arnhem  
The Netherlands

For Country of Origin See Product Label

CE  
2797

## English

## Intended Use

Suture Aid Booties are intended to protect, tag and help locate sutures

## Intended User/ Patient Target Groups

Intended to be used by Healthcare Professionals on patients requiring medical procedures.

## Instructions for Use

Device is intended for limited contact duration, less than or equal to 8 hours.

## Sterilization Instructions:

- Place device into a kit or sterilization tray.
- Using Ethylene Oxide sterilization method, validate according to ISO 11135, sterilize according to the following parameters:

Cycle Parameter	Specification	Duration
<b>Preconditioning</b>		
Temperature	110°F +/- 10°F	
Humidity	65% +/- 15% RH	10.35 Hrs - 72 Hrs
<b>Chamber</b>		
Vacuum	21" +/- 1.0 inHgA	
Humidification/Dwell	22.5" +/- 1.0 HgA	60 +30/-0 minutes
Gas Inject/Dwell	10.3" +/- 0.5 Hg rise	360 +5/-0 minutes
Jacket Temperature	135" +/- 5" F	
<b>Aeration*</b>		
Temperature	120" +/- 10" F	24 Hrs - 72 Hrs Max

\*Proper aeration is essential to eliminating the risk of residuals to meet ISO 10993-7 requirements.

## Warning

Reuse of device could result in infection/contamination and/or device failure which could lead to patient harm.

## Disposal

Dispose or recycle general surgical devices per facility protocol and all applicable federal, state, regional, and/or local laws and regulations.

## Notice to Users and/or Patients in EU

Any serious incident that has occurred in relation to the device, should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

**⚠️ WARNING:** This product can expose you to Diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov). (Excluding X-Tra 053030PBX)

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#### AB'deki Kullanıcılara ve/veya Hastalara Yönelik Bildirim
Emergo Consulting (UK) Limited