Symmetry VersaGrip™ Laparoscopic Instrument

Symmetry Surgical Inc.
3034 Owen Drive
Antioch, TN 37013 USA
 1-800-251-3000
FAX: 1-615-964-5566
www.symmetriesurgical.com

EC REP

Symmetry Surgical GmbH
78532 Tuttlingen, Germany
Maybachstraße 10

+49 7461 96490
Fax +49 7461 77921

€ 2797

LCN-ERG100 Rev E

Revised 08/2020

INSTRUCTIONS FOR USE

Symmetry VersaGrip™ Laparoscopic Instrument


Indications for Use

Laparoscopic surgical manual instruments are designed to be used endoscopically through cannula to perform cutting, grasping, dissecting, retracting, and manipulating functions.

Product Description

Symmetry VersaGrip™ laparoscopic surgical manual instruments consist of a handle and a 5 mm shaft, which hold various 5 mm inserts. The shaft comes in two lengths: 33 cm and 45 cm. One handle can be used in combination with a variety of different 5 mm inserts. Instruments have monopolar posts which allow them to be used for electrosurgery.

The Symmetry VersaGrip™ laparoscopic handles have an adjustable grip. They may be used in one of the three handle positions: pistol grip, hybrid, and inline. To change the handle position, retract the position adjustment button and lock the grip in the chosen position. The handle comes in three sizes: small, medium, and large.

Warnings and Precautions

Warnings

To avoid potential serious injury to the user and the patient and/or damage to this device, the user must obey the following warnings:

1. Before use, ensure that there are no breaks, chips, cracks, scratches, tears, or missing/loose components on the shaft insulation, handle, or housing. Do not use the instrument if any of these defects are present as use of a damaged instrument could cause unintended electrosurgical burns and life-threatening complications. If damage has occurred, discontinue use and immediately return the instrument for repair or replacement.
2. The selected output setting of the ESU generator, when used with Laparoscopic Instruments, should be set in accordance with the clinician's experience, training, or by reference to appropriate clinical references.
3.  The ESU generator may emit non-ionizing radiation. Carefully read all instructions for the generator. Observe all warnings and precautions noted throughout those instructions. Failure to do so may result in complications and injury to the operator or patient.
4. Do not use in patients who have electronic implants, including but not limited to cardiac pacemakers without first consulting a qualified professional (e.g., cardiologist). A possible hazard exists because interference with the action of the electronic implant may occur or the implant may be damaged.
5. Due to concerns about the risk of potential carcinogenic and infectious diseases resulting from electrosurgical by-products (such as tissue smoke plume and aerosols), protective eyewear, filtration masks, and effective smoke evacuation equipment should be used in both open and laparoscopic procedures.
6. Use caution when handling the instruments as the insert tip may be sharp.

7. Δ Do not exceed 7000 Volts Peak-to-Peak.
8. To prevent the unintentional flow of electricity, activate the electrode when only the end of the insert is in contact with tissue. Keep active electrodes isolated from the patient when not in use.
9. To avoid inadvertent shock, place surgical electrode cables in such a way that any contact with the patient or other lead is avoided.
10. Ensure that the patient does not come into contact with metal objects that are earthed or that have an appreciable capacitance to earth, as this may result in high current densities and unwanted burns.
11. To prevent unwanted current flow, avoid skin-to-skin contact with the patient, for example, between the arms and body of the patient.
12. Do not use if the electrode is loose or damaged as doing so may cause loose parts to enter the surgical field.
13. Connect adaptors and accessories to the electrosurgical unit only when the unit is off. Failure to do so may result in an injury or electrical shock to the patient or operating room personnel.
14. Interference from high frequency surgical equipment may adversely influence operation of other electronic equipment.
15. Failure of high frequency surgical equipment could result in an unintended increase of output power.
16. Do not use if plugs are loose or improperly attached.
17. Do not use if the post on the electrode is incompletely covered by the monopolar cord.
18. Do not use if there is excessive buildup of debris on the tip during use.
19. Do not use if explosive gas concentrations are present in the area of use.
20. Do not place instruments near or in contact with flammable materials (such as gauze or surgical drapes). Instruments that are activated or hot from use may cause a fire.
21. Do not use with hybrid trocar systems, i.e., a combination of metal and plastic as doing so may result in alternate site burns due to capacitive coupling. Use only all-metal or all-plastic trocar systems.
22. Do not activate the instrument when not in contact with target tissue, as this may cause injuries due to capacitive coupling.
23. Aspirate fluid from the area before activating the instrument. Conductive fluids (e.g., blood or saline) in direct contact with or in close proximity to an active electrode may carry electrical current or heat away from target tissues, which may cause unintended burns to the patient.
24. Do not bend or pry the instrument, and avoid contacting instrument with metal objects as doing so may cause the instruments to break and leave pieces in the surgical field which may be difficult to remove. This may also result in electrical shock and/or the unintentional flow of electricity.
25. Do not apply hard rotation or "torqueing" of the instrument with the jaw under load as this may cause the jaw to break.
26. The surface of the active electrode may remain hot enough to cause burns after the RF current is deactivated.
27. When not using instruments, place them in a clean, dry, highly visible area not in contact with the patient. Inadvertent contact with the patient may result in burns.
28. Test the instrument function prior to use. If there is any sign of malfunction, the instrument should not be used.
29. Do not attempt to repair or replace any part of the instrument. Doing so could result in serious injury to the patient and/or user resulting in life-threatening injuries and complications.

Precautions:

1. Verify compatibility of all equipment and ensure grounding is not compromised.
2. Keep the active electrodes clean. Build-up of eschar may reduce the instrument's effectiveness. Do not activate the instrument while cleaning as injury to operating room personnel may result.
3. To avoid interference, always place all physiological monitoring electrodes as far as possible from the surgical electrodes. Needle monitoring electrodes are not recommended. Monitoring systems incorporating high frequency current-limiting devices are recommended.
4. Repair of a Symmetry product by any agent other than Symmetry's authorized personnel may negate any warranties, cause the instrument to fail and/or cause the jaw to break.
5. R_{ONLY} Federal law (United States of America) restricts this device to sale by or on the order of a physician.

Device Compatibility for Monopolar Surgery
Neutral electrodes

The use and proper placement of a neutral electrode (ground plate) is a key practice to ensure safe and effective monopolar surgery, particularly in the prevention of burns. Consult the neutral electrode user manual for placement and usage instructions.

High-frequency electrosurgical cables

High-frequency (HF) electrosurgical cables connect electrosurgical generators to laparoscopic manual

instruments. Do not use instruments and cables unless they are designed to fit properly. Consult the high frequency electrosurgical cable's user manual for additional instructions.

Recommended HF Electrosurgical Cables Such As:	
50-2111	Reusable Monopolar Cable
40-2110	Disposable Monopolar Cable


Electrosurgical generators

Use a compatible electrosurgical generator for monopolar surgery. Do not exceed 7000 volts peak-to-peak. Consult the generator user manual for additional instructions.

Instructions for Use

Attach one end of a compatible sterile high frequency electrosurgical cable to the laparoscopic manual instrument.

Size	Color	Approximate Glove Size	Working Length
Small	Purple	< 6.5	33cm
Small	Purple	< 6.5	45cm
Medium	Green	7-8	33cm
Medium	Green	7-8	45cm
Large	Blue	> 8.5	33cm
Large	Blue	> 8.5	45cm

- (A) Flush Port
- (B) Rotation Knob
- (C) Electrosurgical Post
- (D) Position Adjustment Button
- (E) Handle Grip
- (F) Trigger
- (G) Insert Release Button
- (H) Shaft
- (I) Insert

Repositioning

1. Retract the position adjustment button (D).
2. Rotate handle to desired position.
3. Release the position adjustment button (D).

Assembly

1. With handle in the pistol-grip position, depress the insert release button (G), move trigger (F) fully forward into the eject position.
2. While holding the trigger (F) forward, let go of insert release button (G) and slide insert (I) into distal end of shaft.
3. Allow insert (I) to push trigger (F) back while aligning up insert pegs with slots of shaft (H).
4. Once insert release button (G) engages (pops up), pull trigger (F) closed until inserts clicks into place.

Disassembly

1. Retract the position adjustment button (D).
2. Rotate the handle (E) into the pistol-grip position.
3. Depress the insert release button (G).
4. Pull trigger (F) fully open until insert (I) snaps out.
5. Pull insert (I) out at distal end of shaft (H).

Reprocessing Instructions

These reprocessing instructions are also for inserts compatible with the monopolar handles.

Warnings

- Wear appropriate protective equipment: gloves, eye protection, etc.
- These instruments are supplied non-sterile and must be cleaned and sterilized prior to the first and every subsequent use.
- Do not stack trays in the sterilizer as this may compromise the sterility of the products.

Cautions

- Do not use brushes or pads with metal or abrasive tips during manual cleaning, as permanent scoring or damage to the instrument could result.

- Allow the instruments to air cool following steam sterilization. Rapid cooling or "quenching" the instruments in a liquid will damage the instrument and void the warranty.
- Immediate-use steam sterilization (IUSS) should only be used in carefully selected clinical situations (e.g., an instrument needed for a case falls on the floor and no replacement instrument is available). IUSS should not be used due to inadequate inventory.
- Ensure no holes are blocked and tips are in the open position.

Limitations on Reprocessing

- Do not cross-sterilize the instruments. Using multiple sterilization methods may significantly reduce the performance of the instruments.
- Do not leave the instruments in solutions longer than necessary, as doing so may accelerate normal product aging wear and tear.
- Proper processing has a minimal effect on these instruments. End of life is normally determined by wear and damage due to use.
- Damage resulting from improper processing will not be covered by the warranty.

Cleaning Instructions

- Wipe excess soil from the instruments.
- Transport instruments in a tray to avoid damage.
- The instruments should be reprocessed as soon as possible following use.
- If an automated reprocessing method is used, rinse the channels with 50 mL of sterile distilled water immediately after use.

Preparation for Cleaning

1. Disassemble the instruments per the "Disassembly" section included in this manual.
2. Prepare an enzymatic detergent¹ per the manufacturer's recommendations.
3. Wipe the instruments with the detergent, using a clean cloth.
4. Immerse the instruments in the detergent. Using a syringe, inject any inside regions of the instruments with 50 mL of the detergent to ensure all parts are reached.
5. Soak the instruments in the detergent for 15 minutes.

Manual Cleaning

1. Brush

- Prepare a fresh solution of the enzymatic detergent¹ per the manufacturer's recommendations.
- Thoroughly brush the exterior of the instrument with a soft-bristled brush², focusing on any mated or rough surfaces.
- Using a syringe, inject any lumen or mated surface 5 times with 50 mL of the detergent.
- Brush any lumens 5 times from each end, using an appropriate bottle brush³.
- NOTE: An appropriate brush for the interior of the instrument is a bottle brush that has the same diameter as the lumen.
- Brush any movable parts in all extreme positions.

2. Rinse

- Rinse the instruments with RO/DI water⁴ at ambient temperature until all detergent residue is removed.
- Flush any lumens or mated surfaces 5 times.
- Rinse the devices for 1 minute until all detergent residues are removed, then continue rinsing for 30 seconds.
- Drain excess water from the instruments and dry them using a clean cloth or pressurized air.
- Visually inspect the instrument for cleanliness, paying close attention to hard-to-reach areas. If visible soil remains, repeat steps 1 and 2.

3. Soak

- Prepare a non-enzymatic detergent⁵ according to the manufacturer's recommendations.
- Fully immerse the instruments and use a syringe to inject any lumens and mated surfaces with 50 mL of the detergent.
- Soak the instrument for 15 minutes.

4. Brush

- Thoroughly brush the exterior of the instrument using a soft-bristled brush².
- Using a syringe, inject any lumen or mated surface 5 times with 50 mL of the detergent.
- Brush any lumens 5 times from each end, using an appropriate bottle brush³.
- NOTE: An appropriate brush for the interior of the instrument is a bottle brush that has the same diameter as the lumen.
- Actuate the instrument, brushing around any movable parts in all extreme positions.

5. Rinse

- Thoroughly rinse the instruments with RO/DI water⁴ until all detergent residue is removed.
- Flush any lumens or crevices 5 times.
- Rinse the devices for 1 minute or until all detergent residues are removed, then continue rinsing for 30 seconds.
- Drain excess water from the instruments by holding them at an incline, and dry them using a clean cloth or pressurized air.

Automated Cleaning

1. Prepare a fresh solution of the enzymatic detergent¹ per the manufacturer's recommendations.
2. Inject each lumen and mated surface with 50 mL prepared detergent.
3. Brush the lumen using an appropriate bottle brush³ 5 times from each end.
NOTE: An appropriate brush for the interior of the instrument is a bottle brush that has the same diameter as the lumen.
4. Rinse the instrument with treated RO/DI water⁴ at ambient temperature until there is no visible detergent residue. Rinse the devices for 2 minutes until all detergent residues are removed, then continue rinsing for 30 seconds.
5. Place the disassembled instrument into the washer on an incline to aid in drainage.
6. Program the washer using the following parameters:

Pre-Wash	
Recirculation Time	2 Minutes
Water Temperature	Cold
Detergent Type	N/A

Enzyme Wash	
Recirculation Time	2 Minutes
Water Temperature	Hot
Detergent Type	Enzymatic ¹

Wash 1	
Recirculation Time	2 Minutes
Water Temperature	66°C
Detergent Type	Non-Enzymatic ⁵

Rinse 1	
Recirculation Time	2 Minutes
Water Temperature	Hot
Detergent Type	N/A

Dry	
Recirculation Time	7 Minutes
Temperature	115°C
Detergent Type	N/A

7. Use pressurized air to aid in drying time, if necessary.
8. Inspect each instrument for visual cleanliness, pay close attention to hard-to-reach areas.
9. If visible soil remains, repeat steps 1-8.

¹ The devices were validated with the Enzymatic Detergent: Enzol® Enzymatic Detergent

² The devices were validated with an M16 soft-bristle brush

³ The devices were validated with a 5 mm x 24" bottle brush

⁴ The devices were validated with reverse osmosis (RO/DI) water

⁵ The devices were validated with the Non-Enzymatic Detergent: Prolystica 2x Concentrate Neutral Detergent

Disinfection (Optional)

Warning: Not to be used in place of sterilization.

Detailed Manual High-Level Disinfection Instructions:

1. Disassemble the instruments per the "Disassembly" section included in this manual.
2. Disinfect the instruments by submerging them in the disinfecting solution listed below. Flush lumens and crevices using a syringe with the disinfecting solution that has the following active ingredient: ≥ 2.4% glutaraldehyde⁶ with a minimum soak time of 45 minutes at 25°C
3. Thoroughly rinse and flush all parts and lumens with running, demineralized water to remove the disinfectant.
4. Dry all parts with a lint-free towel immediately after rinsing.

⁶ The devices were validated with Cidex® Activated Dialdehyde Solution

Detailed Thermal Low-Level Disinfection Instructions:

1. Disassemble the instruments per the "Disassembly" section included in this manual.
2. Place the disassembled instruments into the washer disinfector.
3. Program the washer disinfector using the following parameters:

Thermal Rinse	
Recirculation Time	1 Minute
Temperature	90°C

Inspection/Lubrication

Carefully inspect each device for damage and wear. Check all moving parts to confirm the device smoothly operates throughout the intended range of motion.

After cleaning and before sterilization, instruments with moving parts (e.g. hinges, connectors, box-locks, sliding, or rotating parts) should be lubricated with instrument milk or a

water-soluble lubricant intended for medical device applications and compatible with the sterilization modality. Always follow the lubricant manufacturer's instructions for dilution, shelf life, and application method.

Sterilization

Warnings – When using a sterilization tray:

- Consult the tray manufacturer's instructions to confirm compatibility with Symmetry VersaGrip™ laparoscopic instruments (including but not limited to the lumen length and diameter).
- ONLY use pre-vacuum steam sterilization cycles. Symmetry VersaGrip™ laparoscopic instruments are not validated for other sterilization methods.

1. Inspect all components for cleanliness. If fluid or tissue buildup is present, repeat the above cleaning and disinfection procedures.
2. Inspect for burns, cuts, and abrasions in the electrical insulation on the insert and/or the handle for instruments equipped with electrosurgical capabilities as these defects could cause unintended electrosurgical burns and life-threatening complications.
3. Inspect all jaws and cutting edges for damage.
4. Ensure instrument jaws align properly and open and close completely.
5. Instrument may be sterilized assembled as long as the cleaning process prior to sterilization was completed as a dissembled unit.
6. Instruments may be placed in a specialized tray or containment device for protection. Follow the instructions provided by the manufacturer.
7. Perform sterilization per the following specifications. Failure to do so will compromise sterility.

Steam Sterilization

Warnings: Drying time depends on several variables, including altitude, humidity, type of wrap, preconditioning, size of chamber, mass of load, material of load, and placement in chamber. Users must verify that drying time set in their autoclave yields dry surgical equipment.

NOTE: Open all flush ports for sterilization or incomplete sterilization may result.

Pre-Vacuum (US)	
Wrapping	Double ⁷
Temperature	132°C (270°F)
Sterilization Time	4 minutes
Dry Time	30 minutes

Pre-Vacuum (Europe)	
Wrapping	Double ⁷
Temperature	134°C (274°F)
Sterilization Time	3 minutes
Dry Time	30 minutes

⁷ Sterilization was validated using two layers of 1-ply KimGuard KC600 polypropylene wrap.

Immediate-Use

 **Warning:** Immediate-use steam sterilization is intended only for emergency situations. In the event that immediate-use reprocessing is required, the following instructions should be used.

"Immediate-Use" Pre-Vacuum	
Temperature	132°C (270°F)
Sterilization Time	4 minutes

Warranty

Symmetry Surgical Inc. guarantees every surgical device bearing the Symmetry VersaGrip™ brand name to be free from defects in workmanship and materials when used pursuant to its instructions for use and for its intended surgical purpose for the life of the product. Any other express or implied warranties, including warranties of merchantability or fitness for a particular purpose, are hereby disclaimed and there are no warranties that extend beyond those set forth herein. Abuse or misuse of the product or failure to comply with the instructions for use shall void this warranty. Any Symmetry VersaGrip™ branded device that fails to meet the foregoing warranty will be replaced or repaired, at Symmetry's election, at no charge.

The device's electrical insulation has been verified to withstand 250 reprocessing cycles, although the service life in clinical practice will depend on the intraoperative usage and actual reprocessing practices. Suitability of this medical device for any surgical procedure should be determined by the user in conformance with the manufacturer's instructions for use.

Storage

Store packaged and sterilized instruments in a dry, clean, and dust-free environment at room temperatures (2 – 25°C).

Service/Disposal

Adhere to national regulations when disposing of or recycling the product, its components and its packaging. For service outside the United States, contact your local Symmetry Surgical representative. Inside the United States, call Symmetry Surgical Customer Service: 1-800-251-3000.

Labeling Symbols

	Authorized representative in the European Community
	Manufacturer
	Date of Manufacture
	Catalog Number
	Batch Code
	Quantity
	Consult Instructions for Use
	Caution
	Conformance to 93/42/EEC
	Dangerous Voltage
	Non-Sterile
	Federal law restricts this device to sale by or on the order of a physician
	Emits non-ionizing radiation

 In USA Contact
Symmetry Surgical Inc.
3034 Owen Drive
Antioch, TN 37013 USA
Phone: 1-800-251-3000
Fax: 1-615-964-5566
www.symmetrysurgical.com

 In Europe Contact
Symmetry Surgical GmbH
Maybachstraße 10
78532 Tuttlingen, Germany
Phone: +49-7461-96490
Fax: +49-7461-77921

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