# Bovie<sup>®</sup>

**Bovie Medical Corporation**  
5115 Ulmerton Road  
Clearwater, Florida 33760 USA

Emergo Europe  
Westervoortsedijk 60  
6827 AT Arnhem  
The Netherlands  
SRN #: NL-AR-000000116  
+31.70.345.8570  
EmergoEurope@ul.com

CE 0413

# Bovie®

## Arthroscopic Electrodes

**REF** AR00, AR01, AR02, AR03


**Rx ONLY**

**STORAGE**  
10 to 75 % non-condensing  
**OPERATE**  
30 to 75 % non-condensing  
**TRANSPORT**  
10 to 100 % including condensation

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---

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## **EN Arthroscopic Electrodes**

Bovie Medical Corporation electrodes have a standard 2.36 mm (3/32") diameter stainless steel shaft and plastic insulator, and fit Bovie® and most other electrosurgical pencils.

### **Intended Use:**

Bovie electrodes are used for cutting and coagulation of tissue.

at least one minute before and during use of electrosurgery.

### **WARNINGS:**

This device has been specifically designed for use in arthroscopy. Do not use for other procedures.

**WARNING:** This product is sterilized utilizing Ethylene Oxide. These products can expose you to chemicals including Ethylene Oxide (EtO), which is known to the state of California to cause cancer and/or birth defects or other reproductive harm. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

**Danger – Explosion Hazard:** Do not use in the presence of flammable anesthetics.

**Fire Hazard:** Both oxygen (O<sub>2</sub>) and nitrous oxide (N<sub>2</sub>O) support combustion. Avoid O<sub>2</sub> and N<sub>2</sub>O enriched atmospheres. Enriched atmospheres may result in fires and burns to patients or surgical personnel.

The metal shaft of the electrode must fit completely and securely into the pencil. If metal is visible, do not use the electrode. Exposed metal may result in electric shock or burn to the patient or surgical personnel.

Confirm proper electrosurgical settings prior to and during a procedure. Use the lowest power settings to achieve the desired effect. If increased power settings are requested, check the patient return electrode and all accessory connections before major power setting adjustments.

**Fire/Explosion Hazard:** The following substances contribute to increased fire and explosion hazards in the operating room:

- Oxygen enriched environments.
- Oxidizing agents such as nitrous oxide (N<sub>2</sub>O) atmospheres.
  - \* Verify all anesthesia circuit connections are leak free before and during use of electrosurgery.
  - \* Verify endotracheal tubes are leak free and that the cuff seals properly to prevent oxygen leaks.
  - \* If an uncuffed tube is in use, pack the throat with wet sponges around the uncuffed tube.
  - \* If possible, stop supplemental oxygen

- \* Activate the electrosurgical unit only after vapors from skin prep solutions and tinctures have dissipated.
- Naturally occurring flammable gases (such as methane) that may accumulate in body cavities.

**Fire Hazard:** Always place the active electrode in a clean, dry, insulated safety holster when not in use.

- Electrosurgical accessories that are activated or hot from use can cause unintended burns to the patient or surgical personnel.
- Electrosurgical accessories may cause fire or burn if placed close to or in contact with flammable materials such as gauze or surgical drapes. Place longer electrodes such as extended electrodes away from the patient and drapes.

**Fire Hazard:** The sparking and heating associated with electrosurgery can provide an ignition source.

### **Observe fire precautions at all times:**

- When using electrosurgery in the presence of gases or flammable substances, prevent pooling of fluids and the accumulation of gases under surgical drapes. Tent drapes to allow vapors to mix with room air before using electrosurgery.
- Tissue buildup (eschar) on the tip of an active electrode poses a fire hazard, especially in oxygen enriched environments such as in throat or mouth procedures. With sufficient heating, eschar can become a glowing ember and pose a fire hazard both as an ignition source and as a fuel. Keep the electrode clean and free of all debris.
- Facial and other body hair is flammable. Water soluble surgical lubricating jelly may be used to cover hair close to the surgical site to decrease flammability.

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This medical device cannot be effectively cleaned and/or sterilized by the user and therefore cannot be safely reused. It is intended for single use only. Any attempt by the user to clean and resterilize this device may result in bioincompatibility, infection, or other risks of device failure to the patient.

## WARNINGS:

Conductive fluids (e.g., blood or saline) in direct contact with an active electrode or in close proximity to any active accessory may carry electrical current and cause unintended burns to the patient. This can happen as a result of either direct coupling with the active electrode or capacitive coupling between the active electrode and the external surface of the electrode insulation. Therefore, to prevent unintended burns in the presence of conductive fluids:

- Always keep the external surface of the active electrode away from adjacent tissue while activating the electrosurgical generator.
- Clear conductive fluid from the electrode before activating the electrosurgical pencil.

Some surgeons may elect to buzz the hemostat during surgical procedures. It is not recommended, and the hazards of such a practice probably cannot be eliminated. Burns to the surgeon's hands may result. To minimize the risk, take these precautions:

- Do not "buzz the hemostat" with a Needle electrode.
- Do not lean on the patient, the table, or the retractors while buzzing the hemostat.
- Activate cut rather than coag. Cut has a lower voltage than coag.
- Use the lowest power setting possible for the minimum time necessary to achieve hemostasis.
- Activate the generator after the accessory makes contact with the hemostat. Do not arc to the hemostat.
- Firmly grasp as much of the hemostat as possible before activating the generator. This disperses the current over a larger area and minimizes current concentration at the finger tips.
- "Buzz the hemostat" below hand level (as close as possible to the patient) to reduce the opportunity for current to follow alternate paths through the surgeon's hands.
- When using a coated or nonstick blade electrode, place the edge of the electrode against the hemostat or other metal instrument.

## CAUTION:

Always use the lowest power setting that achieves the desired surgical effect. Use the active electrode for the minimum time necessary in order to reduce the possibility of unintended burn injury.

Do not activate electrodes while in contact with or near other instruments, including cannulae. Localized burns to the patient or physician may result.

Before use, examine the electrosurgical unit and accessories for defects. Do not use cables or accessories with damaged (cracked, burned, or taped) insulation or connectors. Localized burns to the patient or physician may result.

Do not modify or add to the insulation of active electrodes.

Device is rated for 3kVpeak.

Activate the electrosurgical unit only when you are ready to deliver electrosurgical current and the active tip is in view and near target tissue (especially if looking through an an endoscope).

Deactivate the electrosurgical unit before the tip leaves the surgical site. The electrodes are intended for single use only. Discard after use in a puncture resistant container to prevent injury. These electrodes are not designed to withstand sterilization. DO NOT RESTERILIZE.

## IMPORTANT

Clean the electrode often with moist gauze or other material appropriate for the electrode type.

Bovie Medical Catalog No. ESSP Scratch Pad is a sterile single use cleaning aid for non loop electrodes and is useful in removing eschar during electrosurgical procedures.

## OPERATING PRINCIPAL

### Site of use

- No Restrictions

### Site Condition

- Aseptic

### Intended User Profile

- Education
  - o Restricted for use by or on the order of a physician
- Knowledge
  - o Understands electrosurgery and electrosurgical techniques
  - o Reads and understands the supplied Instruction for Use or Device Labeling
  - o Understands hygiene
- Experience

- o Some training on techniques or training under surveillance/supervision
- Permissible impairments
  - o Mild reading vision impairment or corrected vision to 20/20

#### Intended Conditions for Use

- Doctor office, Surgery center or Hospital, intended for professional use only
- Lit, clear and unobstructed view at point of use
- Ambient Temperature 10° C to 40° C

#### Frequency of Use

- Single Use - Single patient multiple activation

#### Disposal

- Biohazard – Dispose of electrode according to facility procedure

### Install the Electrode

1. Ensure the pencil is not connected to the generator, or the generator is Off or in the Standby mode, if available.
2. Grasp the insulating sleeve on the electrode. Remove tip protector. Insert the electrode into the pencil.
3. Ensure the electrode is fully inserted into the pencil. The shank and insulating sleeve should fit securely into the active accessory. If the shank and/or insulating sleeve does not fit, or the insulation will not insert 3.2 mm ( $\frac{1}{8}$  in.), do not use this electrode/accessory combination.
4. Refer to the generator instruction manual for proper procedures for connecting the active accessory to the generator.
5. A tip protector covers the end of some electrodes. If a tip protector is present, remove and discard it before use.  
Do not use it for insulation purposes.

### Maximum Power

#### Notice:

*Exceeding maximum power settings may result in patient injury or product damage.*

Bovie Medical Corporation recommends the following maximum power settings:

<u>Mode</u>	<u>Power Setting</u>
Coag	30 W
Pure Cut or Blend	35 W

*Sterility is guaranteed unless the package is opened or damaged.*

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