<div> <div> <div>For product information, contact:</div> </div> <div> <div>Yerligere produktoplysninger fås ved at kontakte:</div> </div> </div> <div> <div><div>Neeem voor productinformatie contact op met:</div></div> </div> <div> <div>Pour obtenir des informations sur le produit, veuillez contacter :</div> </div> <div> <div>Für Produktinformationen kontaktieren Sie:</div> </div> <div> <div> Για πληροφορίες σχετικά με το προϊόν, επικοινωνήστε με:</div> </div> <div> <div>Per informazioni sul prodotto, contattate:</div> </div> <div> <div>製品情報、お問い合わせ先:</div> </div> <div> <div>Para obter informações sobre o produto, contacte:</div> </div> <div> <div>Para obtener información del producto, póngase en contacto con:</div> </div> <div> <div>For produktinformaon, kontakta:</div> </div>	<div> <div> <div></div> </div> <div> <div></div> </div> </div>
---	--

<div> <div> <div>3230 Commerce Center Place</div> <div>Louisville, Kentucky, 40211 USA</div> <div>(502) 813-5022 <- (800) 297-6344</div> </div> </div> <div> <div>See our products at:</div> <div>Se vores produkter på:</div> <div>Ga voor onze producten naar:</div> <div>Découvrez nos produits sur:</div> <div>Siehe unsere Produkte auf:</div> <div>Δείτε τα προϊόντα μας στην ιστοσελίδα:</div> <div>È possibile vedere i nostri prodotti all'indirizzo web:</div> <div>当社の製品は以下にてご覧いただけます</div> <div>Consulte os nossos produtos em:</div> <div>Consulte nuestros productos en:</div> <div>Se våra produkter på:</div> <div>http://www.symmetryusa.com</div> </div>	<div> <div> <div></div> </div> <div> <div></div> </div> </div>
--	--

CE 1639

EC	REP
----	-----

Symmetry Surgical GmbH
Maybachstraße 10, 78532 Tuttlingen, Germany
+49 7461 96490
Made in U.S.A.
Fremstillet i U.S.A.
Vervaardigd in de VS
Fabriqué aux États-Unis
In den USA hergestellt
Κατασκευασμένο στις ΗΠΑ
Producto negli USA
(作成：印刷(米国))
Fabricado nos EUA
Fabricado en EE.UU.
Tilberkerd i USA

RE ONLY

OLSEN MEDICAL® Single Use Bayonet Bipolar Irrigating Forceps and Precision Single-Use Non-Stick Bipolar Irrigating Forceps
OLSEN MEDICAL® Bipolarøj bajonetskylletang til engangsbrug og Precision klabefri bipolar skylletang til engangsbrug
Bipolaire bajonetspoeltang voor eenmalig gebruik en bipolaire non-stick precisie-bajonetspoeltang voor eenmalig gebruik van OLSEN MEDICAL®
Pinces bipolaires d'irrigation à baïonnette à usage unique et pinces bipolaires d'irrigation de précision à revêtement non-adhésif et à usage unique d'OLSEN MEDICAL®

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ENGLISH

Description:
The OLSEN MEDICAL® Single Use Bayonet Bipolar Irrigating Forceps are packaged sterile, unless otherwise noted. It consists of an inlet pathway running into a tube attached to the left hand side of the bayonet forceps. This gives the physician a visual of the fluid pathway which facilitates the local fluid irrigation through the forceps. The device is for Single Use Only and is designed for use with bipolar electro-surgical units.

Indications:
The OLSEN MEDICAL® Single Use Bayonet Bipolar Irrigating Forceps and Precision Single-Use Non-Stick Bipolar Irrigating Forceps (Catalog No. 20-0xxx) are single use products sold sterile and are intended for use in electro-surgery for coagulation and irrigation of tissue. This device is intended for use with the OLSEN MEDICAL® Integrated Irrigation Tubing and Bipolar Cord Set or similar design of Bipolar Cord and Irrigation Tubing. The Single Use Bayonet Bipolar Irrigating Forceps are intended for use with a maximum voltage of 500 Volts (Vpeak).

Handling:
1. Inspect package for damage before use. Discard if damaged or seal is broken.
2. Inspect device after opening in a sterile field. Discard if damaged.
3. Remove and discard the tip protector on all forceps before use.
4. Device use and application to be determined by physician.

Storage:
These devices should be stored in their original labeled boxes or in protected storage to prevent damage during handling. Take precautions to ensure that the sterile devices are stored in areas where it is completely dry, and protected from extremes of temperature and humidity.

Sterility:
Sterility is guaranteed up to the expiration date unless package is opened or seal is broken.
Sterilization:
Single Use Only (intended to be used on one patient during one procedure). Dispose after use. Do not resterilize.

Caution:
1. Inspect every device before use. Damaged electro-surgical devices may cause electrical burns to the patient or physician.
2. Use a wet sponge to remove eschar build up from the forceps tips.
3. Active electro-surgical devices must be placed away from the patient when not in use. A holster is recommended.
4. Consult RF generator Owner’s Manual. Use only in the active RF generator output connections.
5. Electromagnetic compatibility (EMC) – Olsen Medical devices are considered medical electrical equipment. Medical electrical equipment may require precautions regarding EMC and may need to be installed or put into service according to the EMC information provided in the operator’s manual of the generator.
6. Do not use electro-surgical equipment unless properly trained to use it in the specific procedure being undertaken. Use by physicians without such training may result in serious, unintended patient injury
7. Do not perform electro-surgery in the presence of flammable anesthetics or other flammable gases, fluids or objects, or in the presence of oxidizing agents, as fire could result.
8. Avoid unnecessary or prolonged activation. Activate only when touching or immediately next to the intended point of contact.
9. Keep the voltage/power settings at the lowest possible level to achieve the desired effect.
10. After use, safely discard in accordance with established procedures for biohazardous waste.

Warning:
Implosion hazard of any element in electro-surgical systems may cause serious burn, shock, and fire or explosive hazard. The use of the device with safely functioning and compatible handpieces, cords, electro-surgical generators and accessories is the user’s responsibility. Read and understand all precautions and instructions before attempting to use any electro-surgical devices. This product is for SINGLE USE ONLY. DO NOT RESTERILIZE. Use aseptic technique in all phases of handling. Olsen Medical Single Use devices have not been designed to undergo or withstand any form of alteration, such as disassembly, cleaning or resterilization, after a single patient use. Reuse can potentially compromise device performance and in any usage beyond the design intent of this single use device may result in unpredictable loss of functionality.

Warning:
This product is sterilized utilizing Ethylene Oxide. These products can expose you to chemicals including Ethylene Oxide (EtO), which is known to the state of California to cause cancer and/or birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

Surgical Assembly Instructions:
1. Ensure that the Irrigation Tubing and Bipolar Cord set is prepared for connection to the instrument. If the cord set is included in a kit make sure that the cord is prepared prior to connecting the forceps.
2. Remove the Single Use Bayonet Bipolar Irrigating Forceps from the peel pouch using aseptic technique.
3. Connect the twin pin connector of the OLSEN MEDICAL® Single Use Bayonet Bipolar Irrigating Forceps into the female socket connector of the OLSEN MEDICAL® Integrated Irrigation Tubing and Bipolar Cord Set or of a similar Olsen design. Connect the set’s male LUER-LOK™ connector to the female irrigation port on the forceps.

Warranty:
All products are warranted to be free from defects in workmanship and materials when used for its intended purpose. Any product that proves defective in workmanship or materials will be repaired or replaced at our discretion free of charge. Repair, modification, or alteration to any product that is performed by anyone not authorized by Olsen Medical® will result in the immediate loss of warranty.
DANSK

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