Vessel Loops

aspensurgical.

STERILE EO R_x ONLY MD

EC REP

Emergo Europe
Prinsessegracht 20
2514 AP The Hague
The Netherlands

Made in USA

CE
2797

Aspen Surgical Products, Inc.
6945 Southbelt Drive SE
Caledonia, MI 49316 USA
Phone 616-698-7100
Toll-Free 888-364-7004
aspensurgical.com

English

Intended use:

Vessel loops are intended to be used to occlude, retract, and identify arteries, veins, tendons, and nerves in surgical procedures.

Intended User/ Patient Target Groups:

Intended to be used by Healthcare Professionals on patients requiring medical procedures.

Directions for use:

- After removing loops from package, the surgeon may do as follows:
 - Retract-wrap around vessel. Clamp loop to drape or hang outside wound.
 - Occluding-wrap loop around vessel twice. Clamp adjacent to vessel to secure.
 - Identify-wrap loop around vessel, nerve and ureter or intestine.
- Dispose or recycle general surgical devices per facility protocol and all applicable federal, state, regional, and/or local laws and regulations.

Precautions:

Device is single use only. Device is intended for limited contact duration, less than or equal to 8 hours. Reuse of device could result in infection/contamination and/or device failure which could lead to patient harm. Do not use product if product damaged or particulate is found in sterile packaging

Notice to Users and/or Patients in the UE:

Any serious incident that has occurred in relation to the device, should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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