Bovie



Laparoscopic Electrodes with 4mm (.157") Adaptor

REF

LB02, LC04, LH01, LH08, LJ07, LL06, LN03, LS05

R_x ONLY

EN Not Made With Natural Rubber Latex
DA Ikke fremstillet med naturlig gummilatex
NL Niet vervaardigd uit latex van natuurlijk rubber
FI Valmistuksessa ei ole käyretty luonnonkumilateksia
FR Exempt de latex de caoutchouc naturel
DE Ohne Naturkautschuklatex hergestellt
IT Non realizzato con lattice di gomma naturale
PT Não fabricado com látex de borracha natural
ES No está fabricado con látex de caucho natural
SV Ej tillverkad med naturgummi
EL Δεν είναι κατασκευασμένο από φυσικό καουτσούκ λάτεξ
CN 非天然橡膠乳膠製成
JA 天然ゴムラテックスで作られていません

CE 0413

Bovie

Bovie Medical Corporation
5115 Ulmertown Road
Clearwater, Florida 33760
USA

EC REP

Emergo Europe
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

MC-18846 Rev. 10
2023-05-30

EN Laparoscopic Electrodes

Bovie Medical Corporation Laparoscopic Electrodes have a standard 2.36mm (³/₃₂" diameter shaft and fit Bovie®, and most other, electro-surgical pencils. The 4mm (.157") Adaptor provided with the Laparoscopic Electrode is used to adapt the electrode with an Electrosurgical Cable with a 4mm socket. The device is rated for 3.5 kVpeak.

Indications for Use/ Medical Purpose

Bovie electrodes are used for cutting and coagulation of tissue.

WARNINGS:

Danger – Explosion Hazard: Do not use electrosurgery in the presence of flammable anesthetics.

Fire Hazard: Both oxygen (O₂) and Nitrous Oxide (N₂O) support combustion. Avoid O₂ and N₂O enriched atmospheres. Enriched atmospheres may result in fires and burns to patients or surgical personnel.

The metal shaft of the electrode must fit completely and securely into either the pencil or adaptor. If the adaptor is used, the adaptor must fit completely and securely on the cable. If metal is visible, do not use the electrode. Exposed metal may result in electric shock or burn to the patient or surgical personnel.

Confirm proper electrosurgical settings prior to and during a procedure. Use the lowest power settings to achieve the desired effect. If increased power settings are requested, check the patient return electrode and all accessory connections before major power setting adjustments.

Fire Hazard: Always place the active Electrode in a clean, dry, insulated safety holster when not in use.

Electrosurgical accessories that are activated or hot from use can cause unintended burns to the patient or surgical personnel.

Electrosurgical accessories may cause fire or burn if placed close to or in contact with flammable materials such as gauze or surgical drapes. Place longer electrodes such as extended electrodes away from the patient and drapes.

Fire/Explosion Hazard: The following substances contribute to increased fire and explosion hazards in the operating room:

- Oxygen enriched environments
- Oxidizing agents such as nitrous oxide (N₂O) atmospheres
 - * Verify all anesthesia circuit connections are leak free before and during use of electrosurgery.
 - * Verify endotracheal tubes are leak free and that the cuff seals properly to prevent oxygen leaks.
 - * If an uncuffed tube is in use, pack the throat with wet sponges around the uncuffed tube.
 - * If possible, stop supplemental oxygen at least one minute before and during use of electrosurgery.
 - * Alcohol-based skin prepping agents and tinctures.
 - * Activate the electrosurgical unit only after vapors from skin prep solutions and tinctures have dissipated.
 - * Naturally occurring flammable gases (such as methane) that may accumulate in body cavities.

Fire Hazard: The sparking and heating associated with electrosurgery can provide an ignition source.

Observe fire precautions at all times:

- When using electrosurgery in the presence of gases or flammable substances, prevent pooling of fluids and the accumulation of gases under surgical drapes. Tent drapes to allow vapors to mix with room air before using electrosurgery.
- Tissue buildup (eschar) on the tip of an active electrode poses a fire hazard, especially in oxygen enriched environments. With sufficient heating, eschar can become a glowing ember and pose a fire hazard both as an ignition source and as a fuel. Keep the electrode clean and free of all debris.
- Facial and other body hair is flammable. Water soluble surgical lubricating jelly may be used to cover hair close to the surgical site to decrease flammability.

This medical device cannot be effectively cleaned and/or sterilized by the user and therefore cannot be safely reused. It is intended for single use only. Any attempt by the user to clean and resterilize this device may result in bioincompatibility, infection to the patient, or other risks of device failure.

Conductive fluids (e.g., blood or saline) in direct contact with an active electrode or in close proximity to any active accessory may disperse electrical current and cause unintended burns to the patient. This can happen as a result of either direct coupling with the active electrode or capacitive coupling between the active electrode and the external surface of the electrode insulation. Therefore, to prevent unintended burns in the presence of conductive fluids:

- Always keep the external surface of the active electrode away from adjacent tissue while activating the electrosurgical generator.
- Clear conductive fluid from the electrode before activating the electrosurgical pencil.

Any active accessory may disperse electrical current and cause unintended burns to the patient. This can happen as a result of either direct coupling (active electrode is accidentally activated near another metal instrument causing the secondary instrument to become energized) with the active electrode or capacitive coupling (capacitive coupling transfers current to non-target tissue causing stray electrosurgical burns) between the active electrode and the external surface of the electrode insulation. Therefore, to prevent unintended burns to the patient:

- Avoid accidentally activating the generator while the active electrode is near another metal instrument.
- Do not activate the electrode while in contact with other instruments.

WARNING: This product is sterilized utilizing Ethylene Oxide. These products can expose you to chemicals including Ethylene Oxide (EtO), which is known to the state of California to cause cancer and/or birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

CAUTIONS:

Federal (USA) law restricts this device to sale by or on the order of a physician.

Always use the lowest power setting that achieves the desired surgical effect. Use the active electrode for the minimum time necessary in order to reduce the possibility of unintended burn injury.

Do not contact or arc the activated electrode to metal surfaces. Blade or needle breakage may result.

Confirm that the active electrode is properly connected to the active accessory, and that the accessory is firmly connected to the correct power output receptacle on the electrosurgical generator.

Before use, examine the electrosurgical unit and accessories for defects. Do not use cables or accessories with damaged (cracked, burned, or taped) insulation or connectors. Localized burns to the patient or physician may result.

Do not modify or add to the insulation of active electrodes.

Activate the electrosurgical unit only when you are ready to deliver electrosurgical current and the active tip is in view and near target tissue. Deactivate the electrosurgical unit before the tip leaves the surgical site.

The electrodes are intended for single use only. Discard after use in a puncture resistant container to prevent injury. These electrodes are not designed to withstand resterilization. Do not resterilize.

For monopolar electrosurgical procedures, ensure that the patient return electrode is appropriately selected, properly applied to the patient, and connected to the electrosurgical generator.

When applying the patient return electrode, ensure that the entire pad adheres to the patient. Do not reuse single use patient return electrodes. Do not use the return electrode if the packaging is damaged, or if the gel / adhesive material has dried.

If the patient has moved or is repositioned after the patient return electrode is applied, ensure that the return electrode firmly contacts the skin and that all connections are intact.

Avoid spilling any solutions on the return electrode. If this occurs, remove the return electrode, dry the skin, and apply a new return electrode. A moist environment may cause the return electrode to peel off.

CAUTIONS: For Coated Electrodes Only This electrode has a coating to reduce the sticking of eschar. Damage to the coating may occur if a scratch pad or other abrasive or sharp object is used to clean the electrode.

Using the coated electrode at a high power setting may cause damage to the coating.

NOTICE:

Do not modify the electrode tip. Modifications to the tip may result in tip breakage or other damage.

IMPORTANT:

Clean the electrode often with moist gauze or other material appropriate for the electrode type.

Maximum Power Settings


- High power settings may cause degradation of the electrode and reduce or impair performance. The manufacturer recommends that the lowest possible power settings be used to achieve the desired tissue effect.
- Power settings guidelines may vary due to differences in surgical techniques, patients, electrodes, and surgical set-up. If you are not familiar with the system or the proper setting is not known, use a lower than normal setting until the desired clinical effect is achieved.
- Confirm proper electrocautery generator settings before proceeding with surgery. Use the lowest possible setting to achieve the desired effect. (This is important due to the potential of inadvertent burning at high voltage.)

Installing the Electrode with an Adaptor

The Bovie® Laparoscopic Electrode has a 4mm (.157") Adaptor attached for use with an electrocautery cable with a 4mm socket.

1. When installing the electrode with adaptor, ensure the adaptor is completely and securely inserted on the electrode. If metal is visible, do not use the electrode. Exposed metal may result in electrical shock or burn to the patient or surgical staff.
2. Ensure the electrocautery cable is not connected to the generator, or the generator is Off or in the Standby mode, if available.
3. Grasp the insulating sleeve on the electrode adaptor and insert adaptor into the cable socket.
4. Ensure the adaptor is fully inserted into the cable socket. If metal is visible, do not use the adaptor. Exposed metal may result in electrical shock or burn to the patient or surgical staff.
5. Refer to the generator instruction manual for proper procedures for connecting the active accessory to the generator.

Installing the Electrode without an Adaptor

The Bovie® Laparoscopic Electrode has a 4mm (.157") Adaptor attached for use with an electrocautery cable with a 4mm socket. If the electrode is to be used with a standard electrocautery pencil, remove the adaptor and discard.

1. Ensure the pencil is not connected to the generator, or the generator is Off or in the Standby mode, if available.
2. Grasp the insulating sleeve of the electrode and insert the electrode into the pencil.
3. Ensure the electrode is fully inserted into the pencil. The shank and insulating sleeve should fit securely into the pencil. If the shank and/or insulating sleeve does not fit, or insulation will not insert 3.2mm (1/8"), do not use electrode/pencil combination.
4. Refer to the generator instruction manual for proper procedures for connecting the active accessory to the generator.

Operating Principal

Site of use

- No Restrictions

Site Condition

- Aseptic

Intended User Profile

- Education
 - o Restricted for use by or on the order of a physician
- Knowledge
 - o Understands electrocautery and electrocautery techniques
 - o Reads and understands the supplied Instruction for Use or Device Labeling
 - o Understands hygiene
- Experience
 - o Some training on techniques or training under surveillance/supervision
- Permissible impairments
 - o Mild reading vision impairment or corrected vision to 20/20

Intended Conditions for Use


- Doctor office, Surgery center or Hospital, intended for professional use only
- Lit, clear and unobstructed view at point of use
- Ambient Temperature 10° C to 40° C

Frequency of Use

- Single Use - Single patient multiple activation

Disposal

- Biohazard – Dispose of electrode according to facility procedures

Storage and Handling:

Device must be stored in a clean, cool and dry place. It must be handled with care to avoid damage to the packaging and its content during transportation and storage. Instrument must not be in contact with acids or other corrosive liquids.

Sterility is guaranteed unless the package is opened or damaged.

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