# AlerTox•Lactosa

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**International:**

**Hygiena International**

8, Woodshots Meadow  
Watford, Hertfordshire  
WD18 8YU, UK  
+44 (0)1923-818821

**Hygiena (Shanghai) Trading Co., Ltd.**

Rm.7K, No.518, Shangcheng Rd.  
Pudong New District  
Shanghai, China 200120  
+86 21-5132-1081, +86 21-5132-1077,  
+86 21-5132-1078

**Hygiena Diagnóstica España S.L.**

P. I. Parque Plata, Calle Cañada Real 31-35,  
41900, Camas, Sevilla, Spain  
+34 954-08-1276

[www.hygiena.com](http://www.hygiena.com)

[enquiries@hygiena.com](mailto:enquiries@hygiena.com)

# AlerTox•Lactose

## Samples with low carbohydrate content

Thin-layer chromatography rapid test for the detection of lactose in foods and beverages with a low content (less than 8%) of carbohydrates, rinsing water and on surfaces.

## 1. Test basis

The test is based on the principles of thin-layer chromatography: separation and detection. They are applied to separate and detect lactose in foods and beverages with a low content (less than 8%) of carbohydrates, rinsing water and on surfaces. It consists of three phases common to all sample types: matrix homogenization and lactose extraction, compound separation and visualization.

## 2. Supplied materials (10 tests)

Reference	Component	Units
ASY3179	Extraction Solution A (5 mL)	10
ASY3168	Filter paper	10
BOT3000	Collection tube	10
ASY3186	Tube containing AlerTox Lactose separation plates	1 (10 plates)
ASY3166	Separation chamber	1
ASY3181	Separation Solution (4 mL)	10
ASY3177	Visualization Solution A (1.9 mL)	10
ASY3178	Visualization Solution B (1.5 mL)	1
ASY3176	Sprayer for Visualization Solutions A and B	10
ASY3171	Lactose Standard (1.5 mL)	1

## 3. Necessary material not supplied

- AlerTox Heater (MCH3001)
- Scale (precision 0.1 g)
- Funnel
- Homogenization system for solids (grinder, mortar...)
- Micropipettes (2-20 µL, 20-200 µL and 100-1000 µL) and disposable tips
- Disposable gloves
- Mask
- Cotton swabs (for surface analysis)
- Timer

## 4. Optional material

- Vortex mixer
- Extraction hood

**NOTE:** If you wish to analyze more than one sample at the same time, you will need to have as many separation chambers (ASY3166) as samples you have.

## 5. Storage conditions

AlerTox Lactose must be stored at 15 – 25 °C (59 – 77 °F). **Never store AlerTox Lactose at temperatures above 30 °C (86 °F) or below 2 °C (36 °F).**

## 6. Detection limits

AlerTox Lactose is a qualitative test that allows the detection of lactose in a sample, at the established working dilutions. The limit of detection of AlerTox Lactose is 100 ppm of lactose.

## 7. Precautions

- Carefully read this manual and the Material Safety Data Sheet (MSDS) before performing the test.
- This kit is designed for professional use only.
- Personal protection equipment should be worn during the separation and visualization phases, as indicated in the MSDS.
- Maintain optimal hygienic conditions during the test to avoid contaminations.
- Do not use materials/reagents from different kits, nor use any of them beyond their expiration date.

## 8. Instructions

**WARNING!** It is advisable to wear a mask during the whole procedure, and, if possible, manipulate the components of this kit under an extraction hood.

**WARNING!** The AlerTox Heater has no standby mode. This means its surface will be hot (110 °C/230 °F) during the entire test procedure.

- 8.1.** Switch on the AlerTox Heater and set it to 110 °C (230 °F). Allow AlerTox Heater to reach thermal stability before placing any AlerTox separation plate on its surface (min. 20 min).
- 8.2.** Pour the content of one of the supplied Separation Solution tubes in the separation chamber and close the chamber. **Incubate for at least 20 minutes before performing the test.**

### 8.3. Sample preparation (foods: liquids and semi-liquids, solids)

#### 8.3.1. Homogenize the sample

- Liquids/semi-liquids: shake the sample to ensure it is completely homogeneous.
- Solids: crush the sample as finely as possible to obtain a homogenous powder.

**8.3.2. Add 2 mL (liquids) or 1 g (solids and semi-liquids) of homogenized sample** to one of the tubes with Extraction Solution and close it.

**8.3.3.** Shake the tube vigorously (manually or using a vortex mixer) for at least 2 minutes.

**8.3.4.** Filter the sample using one of the provided filter papers into one collection tube with yellow cap.

### 8.4. Sample preparation (rinse waters, surfaces)

#### 8.4.1. Surfaces

1. Open one of the tubes with Extraction Solution and moisten the tip of the swab with the solution. Rub the swab in a zigzag pattern against an area of 5 cm<sup>2</sup> (2 cm x 2.5 cm)/ 0.77 in<sup>2</sup> (0.38 in x 2 in)
2. Put the tip of the swab back into the Extraction Solution. Cut the swab so that it fits in the tube with the cap closed and close the tube.
3. Shake the tube vigorously (manually or using a vortex mixer) for at least 2 minutes.
4. Filter the sample using one of the provided filter papers into one collection tube with yellow cap. See Figure 1 for details.

#### 8.4.2. Rinse waters

1. Take a sample directly from the CIP (Cleaning In Place) water circuit in a clean container (not provided).
2. Take 2 mL of the sample, add them to one of the tubes with Extraction Solution and close it.
3. Shake the tube vigorously (manually or using a vortex mixer) for at least 2 minutes.
4. Filter the sample using one of the provided filter papers into one collection tube with yellow cap. See Figure 1 for details.

## 8.5. Test procedure

- 8.5.1.** Make sure the AlerTox Heater has reached thermal stability at 110 °C (230 °F). Open the tube containing the AlerTox Lactose separation plates and take one. Close the tube immediately.
- 8.5.2.** Activation of the AlerTox Lactose separation plate. Put the AlerTox Lactose separation plate on the AlerTox Heater and wait 6 minutes. (*This time can be used to prepare the sample*). After 6 minutes, take the plate from the AlerTox Heater and lay it flat on a clean surface.
- 8.5.3.** Take 5 µL of filtered sample with a micropipette and apply the pipette tip on the sample application point (**S**) (left side of the plate, see Fig. 1). Let the sample flow out of the tip by capillarity and, if needed, press the plunger slightly to slowly release the volume. **Make sure the whole volume has been dispensed.**
- 8.5.4.** Dispense 5 µL of the Lactose Standard on its application point (**LS**) (right side of the plate, see Fig. 1) following the instructions in 8.5.3.
- 8.5.5.** Drying of the AlerTox Lactose Separation plate. Put the AlerTox Lactose Separation plate on the AlerTox Heater and dry it for 45 seconds. Take the AlerTox Lactose separation plate from the AlerTox Heater and lay it flat on a clean surface.




Figure 1. Sample preparation and AlerTox Lactose separation plate loading


- 8.5.6.** Insert the AlerTox Lactose separation plate into the Separation Chamber with Separation Solution so that the lower region of the plate (where the Sample and the Lactose Standard have been dispensed) meets the Solution. The plate must be inserted as straight as possible, ensuring that its lower edge reaches the Solution evenly. Close the chamber and allow the Separation Solution to rise to the upper line by capillarity (approximately 25 minutes).
- 8.5.7.** During this time, add the content of one of the tubes with Visualization Solution A to one of the supplied sprayers. Add 0.1 mL of Visualization Solution B to the same sprayer and close the sprayer and mix by inversion.
- 8.5.8.** Once the Separation Solution has reached the upper line of the AlerTox Lactose separation plate, take the plate out of the chamber and dry on the AlerTox Heater for 6 minutes. Take the plate from the AlerTox Heater and lay it flat on top of a paper towel or any other clean, absorbent material.

**Figure 2. Test procedure**

- 8.5.9.** At this point, use a mask to avoid inhalation. Spray the AlerTox Lactose separation plate with the solution prepared in 8.5.7. until it is **completely and uniformly soaked**.
- 8.5.10.** Immediately place the AlerTox Lactose separation plate on the AlerTox Heater for 10 minutes to develop it. The plate will turn pink.
- 8.5.11.** After the 10 minute development period, remove the AlerTox Lactose plate from the AlerTox Heater and read the result (section 8.5.). Switch off the AlerTox Heater. When the surface has cooled down, wipe it out with a damp cloth.

## 9. Reading and interpretation of the results

**Figure 3. Interpretation of results**

**Positive:** After the chromatography, a horizontal band appears in the sample (S) lane, at the same height as the band corresponding to the Lactose Standard (LS). Other bands/spots that may appear will correspond to other carbohydrates from the analyzed sample. It is important to differentiate between the migration smear and the actual band.

**Negative:** After the chromatography, NO band appears in the sample (S) lane at the same height as the band corresponding to the Lactose Standard (LS). Other bands/spots that may appear will correspond to other carbohydrates from the analyzed sample. It is important to differentiate between the migration smear and the actual band.

**Not valid:** The test is considered invalid if there is a total absence of bands both in the sample (S) and in the Lactose Standard (LS) lanes, an absence of the band corresponding to the Lactose Standard (LS) only, and/or an absence of change of coloration of the plate (from white to pink).

## 10. Matrix validation

AlerTox Lactose can be used to detect lactose in foods and beverages with a low content (less than 8%) of carbohydrates, rinse waters and on surfaces. For samples with a higher carbohydrate content and/or hydrolyzed foods, please contact your supplier.

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**Hygiena Mexico, S.A. de C.V.**