



EnSURE® and SystemSURE Plus® Monitoring System Installation & Operational Qualification

Site:	Instrument Serial Number:
-------	---------------------------

Purpose

This Installation/Operational Qualification has been specifically designed for the qualification of Hygiena’s EnSURE® and SystemSURE Plus® Monitoring Systems.

The Installation Qualification (IQ) and Operational Qualification (OQ) sections verify and document that the EnSURE or SystemSURE Plus Monitoring System has been installed and operates in accordance with the design specifications and manufacturer's recommendations. These verifications will be documented in the attached test report forms. Any discrepancies will be documented and reconciled. A fully trained Hygiena® representative will complete all forms, attachments and documentation at the time of installation.

Hygiena recommends that this document is fully reviewed and understood prior to initializing any of the testing contained therein.

Upon completion of each IQ or OQ test, the documentation should be signed, dated and witnessed. The following signatures are found throughout this document:

Signature	Printed Name	Company/Function



Contents

Installation Qualification

- IQ-01 Document Verification**
- IQ-02 Initial Inspection and Equipment List Verification**
- IQ-03 Operating Environment Requirements**
- IQ-04 Luminometer Setup and Registration**
 - **Installation Qualification Checklist**

Operational Qualification

- OQ-01 Luminometer Startup Verification**
- OQ-02 Calibration Verification**
- OQ-03 Test Device Verification**
- OQ-04 Database Verification**
 - **Operational Qualification Checklist**

Appendix A: Discrepancy Report Form



Installation Qualification

IQ-01 Document Verification

1. Objective

Verify that all the required documentation for the luminometer is available, complete and stored in the appropriate location.

Verify that the instrument calibration documentation is available and valid.

2. Test Procedure

Verify that the Operator Manual and Certificate of Conformity were received.

Record instrument serial number from Certificate of Conformity.

Attach Certificate of Conformity to this form.

3. Instrumentation Documentation

Operator Manual ([link](#)):

Yes

No

4. Calibration Documentation

Instrument Serial Number on Certificate:

Certificate of Conformity:

Yes

No

Certificate of Conformity attached:

Yes

No

5. Comments

6. Test Certification

Documented by: _____ Date: _____

Verified by: _____ Date: _____



IQ-02 Initial Inspection and Equipment List Verification

1. Objective

Verify and document that the correct equipment was received and is undamaged.

2. Test Procedure

Record the luminometer serial number.

Verify that all the equipment ordered was received.

Visually inspect all equipment and verify that it was received undamaged.

3. Equipment Record

Serial Number on Instrument:

- All items on the packing list are present: Yes No
- Packing list (or copy) is attached: Yes No
- Equipment received was the equipment ordered: Yes No
- Equipment received is visually undamaged: Yes No

4. Acceptance Criteria

All documents are present, and all equipment specified on the packing list is present.

5. Test Outcome

Meets criteria for equipment received and initial inspection: Yes No

If "No," refer to Discrepancy Number _____.

6. Comments

7. Test Certification

Documented by: _____ Date: _____

Verified by: _____ Date: _____



IQ-03 Operating Environment Requirements

1. Objective

Review the system environment requirements.

2. Test Procedure

Review the environmental requirements with the user.

3. Environment Requirements

Parameter	Specified
Operating Temperature	5 °C to 45 °C
Relative Humidity Range	20% to 85% non-condensing

4. Acceptance Criteria

The instrument is installed correctly and meets the required environmental conditions.

5. Test Outcome

Reviewed Environment Requirements: Yes No

If "No", refer to Discrepancy Number _____.

1. Comments

2. Test Certification

Documented by: _____ Date: _____

Verified by: _____ Date: _____



IQ-04 Setup and Registration

1. Objective

Set up the EnSURE or SystemSURE *Plus* instrument and SureTrend® software.

2. Test Procedure

Follow the instructions in the *Operator Manual* to complete the setup and instrument registration.

3. Acceptance Criteria

The instrument is set up correctly.

4. Test Outcome

Setup is complete: Yes No

If "No," refer to Discrepancy Number _____.

5. Comments

6. Test Certification

Documented by: _____ Date: _____

Verified by: _____ Date: _____



Installation Qualification Checklist

IQ-01 Document Verification: Passed Failed

IQ-02 Initial Inspection and Equipment List Verification: Passed Failed

IQ-03 Operating Environment Requirements: Passed Failed

IQ-04 Setup and Registration: Passed Failed

Installation Qualification passed? Yes No

If "Yes," then proceed to *Operational Qualification*.

Comments:

Documented by: _____ Date: _____

Verified by: _____ Date: _____



Operational Qualification

OQ-01 Luminometer Start-Up Verification

1. Objective

Verify the operation of the luminometer.

2. Materials and Reagents

Test Kit

Name: _____

Cat. No.: _____ Lot No.: _____ Exp.: _____

3. Test Procedure

Detailed instructions for the use and storage of Hygiena test devices are found in their respective kit inserts and on the website. Follow the instructions on the kit insert to activate the test and measure the results using the luminometer.

Results Screen appeared: Yes No

4. Acceptance Criteria

The results must appear on the instrument screen. If results appear out of range, please refer to *OQ-03 Calibration Verification*.

5. Test Outcome

Meets acceptance criteria for luminometer startup: Yes No
If "No," refer to Discrepancy Number _____.

6. Comments

7. Test Certification

Documented by: _____ Date: _____

Verified by: _____ Date: _____



OQ-02 Calibration Verification

1. Objective

Establish the Instrument Blank RLU value and verify that this value is within the acceptable range.

2. Materials and Reagents

CalCheck LED Calibration Verification Device (Cat No. CAL)
EnSURE or SystemSURE Plus Luminometer

3. Test Procedure

Detailed instructions for the use of Hygiena calibration verification devices are found in their respective kit inserts and on the website. Follow the instructions on the kit insert to perform the calibration verification.

4. Records

Record the results for calibration verification (as applicable) in the spaces provided below.

CalCheck (Cat No. CAL)

Calibration Verification test 1: _____, test 2: _____, test 3: _____.

5. Acceptance Criteria

The Calibration Verification must fall within the accepted range as outlined in the respective kit insert.

6. Test Outcome

Meets acceptance criteria for Calibration Verification: Yes No

7. Comments

8. Test Certification

Documented by: _____ Date: _____

Verified by: _____ Date: _____



OQ-03 Test Device Verification

1. Objective

Establish the competency of the user to test a surface using the luminometer and test device.

2. Materials and Reagents

UltraSnap® or SuperSnap® device
EnSURE or SystemSURE Plus luminometer

3. Test Procedure

Follow the instructions on [Form 1031, In-Service/Competency Checklist](#) and attach when complete.

4. Records

Attach *In-Service Competency Checklist*.

5. Acceptance Criteria

The Checklist is attached and all boxes are marked "Yes."

6. Test Outcome

Meets acceptance criteria for Test Device Verification: Yes No

7. Comments

8. Test Certification

Documented by: _____ Date: _____

Verified by: _____ Date: _____



OQ-04 Database Verification

1. Objective

Verify that data generated by the luminometer is transferred to the SureTrend database.

2. Test Procedure

Open SureTrend software.

Open the Sync application and press 'Sync Now' on the luminometer to upload data.

Verify that the data generated for OQ-3, OQ-4 was accurately transferred to the database.

Print the results from the database and attach them to this form.

Data transferred: Yes No

Results attached: Yes No

3. Acceptance Criteria

The data generated in OQ-3, OQ-4 must appear in the database.

4. Test Outcome

Meets acceptance criteria for Database Verification: Yes No

If "No," refer to Discrepancy Number _____.

5. Comments

6. Test Certification

Documented by: _____ Date: _____

Verified by: _____ Date: _____



Operational Qualification Checklist

- | | | |
|---|---------------------------------|---------------------------------|
| OQ-01 Luminometer Startup Verification: | <input type="checkbox"/> Passed | <input type="checkbox"/> Failed |
| OQ-02 Calibration Verification: | <input type="checkbox"/> Passed | <input type="checkbox"/> Failed |
| OQ-03 Test Device Verification: | <input type="checkbox"/> Passed | <input type="checkbox"/> Failed |
| OQ-04 Database Verification: | <input type="checkbox"/> Passed | <input type="checkbox"/> Failed |
| Operational Qualification passed? | <input type="checkbox"/> Yes | <input type="checkbox"/> No |

Comments:

Documented by: _____ Date: _____

Verified by: _____ Date: _____

End of Document



Appendix A: Discrepancy Report Form

Test ID: _____

Discrepancy Number: _____

Discrepancy	
Discrepancy Workout	<input type="checkbox"/> Permanently resolved. See Resolution below. <input type="checkbox"/> Temporarily resolved. See Resolution and Action Plan below.
Resolution	
Action Plan	

Done By/Date:	
Verified By/Date:	