



One Health Diagnostics™

Authors: Rumeysa Göcen, Anne Rölling, Patrice Chablain
Hygiena® Diagnostics GmbH, Hermannswerder 15, 14473 Potsdam, Germany

30-Minute Screening of Microorganisms for Plant- and Milk-Based Dessert Products with the Hygiena® Innovate™ System

Innovate™



Innovate™
AUTOSAMPLER III



INTRODUCTION:

The conversion of adenosine triphosphate (ATP) from the metabolic processes of living cells to adenosine diphosphate (ADP) can be measured using light detection systems. The Innovate™ Rapid Microbial Screening System is based on this bioluminescence-based method for fast and reliable detection of microbial contaminations. After product incubation for the enrichment of microorganisms in the package, it requires only 30 minutes with five minutes of hands-on time to decide if products can be released.

In this study, six different dessert products were inoculated with low-level spiking, and measured after 24, 48 and 72 hours of product incubation at 30 °C with the Innovate System and plated in parallel onto agar. The values, measured in RLU (Relative Light Units), show a clear difference between uncontaminated and contaminated products and demonstrate that the Innovate System is an alternative sterility control option not only for liquid products, but also for viscous dessert products.

In comparison to the plate count method, which needs an additional incubation step for at least two days, this method can detect contamination in dessert products within 48 to 72 hours, regardless of whether it is an ESL, UHT, or a dairy or plant-based dessert.

PURPOSE:

The objective of this study was to show that the 30-minute Innovate-based sterility test for milk- and plant-based desserts, produces qualitative results comparable to the plate count-based reference method after 72 hours incubation.

REGISTERED TRADEMARKS / GLOBAL CERTIFICATIONS:

Hygiena® is a registered trademark of Hygiena. Innovate™ and RapiScreen™ are trademarks of Hygiena. The RapiScreen Dairy Kit is AOAC certified #092301.

RapiScreen™ Dairy Kit



METHOD:

Microbial cultures from six gram-positive and gram-negative bacteria, spore-formers and yeasts were enriched in corresponding enrichment broth. Defined inoculum concentrations from 1 to 20 CFUs/Sample were spiked directly into six different dessert product packs. Spiked and non-inoculated products were incubated up to three days at 30 °C. After product incubation, aliquots of the product were taken for parallel testing using the Innovate System and the plate counting method. Due to the viscosity of the samples, 33.5 µL of α-amylase was added to 3 g of sample and shaken at 110 rpm for 10 minutes to simplify pipetting. Tryptic soy agar, Sabouraud dextrose agar and Columbia Agar with Sheep Blood were used for plate-count testing. Innovate System measurements were conducted according to the RapiScreen (RS) Dairy Kit. The standard RS dairy protocol was used for the Innovate System measurement.

In order to distinguish contaminated samples from non-contaminated samples, a threshold was set for each product: non-inoculated products were incubated at 30 °C for 48 hours as a negative control. After 48 hours of incubation, all samples were measured using the Innovate System. Four columns (32 wells) were measured to establish the baseline. The threshold for each product was set according to the baseline multiplied by three. All samples with an RLU value higher than the threshold value were considered contaminated.

Table 1: Enrichment/ Measurement Steps

Enrichment	Inoculum Concentration	1 - 20 CFUs/Sample
Time		24 - 72 hours
Temperature:		30 ± 2 °C
Kit		RapiScreen™ Dairy Kit
Measurement with Innovate		
Sample Volume		50 µL
Protocol		Standard dairy protocol
Time		30 min
Confirmation		
Method		According to ISO 4833-2:2013 Microbiology of the food chain - Horizontal method for the enumeration of microorganisms - Part 2: Colony count at 30 °C by the surface plating technique
Sample Volume		100 µL
Time		24 - 72 hours

Table 2: Strains

	Milk-Based Products	Plant-Based Products
<i>Bacillus cereus</i>	-	+
<i>Bacillus subtilis</i>	+	-
<i>Saccharomyces cerevisiae</i>	+	+
<i>Escherichia coli</i>	+	+
<i>Pseudomonas aeruginosa</i>	+	+
<i>Streptococcus thermophilus</i>	+	+
<i>Clostridium sporogenes</i>	+	+

RESULTS:

Table 3-5: Results of milk-based dessert products

[pos = contaminated samples, neg= contamination free samples, pos in orange = CFU/mL < 10⁴, pos in red = CFU/mL > 10⁴]

Cream Chocolate								
Threshold: 30 RLU								
Spore/Organism	Inoculation Level (CFUs/sample)	Incubation Temperature [°C]	24 hrs		48 hrs		72 hrs	
			Innovate System Mean RLU	ISO (pos/neg)	Innovate System Mean RLU	ISO (pos/neg)	Innovate System Mean RLU	ISO (pos/neg)
<i>Bacillus subtilis</i>	11	30 °C	4,796 7,300	pos pos	23,544 19,828	pos pos	14,641 10,857	pos pos
<i>Saccharomyces cerevisiae</i>	5.6	30 °C	15 17	pos pos	35,920 38,305	pos pos	burst	
<i>Escherichia coli</i>	15	30 °C	2,072 54	pos pos	3,864 2,720	pos pos	1,331 775	pos pos
<i>Pseudomonas aeruginosa</i>	7	30 °C	22 10	pos pos	10,000 10,225	pos pos	9,373 10,593	pos pos
<i>Streptococcus thermophilus</i>	20	30 °C	83 11	pos pos	20,604 21,957	pos pos	7,895 7,294	pos pos
<i>Clostridium sporogenes</i>	12	30 °C	8 7	neg neg	8 9	neg neg	10 7	neg neg

Cream Vanilla								
Threshold: 42 RLU								
Spore/Organism	Inoculation Level (CFUs/sample)	Incubation Temperature [°C]	24 hrs		48 hrs		72 hrs	
			Innovate System Mean RLU	ISO (pos/neg)	Innovate System Mean RLU	ISO (pos/neg)	Innovate System Mean RLU	ISO (pos/neg)
<i>Bacillus subtilis</i>	11	30 °C	43,658 43,209	pos pos	80,056 83,225	pos pos	76,155 61,573	pos pos
<i>Saccharomyces cerevisiae</i>	5.6	30 °C	62 56	pos pos	435,069 384,967	pos pos	burst	
<i>Escherichia coli</i>	15	30 °C	10,884 12,452	pos pos	57,533 62,761	pos pos	21,306 24,326	pos pos
<i>Pseudomonas aeruginosa</i>	7	30 °C	557 399	pos pos	44,795 55,240	pos pos	39,106 46,737	pos pos
<i>Streptococcus thermophilus</i>	20	30 °C	120 79	pos pos	73,959 78,567	pos pos	55,454 54,786	pos pos
<i>Clostridium sporogenes</i>	12	30 °C	11 12	neg neg	13 12	neg neg	11 10	neg neg

Mousse Chocolate								
Threshold: 36 RLU								
Spore/Organism	Inoculation Level (CFUs/sample)	Incubation Temperature [°C]	24 hrs		48 hrs		72 hrs	
			Innovate System Mean RLU	ISO (pos/neg)	Innovate System Mean RLU	ISO (pos/neg)	Innovate System Mean RLU	ISO (pos/neg)
<i>Bacillus subtilis</i>	11	30 °C	244 428	pos pos	17,582 13,695	pos pos	2,708 437	pos pos
<i>Saccharomyces cerevisiae</i>	5.6	30 °C	13 12	pos pos	12,747 13,751	pos pos	Burst	
<i>Escherichia coli</i>	15	30 °C	6 18	pos pos	970 944	pos pos	1,920 1,588	pos pos
<i>Pseudomonas aeruginosa</i>	7	30 °C	12 12	pos pos	1,214 1,042	pos pos	1,337 1,645	pos pos
<i>Streptococcus thermophilus</i>	20	30 °C	77 20	pos pos	3,422 4,104	pos pos	1,716 2,303	pos pos
<i>Clostridium sporogenes</i>	12	30 °C	5 8	neg neg	8 11	neg neg	6 6	neg neg

Table 6-8: Results of milk-based protein-rich dessert products

[pos = contaminated samples, neg= contamination free samples, pos in orange = CFU/mL < 10⁴, pos in red = CFU/mL > 10⁴]

Protein Pudding Hazelnut chocolate								
Threshold: 33 RLU								
Spore/Organism	Inoculation Level (CFUs/sample)	Incubation Temperature [°C]	24 hrs		48 hrs		72 hrs	
			Innovate System Mean RLU	ISO (pos/neg)	Innovate System Mean RLU	ISO (pos/neg)	Innovate System Mean RLU	ISO (pos/neg)
<i>Bacillus subtilis</i>	11	30 °C	6,718 5,400	pos pos	18,293 21,950	pos pos	40,381 43,586	pos pos
<i>Saccharomyces cerevisiae</i>	5.6	30 °C	10 17	pos pos	173 196	pos pos	63 73	pos pos
<i>Escherichia coli</i>	15	30 °C	17 3,977	pos pos	129 66	pos pos	48 45	pos pos
<i>Pseudomonas aeruginosa</i>	7	30 °C	116 425	pos pos	2,325 1,125	pos pos	2,579 22,021	pos pos
<i>Streptococcus thermophilus</i>	20	30 °C	529 257	pos pos	10,244 90,745	pos pos	2,199 2,065	pos pos
<i>Clostridium sporogenes</i>	12	30 °C	6 8	neg neg	7 8	neg neg	7 8	neg neg

Protein Pudding Vanille								
Threshold: 51 RLU								
Spore/Organism	Inoculation Level (CFUs/sample)	Incubation Temperature [°C]	24 hrs		48 hrs		72 hrs	
			Innovate System Mean RLU	ISO (pos/neg)	Innovate System Mean RLU	ISO (pos/neg)	Innovate System Mean RLU	ISO (pos/neg)
<i>Bacillus subtilis</i>	11	30 °C	18,248 12,053	pos pos	9,381 8,029	pos pos	14,302 12,150	pos pos
<i>Saccharomyces cerevisiae</i>	5.6	30 °C	17 21	pos pos	76 24	pos pos	1,243 78	pos pos
<i>Escherichia coli</i>	15	30 °C	444 353	pos pos	1,268 4,680	pos pos	8,764 180	pos pos
<i>Pseudomonas aeruginosa</i>	7	30 °C	1,107 117	pos pos	4,404 6,184	pos pos	10,193 9723	pos pos
<i>Streptococcus thermophilus</i>	20	30 °C	135 138	pos pos	32,967 32,212	pos pos	24,797 27,432	pos pos
<i>Clostridium sporogenes</i>	12	30 °C	13 17	neg neg	13 6	neg neg	12 11	neg neg

Table 8: Results of the plant-based dessert products

[pos = contaminated samples, neg= contamination free samples, pos in orange = CFU/mL < 10⁴, pos in red = CFU/mL > 10⁴]

Dark Chocolate (plant-based)								
Threshold: 30 RLU								
Spore/Organism	Inoculation Level (CFUs/sample)	Incubation Temperature [°C]	24 hrs		48 hrs		72 hrs	
			Innovate System Mean RLU	ISO (pos/neg)	Innovate System Mean RLU	ISO (pos/neg)	Innovate System Mean RLU	ISO (pos/neg)
<i>Bacillus cereus</i>	6.6	30 °C	3 7	pos neg	16,326 4	pos neg	11,006 4	pos neg
<i>Saccharomyces cerevisiae</i>	8	30 °C	3 11	pos pos	13 10	pos pos	7,130 195	pos pos
<i>Escherichia coli</i>	8.2	30 °C	439 916	pos pos	1,870 1,693	pos pos	1,537 1,794	pos pos
<i>Pseudomonas aeruginosa</i>	3.9	30 °C	12 11	pos pos	9 9	pos pos	323 433	pos pos
<i>Streptococcus thermophilus</i>	6.5	30 °C	8 9	neg neg	3,005 127	pos pos	4,356 523	pos pos
<i>Clostridium sporogenes</i>	12	30 °C	8 8	neg neg	838 880	pos pos	624 955	pos pos

After 48 hours of static incubation of the inoculated dessert product packs at 30 °C, all products were able to be analyzed correctly. The Innovate System results corresponded 100% to the plate-count measurements, which required a second plate incubation after enrichment in the product.

Clostridium sporogenes was not detected in milk-based products either with the Innovate System nor by plate growth.

For better pipetting, measurements were also taken with and without an amylase treatment. The amylase had no effect on the measurements, but it led to a comfortable handling of the samples. (Data not shown in the tables).

The 100% agreement of the results after 48 hours of incubation between ISO and the alternative Innovate System method were independent of plant-based or milk-based production, extended-shelf-life (ESL) or ultra-heat-treated (UHT) products.



SIGNIFICANCE:

ATP detection with the Innovate System RapiScreen Dairy system is a rapid, accurate and safe 30-minute testing for processed dessert products. This demonstrates that the Innovate system is suitable for commercial sterility testing of challenging dessert matrices, even chocolate cream products.

Without a time-consuming second plate incubation, the Innovate System generates comparable results to the reference method and shortens product release by more than 3 days, especially for ESL dessert products.