



Development of a Multiplex Real-time PCR Assay for the Detection of Highly Pathogenic *Salmonella enterica* (HPS) in Beef, Poultry and Pork

Hans-Henno Dörries¹, Dayna M. Harhay², Kerry D. Brader², Monali Gandhi³, Ivonne Remus-Dörries¹, Cordt Grönwald¹, Patrice Chablain¹, Rebecca L. Olsen³, and Stacy K. Stollenberg³.
¹Hygiena[®] Diagnostics GmbH, Hermannswälder 17, Potsdam, Germany 17743
²USDA-ARS Meat Animal Research Center, Clay Center, Nebraska, USA 68933
³Hygiena[®], 941 Avenida Acea Carrillo, California, USA 93012

BAX[®] System Q7

BAX[®] System X5

foodproof[®]

microproof[®]

INTRODUCTION:

Salmonella remains a significant health concern causing millions of cases of foodborne illness annually. The CDC estimates *Salmonella* bacteria causes 1.35 million infections, more than 26,000 hospitalizations, and over 400 deaths in the United States each year. Food is the source of a majority of these illnesses.

While there are established performance standards set by the USDA to reduce the prevalence of *Salmonella* in poultry, there has been no significant reduction in the number of human infections from *Salmonella* in poultry products.

The USDA and the food industry work to find a balance to continue to offer wholesome, nutritious meat products while reducing the risk of illness from *Salmonella* for humans. The established performance standards measure the prevalence of *Salmonella*, but the risk of human illness is associated with the amount of *Salmonella* in a product and its virulence.

There is a need in the poultry and beef industry for a rapid method to identify the risk of *Salmonella* in food samples based on virulence or pathogenicity, instead of prevalence alone.

PURPOSE:

1. Development of a real-time PCR assay to detect *Salmonella* in relevant food samples based on genes associated with pathogenicity as researched by USDA.
2. Evaluate the performance of the assay in detecting the pathogenicity of samples previously characterized by USDA as highly pathogenic *Salmonella* (HPS).
3. Demonstrate the inclusivity and exclusivity of the assay.

REGISTERED TRADEMARKS

BAX[®] is a registered trademark of Hygiena[®]
 naica[®] is a registered trademark of Stilla Technologies

METHODS:

The USDA Agriculture Research Service (USDA-ARS) conducted a study to characterize differences in virulence gene content among *Salmonella* serovars most commonly linked to outbreaks and human illness, and serovars more commonly found as commensals associated with food animals. Primers and probes for the PCR assay were selected based on the virulence genes and targets researched by USDA-ARS and shared via a CRADA. A gel-based PCR assay was developed by USDA-ARS, and a real-time lyophilized PCR assay was developed by Hygiena[®].

The lyophilized PCR assay was developed to detect three HPS targets and *Salmonella* spp. In order to check the results regarding PCR inhibition, an internal control was added as a fifth target. The assay was tested on FSIS *Salmonella*-positive enrichments and isolates prepared by USDA from each of the four major commodities: chicken, turkey, pork, and beef.


In this study, 30 turkey samples, 36 chicken samples, 33 pork samples, and 30 beef samples were evaluated using the newly developed real-time PCR assay.

Pure, RNA-free DNA of *Salmonella* Enteritidis and *Salmonella* Dublin was quantified by droplet dPCR (dPCR; Stilla, naica[®] instrument) and used in a POD approach to estimate the absolute LOD_{95%}. DNA concentrations of 25, 5, 2.5, 0.25, and 0.025 copies/reaction were tested in 32 replicates each.

RESULTS:

Sensitivity (Absolute Limit of Detection (LOD_{95%}))

Target	LOD _{95%}
HPS-B	3.08
HPS-A	3.42
HPS-X	3.82
<i>Salmonella</i> spp.	3.66



S. Enteritidis	Number of positive vs total number of samples; DNA Concentration [copies/μL; (copies/reaction)] (CV [%])			
Target	1 (25)	0.2 (5)	0.1 (2.5)	0.01 (0.25)
HPS-B	16/16; 32.29 (0.41)	16/16; 34.59 (1.30)	13/16; 33.82 (1.04)	4/16; 36.43 (0.51)
HPS-A	16/16; 32.82 (0.23)	16/16; 34.96 (0.83)	15/16; 36.14 (1.34)	2/16; 36.96 (0.10)
HPS-X	16/16; 30.81 (0.27)	16/16; 32.78 (0.94)	13/16; 33.82 (1.04)	3/16; 36.12 (2.53)
<i>Salm. spp.</i>	16/16; 32.02 (0.42)	15/16; 34.76 (1.00)	15/16; 35.34 (0.78)	4/16; 36.37 (0.69)
IPC	16/16; 29.36 (0.34)	16/16; 29.22 (0.36)	16/16; 29.37 (0.29)	16/16; 29.51 (0.50)

S. Dublin	Number of positive vs total number of samples; DNA Concentration [copies/μL; (copies/reaction)] (CV [%])			
Target	1 (25)	0.2 (5)	0.1 (2.5)	0.01 (0.25)
HPS-B	16/16; 30.93 (0.47)	16/16; 32.85 (0.71)	16/16; 33.86 (0.89)	3/16; 35.77 (0.62)
HPS-A	16/16; 31.72 (0.30)	16/16; 33.45 (0.58)	14/16; 34.44 (0.67)	2/16; 34.93 (0.24)
HPS-X	16/16; 29.80 (0.69)	16/16; 32.16 (0.60)	14/16; 32.87 (0.65)	2/16; 32.77 (0.61)
<i>Salm. spp.</i>	16/16; 31.33 (0.62)	16/16; 33.70 (0.77)	13/16; 34.98 (0.86)	4/16; 35.95 (0.82)
IPC	16/16; 30.20 (0.83)	16/16; 30.09 (0.62)	16/16; 30.17 (0.75)	16/16; 30.11 (1.12)

COMPARISON OF USDA ARS GEL-BASED PCR WITH REAL-TIME ASSAY FOR ISOLATES:

	HPS-B	HPS-A	HPS-X	<i>Salmonella</i> spp.
Agreement	30	30	28	30
Deviation	0	0	2	0
Accuracy [%]	100.0	100.0	93.3	100.0

	HPS-B	HPS-A	HPS-X	<i>Salmonella</i> spp.
Agreement	30	29	30	30
Deviation	0	1	0	0
Accuracy [%]	100.0	96.7	100.0	100.0

	HPS-B	HPS-A	HPS-X	<i>Salmonella</i> spp.
Agreement	35	36	36	36
Deviation	1	0	0	0
Accuracy [%]	97.2	100.0	100.0	100.0

	HPS-B	HPS-A	HPS-X	<i>Salmonella</i> spp.
Agreement	33	33	33	33
Deviation	0	0	0	0
Accuracy [%]	100.0	100.0	100.0	100.0

RESULTS:

INCLUSIVITY & EXCLUSIVITY

Specificity studies yielded 100% inclusivity on 204 *Salmonella* serovars (Table 5) and 100% exclusivity on 53 closely related species of *Enterobacteriaceae*.

The data from the specificity study with serovars from the in-house culture collection also show 100% agreement with the USDA-ARS gel-based PCR.

Table 5: Comparison of gel-based PCR vs. real-time PCR assay for specificity study

Serotype	Number tested	Original primer sets (% positive)			HPS real-time PCR kit (% positive)			Serotype	Number tested	Original primer sets (% positive)			HPS real-time PCR kit (% positive)		
		HPS-B	HPS-A	HPS-X	HPS-B	HPS-A	HPS-X			HPS-B	HPS-A	HPS-X	S. spp.		
4,12:12i:-	3	100	100	100	100	100	100	Manhattan	3	0	0	0	0	0	0
Abaketeuba	1	0	0	0	0	0	0	Manila	1	0	100	0	0	100	0
Aberdeen	1	0	100	100	0	100	100	Meleagridis	1	0	100	0	0	100	0
Abony	2	0	100	100	0	100	100	Mississippi	1	100	100	100	100	100	100
Agona	1	0	100	0	0	100	0	Montevideo	7	0	14	0	0	14	0
Alatum	3	0	0	0	0	0	0	Munich	1	0	0	100	0	0	100
Bardo	3	100	100	100	100	100	100	Munster	2	0	0	0	0	0	0
Bareilly	2	0	0	100	0	0	100	Napoli	1	0	0	100	0	0	100
Blockley	4	100	100	100	100	100	100	Newport	6	0	100	100	0	100	100
Bransfordians	2	0	100	100	0	100	100	Oxensburg	1	0	0	0	0	0	0
Brandenburg	3	0	0	0	0	0	0	Paratyphi	4	0	0	0	0	0	0
Brenderup	1	100	100	100	100	100	100	Paratyphi B	9	11	22	33	11	22	33
Cerro	2	0	0	0	0	0	0	Paratyphi C	1	100	100	0	100	100	0
Choleraesuis	1	100	100	100	100	100	100	Reading	3	0	33	33	0	33	33
Concord	1	0	100	100	0	100	100	Rissen	1	0	100	0	0	100	0
Cottbus	1	0	0	100	0	0	100	S. bongori	2	50	100	0	50	100	0
Crossness	1	0	100	100	0	100	100	S. enterica subsp. anthonomus	10	0	0	0	0	0	0
Cubana	1	0	100	0	0	0	0	S. enterica subsp. dissonans	8	33	0	0	33	0	0
Derby	3	0	33	0	0	33	0	S. enterica subsp. enterica	6	33	67	33	33	67	33
Dublin	8	88	100	88	88	100	88	S. enterica subsp. houtense	7	0	0	100	0	0	100
Enteritidis	3	100	100	100	100	100	100	S. enterica subsp. indica	2	50	0	0	50	0	0
Gallinarum	3	100	100	100	100	100	100	S. enterica subsp. salmoe	8	33	33	100	33	13	0
Gallinarum-Pullorum	3	100	100	100	100	100	100	Saintpaul	3	0	100	100	0	100	100
Hadar	5	60	100	100	60	100	100	Sentfihenberg	14	0	100	0	0	100	0
Havana	1	100	0	0	100	0	0	Stanley	1	0	100	0	0	100	0
Heidelberg	6	0	100	100	0	100	100	Stanleyville	1	0	0	0	0	100	0
Hvittingfoss	1	0	100	100	0	100	100	Tennessee	1	0	100	0	0	100	0
Infantis	11	0	100	100	0	100	100	Thompson	2	100	100	100	100	100	100
Johannesburg	1	100	0	0	100	0	0	Typilmurum	1	100	100	100	100	100	100
Kentucky	1	0	100	0	0	100	0	Uganda	1	0	0	100	0	0	100
Krefeld	1	0	100	0	0	100	0	Virchow	2	0	100	100	0	100	100
Lille	1	0	100	0	0	100	0	Waycross	1	0	0	0	0	0	0
Utchfield	1	100	100	100	100	100	100	Waltveden	2	100	0	100	100	0	100
Livingstone	4	0	50	25	0	50	25	Westhampton	1	0	100	0	0	100	0
London	1	0	100	100	0	100	100		204						

SIGNIFICANCE:

The results demonstrate that highly pathogenic *Salmonella* can be quickly detected in meat and poultry using this next-generation method for *Salmonella* detection. This approach to *Salmonella* testing can enable the industry to make risk-based decisions about their product by offering an indication of when product has more infectious potential. It would help achieve the objective of simultaneously reducing the number of illnesses caused by *Salmonella* while continuing to provide consumers with nutritious and affordable meat.

REFERENCES:

1. <https://www.cdc.gov/salmonella/index.html>
2. Response to Questions Posted by the Food Safety and Inspection Service: Enhancing *Salmonella* Control in Poultry Products, Journal of Food Protection, 87, 2024
https://www.fsis.usda.gov/sites/default/files/media_file/documents/NACMCF_Salmonella_2023.pdf

The research presented was conducted by scientists with the USDA-ARS US Meat Animal Research Center (USMARC). A Cooperative Research and Development Agreement (CRADA No. 58-3040-1-009) is established between USDA and Hygiena[®] LLC for the development of the HPS markers into a real-time assay intended for the BAX[®] System.

Disclaimer

Names are necessary to report factually on available data; however, the USDA neither guarantees nor warrants the standard of the product, and the use of the name by USDA implies no approval of the product to the exclusion of others that may also be suitable. USDA is an equal opportunity provider and employer.