

InSite® | *Salmonella*

Hygiena's InSite® *Salmonella* is a quick and easy test for the detection of *Salmonella* species in food processing environments. The test contains a liquid medium formulated with growth enhancers and chromogenic compounds selective for *Salmonella* species.



Benefits:

- Easy to use with no preparation steps required
- Large sponge tip maximizes sample collection for accurate results
- Eliminates costly outside lab fees for an equivalent test
- Color results are easy to interpret—bright yellow indicates presence of *Salmonella* species
- Presumptive positive results in as little as 24 hours
- Negative results in 48 hours
- Self-contained tube eliminates risk of contamination, makes testing easier, faster, and more cost-effective

Sensitivity:

- Detects *Salmonella* species down to <10 CFU, even in the presence of high numbers of competing organisms

Interpretation:

- InSite *Salmonella* features a viewing window with a result chart for easy interpretation. Results are presumptive positive when the Insite *Salmonella* media turns bright yellow.




Results in the Viewing Window of the InSite Device
(This diagram shows the label laid out flat.)

Results	Medium (viewing window)	Swab Tip
Positive	Yellow	Yellow
Negative	Brown	Yellow
Negative	Purple	Yellow
Negative	Purple	Purple

Interpretation of Color Results

Catalog No.	Description	Quantity
IS050	InSite <i>Salmonella</i> Environmental <i>Salmonella</i> Species Test	50

InSite® Salmonella Procedure



See an instructional demo at [HygienaTV](#) or [HygienaTV](#)

Related Products

Digital Dry Block Incubators

Digital incubators are capable of up to 105 °C and flexible for any Hygiena test device. Available in two sizes.




InSite® *Listeria* Environmental *Listeria* Species Test

InSite *Listeria* is a quick and easy test for the detection of *Listeria* species in 24 – 48 hours. The test contains a liquid media formulated with selective agents, growth enhancers, and color changing compounds specific to the growth of *Listeria* species.



MicroSnap® Coliform & *E. coli*

MicroSnap Coliform & *E. coli* are rapid tests for detection and enumeration of Coliform bacteria. Single figure organisms can be detected in the same day.



MicroSnap® Total

MicroSnap Total delivers total viable count (TVC) results in the same day.

MicroSnap® EB

MicroSnap EB delivers results in the same day.

BAX® System Instruments

The BAX System harnesses the power of DNA-based foodborne pathogen detection to help you make product decisions with confidence. An extensive assay portfolio, multiple platforms, and flexible protocols enable food labs of any size to take advantage of this benchmark technology.



Catalog No.	Description	Quantity
INCUBATOR/INCUBATOR2	Digital Dry Block Incubators	1
IL050/IL100	InSite <i>Listeria</i> Environmental <i>Listeria</i> Species Test	50/100
MS2-COLIFORM	MicroSnap Coliform	100
MS2-TOTAL	MicroSnap Total	100
MS2-EB	MicroSnap EB	100
MS2-ECOLI	MicroSnap <i>E. coli</i>	100
MCH-2007	BAX System X5	1
ASY2018	BAX System Q7	1