

Deja Latney<sup>1</sup>, Julie Weller<sup>1</sup>, Savannah Applegate<sup>1</sup>, Jerri Lynn Pickett<sup>2</sup> and **Jacquelyn Adams<sup>2</sup>** 1. Hygiena<sup>®</sup>, 2 Boulden Circle, New Castle, DE 19720, 2. Tyson Foods, 3609 Johnson Rd, Springdale, AR 72762

### **INTRODUCTION:**

Breaded, stuffed raw chicken products have been associated with 14 Salmonella outbreaks since 1998. These products appear to be fully cooked due to a crispy brown exterior, but proper cooking is required to eliminate potential microbial hazards. Since this type of product has repeatably been linked to illness, the U.S. Department of Agriculture's Food Safety and Inspection Service (USDA FSIS) proposed new regulatory action on April 25, 2023, to declare Salmonella an adulterant in breaded stuffed raw chicken products (1).

For many years, prevalence testing for Salmonella in poultry products has provided the industry with important data, but new diagnostic capabilities for quantification are needed to determine how much of the pathogen is present.

#### **PURPOSE:**

The objective of this study was to meet the new proposed standard of 1 CFU/g limit by:

- 1. Developing a linear equation for breaded stuffed raw chicken products
- 2. Verifying the BAX<sup>®</sup> System Real-Time PCR assay for Salmonella quantification (SalQuant<sup>®</sup>)

### **REGISTERED TRADEMARKS:**

BAX<sup>®</sup> is a registered trademark of Hygiena for its line of equipment, reagents and software used to analyze samples for microbial contamination. Hygiena® and SalQuant® are registered trademarks of Hygiena. Tyson<sup>®</sup> is a registered trademark of Tyson Foods, Inc.

### **METHODS:**

Four varieties of frozen, breaded, stuffed raw chicken products were thawed and divided into 375 g test portions. Primary homogenates were made by adding 375 mL of BPW. Then, 15 x 30 mL samples were removed for inoculation. A coldstressed culture of Salmonella was added at varying concentrations in replicates of 3 to reach levels of 1, 10, 100, 1,000 and 10,000 CFU/mL. One additional sample was reserved as a negative control.

Enrichment Following inoculation, samples were enriched in 30 mL of prewarmed (42 °C) BAX MP media with 0.5 mL/L Quant Solution, incubated for 6 hours and then tested in guintuplet using realtime PCR.

At the same time, a 3-tube x 5-dilution MPN was conducted for each inoculation level following the USDA FSIS Appendix 2.05.

# Salmonella Quantification (SalQuant<sup>®</sup>) with Hygiena's **BAX® System for Breaded Stuffed Raw Chicken Products**

# BAX<sup>®</sup> System 7

#### **Sample Preparation**

### **FIGURES:**



Figure 1 (Left). Mean (Salmonella Ct) and Salmonella Ct of Combined Data from All Four Products vs. Inoculated Log CFU/g.

Figure 2 (Right). MPN and SalQuant Comparison per Inoculation Level at 6 Hours of Enrichment of All Four Products, Combined.

### **RESULTS:**

Curve Development (Figure 1): At 6 hours, a curve was constructed for each product type, and then the combined mean was used to generate one equation with R<sup>2</sup> and Log Root Mean Square Error (RMSE).

# BAX<sup>®</sup> System X 5

# foodproof®

# microproof<sup>®</sup>

• R<sup>2</sup> of 0.86 • Log RMSE of 0.57 Enumerable range of 0.0 – 4.0 Log CFU/g

MPN Comparison (Figure 2): There were no statistical differences between MPN values and SalQuant estimations since the error bars for each inoculation level overlap.



# **SIGNIFICANCE:**

These results demonstrate that complete quantification of Salmonella from 0.0 to 4.0 Log CFU/mL can be achieved using the BAX System Real-Time PCR Assay for Salmonella.

## **REFERENCES:**

. USDA Press. April 25, 2023. USDA Proposes Declaring Salmonella an Adulterant in Breaded Stuffed Raw Chicken Products. Press Release No. 0090.23. https://www.usda.gov/media/pressreleases/2023/04/25/usda-proposes-declaring-salmonella-adulterantbreaded-stuffed-raw

