



Innovate System Software Update

Introducing Innovate System Windows 10 Software Release v5.25

Hygiena's Innovate System's v5.25 upgrade features Win10 compatibility, new and improved software features, and the latest user documentation to keep you and your Innovate System up and running.

Innovate System software v5.25 contains all new Hygiena branding. As of this release, Windows 10 Professional, 64-bit is now a fully supported platform for the Innovate System.

Because each release of Innovate Software builds on improvements from previous versions, it is strongly recommended that you always install the most recent version as soon as possible for optimal performance.

Software Updates for Release of v5.25

Functional Update	Old Innovate V5 Protocol	New Innovate V5.25 Protocol
Scan COM Ports	<ol style="list-style-type: none">1. Close Innovate V5 Software2. Navigate to C:\ProgramData\InnovateV5\Cell Scan\Service.exe3. Open service.exe application4. Select 'Scan COM Ports' icon5. Close Service.exe application6. Login to Innovate Software	<ol style="list-style-type: none">1. Login to Innovate Software2. Navigate to 'Instrument' tab3. Select 'Detect Instrument'
Single Installation File	<ol style="list-style-type: none">1. Install Innovate.im V5 software2. Install Innovate Driver3. Once install is completed, navigate to C:\ProgramData\InnovateV5 and rename 'Para' folder 'ParaBackup' folder4. Create new 'Para' folder and paste RS_Dairy, RS_Dairy_Control, and RS_Injector protocols into new Para folder5. Navigate to desktop Innovate software icon and right click6. Select 'Properties'7. Delete 'RS_Dairy_Square.par and replace with RS_Dairy.par in the target location	<ol style="list-style-type: none">1. Install Innovate software and driver with single installation file



Supported Languages	1. English 2. German	1. English 2. German 3. French 4. Spanish 5. Portuguese 6. Italian
Operating System Compatibility	1. Windows XP 2. Windows 7	1. Windows 10

v5.25 Software Installation Instructions

1. (If current Innovate User): Backup historical data and protocols by copying the InnovateV5 folder from C:\ProgramData\InnovateV5 and C:\Program files (x86)\InnovateV5\Database and pasting on a separate flash drive.
2. Right click on Innovate.exe application from a USB or computer download. Click **Run as Administrator** to ensure the software is installed properly.
3. Follow the on-screen prompts carefully. For a standard installation, choose default responses.
4. (If current Innovate user): Once the software is installed successfully, copy&replace the existing InnovateV5 folder with your saved backup folder. Do the same for the database folder.
5. (If current Innovate user): Login to Innovate software and confirm your settings are loaded successfully.