



Boost in-clinic
results with
Dermapen HOME™

Patients will be relying on you
to recommend the pen level and
speed that's right for them.

Dermapen HOME™ is hygienic and safe

Dermapen HOME™ treatments are:

- ✔ Convenient
- ✔ Fast
- ✔ Comfortable

Your patients can adjust the pen level and speed to deliver a customised treatment that suits their skin concern.



Dermapenworld
ORIGINAL TRUSTED
AUSTRALIAN



Dermapenworld™

GLUE AREA
(tear-off pad)

Name: _____ Today's Date: _____

Next Professional Dermapen Treatment Date: _____

Mark the skin concerns which apply



Ageing/sun-damaged skin



Pigmentation



Skin texture




Scarring



Fine lines and wrinkles



Rosacea



Problematic skin



Enlarged pores

Before microneedling

Start by sanitising your hands

Cleanse your face



TRI-PHASE CLEANSER™

OR



CLR FOAM CLEANSER™

Clear away dull or dead skin cells



MICRO DERM EXFOLIANT™

Spritz & remove impurities with gauze (or cotton pads)



CLINIPREP™

During microneedling

Apply Meso-Glide™ serums to the six areas shown overleaf

General skin concerns



HYLA ACTIVE™

Pigmentation



HYLA ACTIVE™ and BRITE LITE™

Problematic skin




HYLA ACTIVE™ and CLR LOTION™

Anti-ageing



HYLA ACTIVE™ and ANTIOXIDANT COCKTAIL™

Your **Dermapen**HOME™ settings



1. Forehead

Pen level 1 2 3 4
 Pen speed 1 2 3 4

2. Eye area

Pen level 1 2 3 4
 Pen speed 1 2 3 4

3. Nose

Pen level 1 2 3 4
 Pen speed 1 2 3 4

4. Cheeks

Pen level 1 2 3 4
 Pen speed 1 2 3 4

5. Upper lip and lips

Pen level 1 2 3 4
 Pen speed 1 2 3 4

6. Chin/jawline

Pen level 1 2 3 4
 Pen speed 1 2 3 4

Immediately after microneedling



**HYLA ACTIVE™ 3D
SCULPTURED MASK**
for 10-30 minutes



VITAMIN RICH REPAIR™



COVER RECOVER™

Your Dp Dermaceuticals™ morning and evening routine

Cleanse		<input type="checkbox"/> CLINIPREP™	<input type="checkbox"/> TRI-PHASE CLEANSER™	<input type="checkbox"/> CLR FOAM CLEANSER™	<input type="checkbox"/> MICRO DERM EXFOLIANT™
		<input type="checkbox"/> CLINIPREP™	<input type="checkbox"/> TRI-PHASE CLEANSER™	<input type="checkbox"/> CLR FOAM CLEANSER™	<input type="checkbox"/> MICRO DERM EXFOLIANT™
Nourish		<input type="checkbox"/> HYLA ACTIVE™	<input type="checkbox"/> ANTIOXIDANT COCKTAIL™	<input type="checkbox"/> EXO-SKIN™	
		<input type="checkbox"/> HYLA ACTIVE™	<input type="checkbox"/> ANTIOXIDANT COCKTAIL™	<input type="checkbox"/> EXO-SKIN™	
Correct		<input type="checkbox"/> VITAMIN RICH REPAIR™	<input type="checkbox"/> RETINAL ACTIVE™	<input type="checkbox"/> BRITE LITE™	<input type="checkbox"/> SSSCAR™
		<input type="checkbox"/> VITAMIN RICH REPAIR™	<input type="checkbox"/> RETINAL ACTIVE™	<input type="checkbox"/> BRITE LITE™	<input type="checkbox"/> SSSCAR™
Moisturise		<input type="checkbox"/> SKIN VENEER™	<input type="checkbox"/> CLR LOTION™		
		<input type="checkbox"/> SKIN VENEER™	<input type="checkbox"/> CLR LOTION™		
Protect		<input type="checkbox"/> COVER RECOVER™	<input type="checkbox"/> BARRIER BODY BUTTER™		
Boosters					

Our Synergy of Solutions

Best results are achieved when Dermapen HOME™ is used holistically. DermapenWorld Treatments™ combine:

- Dermapen 4™ in-clinic procedures
- MG-Collection™ professional serums
- ÜBER PRO & MD™ professional peels
- Dp Dermaceuticals™ EXOSOME™ therapy in-clinic and at-home
- Dp Dermaceuticals™ L.E.D.™ FACE mask for use in-clinic and at-home
- Dp Dermaceuticals™ for in-clinic and at-home procedures
- Weekly maintenance with Dermapen HOME™

Nothing routine



Dermapenworld™

Scan to learn more at
dermapenworld.com/dermapen-home

