BOOKWALTER® WISHBOOK Adjustable Frame

Symmetry Surgical Inc.
3034 Owen Drive,
Antioch, TN 37013 USA
☎ 1-800-251-3000
FAX: 1-615-964-5566
www.symmetrysurgical.com

EC REP

Symmetry Surgical GmbH
Maybachstraße 10
78532 Tuttlingen, Germany
☎ + 49 7461 96490
Fax: + 49 7461 77921

LCN 204046-001/H
© 2004–2015 Symmetry Surgical
♻️ Revised 09/15

ENGLISH

IMPORTANT INFORMATION

Please Read Before Use

BOOKWALTER® WISHBOOK Adjustable Frame

Description

The BOOKWALTER WISHBOOK Adjustable Frame, when used with other components of the BOOKWALTER Retractor system, provides multidirectional exposure of the operative site. The arms of the adjustable frame can be positioned independently, as desired.

The frame is provided nonsterile and must be sterilized before use, including initial use.

Precautions

If using monopolar electrosurgical equipment while using the adjustable frame and other BOOKWALTER retractor components, consult the basic rules of usage concerning electrosurgery (see EN 60601-2-2), as well as the manufacturer's instructions for use regarding specific warnings and conditions.

Be sure to electrically isolate the patient from high frequency leakage currents to the operating table, in accordance with the instructions provided by the manufacturer of the electrosurgical equipment.

Ensure that the patient does not come into contact with metal parts that are earthed or that have an appreciable capacitance to earth (for example: operating table supports), in accordance with the instructions provided by the manufacturer of the electrosurgical equipment.

Avoid small-area contacts with objects having a low impedance to earth at high frequencies, as unwanted burns can result.

Using the Adjustable Frame

The following BOOKWALTER components are also required when using the adjustable frame. Read the instructions for use booklet that accompanies each component.

Component	Reorder Number
1 Table Post	50-4581
1 Post Coupling	50-4554
BOOKWALTER Retractor Blades	Appropriate number and configuration procedure requires
Ratchet Mechanisms	50-4579, 50-4580, 50-4589

1. With the table post secured to the rail of the O.R. table, insert the support arm of the adjustable frame into the post coupling and secure the frame over the patient (Figure 1). Lock the arms in position by pushing the securing handle back (Figure 2a).
2. **NOTE:** For complete closure of the adjustable frame, position the WISHBOOK bridge (50-4687) to span the opening of the adjustable arms (Figure 2b), and lock in position by pushing the securing handle back as instructed in Step 1.
3. To readjust the frame, unlock the securing handle (Figure 3).
4. Assemble the retractor blades and ratchet mechanisms. Position the ratchet mechanisms as desired on the arms. Retract the blades as needed.

Cleaning and Sterilization

Refer to Symmetry Surgical Instructions for Use #204233 for updated information regarding reprocessing the instrument.

See **CAUTIONS** below.

Copies of this instruction sheet are available by requesting IFU #204233 from Symmetry Surgical Customer Service at:
Telephone: 800 251 3000
Fax: +1 615 964 5566
E-mail: customerservice@symmetrysurgical.com

CAUTION: Do not lubricate the BOOKWALTER WISHBOOK Adjustable Frame. Lubricant might interfere with the device's ability to maintain position while under load.

Service and Repair

For service or repair outside the United States, contact your local Symmetry Surgical representative.

Inside the United States, send instruments for service or repair to:

Symmetry Surgical Inc.
3034 Owen Drive,
Antioch, TN 37013 USA

Always include a repair purchase order number and a written description of the problem.

Warranty

Symmetry Surgical warrants that this medical device is free from defects in both materials and workmanship. **Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed. Suitability for use of this medical device for any particular surgical procedure should be determined by the user in conformance with the manufacturer's instructions for use. There are no warranties that extend beyond the description on the face hereof.**

© BOOKWALTER is a registered trademark of Symmetry Surgical

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