



QBD1200 TOC Analyzer Field Service Partnership

Your Hach **QBD1200 TOC Analyzer** Field Service Partnership (FSP) provides all-inclusive parts and one scheduled preventative maintenance visit performed by a Hach Field Service Technician. The FSP also includes all visits authorized by the Hach Technical Support Team and a special priority toll free number that will be included with your Partnership documentation.

During the pre-scheduled site visits, your Hach Field Service Technician will complete:

Verification of Instrument Performance/Maintenance

- Perform limited instrument cleaning.
- Review and evaluate user-programmed parameters.
- Evaluate any instrument error or warning conditions.
- Verify instrument operating voltages.
- Annually replace ozone destructor, tubing, syringe and UV reactor.
- Calibrate to factory-recommended specifications.
- Verify software version and update as necessary.

Repairs

- Perform required repair service including parts and labor.
- Includes sending unit to the factory at no additional charge if unable to repair in the field. This instrument will go to the head of the repair queue.
- Abuse, neglect, and acts of God not covered.

Reporting/Certificate of Performance

- Provide Hach Field Service Report with complete documentation of service performed and measurements/readings.
- Issue Certificate of Instrument Performance for each instrument that successfully passes final testing.

Training

- Provide basic end-user training on general instrument operation and maintenance (advance notice required).