
TOPCON HARMONY®

Clinical Data Management tailored for you.

With customizable Clinical Views to help you make a faster diagnosis or treatment decision.

Age-related macular degeneration (AMD)

Harmony is a comprehensive data management and communication platform for eye care. Its customizable Clinical View layout is a powerful tool for eye care specialists to efficiently access, review, and document patients' diagnostic data at a glance. The healthcare professional is in complete control and can arrange the screen layout based on the disease entity, in this case, AMD, to save time and increase efficiency.



The timeline feature shows the examinations performed over time to provide a complete picture of the patient's journey

Magic Lens tool **enables analysis of multiple imaging modalities** in an isolated area of interest



TOPCON HARMONY®

Clinical Data Management tailored for you.


With customizable Clinical Views to help you make a faster diagnosis or treatment decision.

Glaucoma

Harmony provides eye care professionals with convenient navigation and analysis tools for parallel image review so they can easily monitor progression over time and quickly identify a significant change in complex diseases such as the Glaucoma case shown below.

The overlay tool allows comparison of examination items throughout the course of the patient's treatment

Side-by-side comparison of all available OCT and visual field reports



The Slider tool enables a dynamic 'fly through' option to visually inspect patterns of change

