



See all of the features that make Hygiena's ATP Cleaning Verification System the preferred choice of hospitals and healthcare facilities worldwide.

SystemSURE *Plus* + UltraSnap[®] ATP CLEANING VERIFICATION SYSTEM SPECIFICATIONS

Connects to PC via USB.

Pass, Caution, or Fail result symbols are easy to interpret for any user.

100% recyclable plastic.

Easy 3-step activation requires minimal training.

200 programmable user ID's to track staff performance.

Powered by 2 AA batteries for 8+ months of uninterrupted testing.

Online training videos and instructional guides available.

Organize test locations into plans (i.e., departments, floors, buildings, etc).

5,000 programmable test locations.

Self-checks calibration at startup.

Advanced sensor technology is durable and stays in calibration for years.

Lightweight, hand-held design is convenient for one hand use.

Patented device design and unique liquid-stable chemistry provides consistent, reliable results.

15 month shelf life at refrigerated temperatures, 4 weeks at room temperature.

Pre-moistened swab bud breaks through biofilm for maximum sample collection.

Catalog No:
SS3H SystemSURE Plus & SureTrend Software
US2020 UltraSnap Surface ATP Swabs (100)

RLU
0

Prog:1
Bed Rail
Plan:1
Patient Room
13:58 2/ 4/2013

SystemSURE *Plus*

hygiena
UltraSnap[®]

WWW.HYGIENA.COM

INFO@HYGIENA.COM



- ✓ **TRACK TEST RESULTS**
- ✓ **QUICKLY IDENTIFY PROBLEM AREAS**
- ✓ **ANALYZE CLEANING EFFECTIVENESS**
- ✓ **AUTOMATICALLY GENERATE REPORTS**
- ✓ **FACILITATE CONTINUOUS IMPROVEMENT**

BENEFITS

- Free upgrades
- No licensing fees
- Unlimited technical support
- Fraud-proof data protection
- Localized data storage
- Password protection
- PC/Windows Software

HOSPITAL HYGIENE IN YOUR INBOX

SureTrend automatically e-mails charts and reports to as many recipients as you wish.

40+ PRE-SET REPORTS

Customizing charts and reports is only a few clicks away. With dozens of pre-set reports, managers can start analyzing data immediately.

GLOBAL LOCATION DIRECTORY

Infection Control and EVS can collaborate to create a master list of locations to be tested in a hospital system and share that list with any user in their network.

NETWORK YOUR GROUP

SureTrend may be installed on a network server to link devices, departments, or even hospitals.

MORE MEANINGFUL HUDDLES

Use the charts and graphs generated by SureTrend as visual aids to illustrate a problem, motivate team members, and work towards a common goal.

TEST PLAN RANDOMIZATION

Randomization offers truly statistical coverage for facilities that wish to monitor 100% of locations over time without testing every location daily.