



## 15211 Film Anti-Fog Surgical Mask

### Intended Use

Intended to be used in environments where low fluid exposure is a risk.



### Complies with ASTM F2100

Fluid Resistance (mmHg): 80  
BFE (%) @ 3.0  $\mu\text{m}$ :  $\geq 95\%$   
PFE (%) @ 0.1  $\mu\text{m}$ :  $\geq 95\%$   
 $\Delta P$  (mm H<sub>2</sub>O/cm<sup>2</sup>):  $< 4$   
Flame Spread: Class 1

### Key Features

- Complies with ASTM F2100 Level 1 requirements
- Anti-fog film inhibits eyewear fogging improving vision
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft longer ties promote secure comfortable fit
- Made in North America
- Not made with natural rubber latex

# Film Anti-Fog Surgical Mask

Intended to be used in environments where low fluid exposure is a risk.



## Technical Specifications

|                   |                             |
|-------------------|-----------------------------|
| Part No.          | 15211                       |
| Description       | Film Anti-Fog Surgical Mask |
| Color             | Green                       |
| Fluid Resistance  | 80 mm Hg                    |
| Filtration – BFE  | ≥99.7                       |
| Filtration – PFE  | ≥99.7                       |
| Delta P           | ≤2.9                        |
| Flammability      | Class 1                     |
| Attachment Design | Vertical Ties               |
| Mask Construction | Sonic Welded                |
| Qty               | 50/Box, 6/Boxes, 300/Case   |

Precept® is a registered trademark of Aspen Surgical Products, Inc.

Aspen Surgical reserves the right to make changes without notice in design, specifications and model.

©2021 Aspen Surgical Products, Inc. ALL RIGHTS RESERVED.  
M-SM-CAT-1051525 rev 1 29-JUN-2021 ENG – US

**Aspen Surgical Products, Inc.**  
phone +1 616.698.7100 toll-free +1 888.364.7004  
[aspensurgical.com](http://aspensurgical.com)

