



Vein Strippers

Symmetry Surgical single-use **Vein Strippers** are designed to provide proven, cost-effective removal of saphenous varicose veins, while eliminating the need to clean and process the vein stripper after each use.

Symmetry Surgical Vein Strippers:

- Sterile Package Includes:
 - Handle (1)
 - 1m Cables (2)
 - 9.5mm Olive (1)
 - 12.85mm Olive (1)
 - 15.4mm Olive (1)
- Latex free
- Radiopaque



CODE	DESCRIPTION	BOX QUANTITY
VE-022	Dermo-Strip Vein Stripper	12 Packs

Contact your Specialist or [visit **symmetrysurgical.com**](http://www.symmetrysurgical.com) for more information.