


# OPTIMAL CONTROL BETTER PRECISION


## BOVIE® MICRODISSECTION NEEDLES

*Delivers precision and accuracy in electrosurgery for a multitude of plastic, ENT, head, neck, oral and maxillofacial procedures*

- Highest grade tungsten wire needle
- Ultra sharp 5 micron tip
- Broad selection of microdissection needles



# MICRODISSECTION NEEDLES



CALL YOUR SYMMETRY SURGICAL SALES REPRESENTATIVE, LEARN MORE AT [SYMMETRYSURGICAL.COM](http://SYMMETRYSURGICAL.COM) OR CALL 1-800-251-3000.

Needle Electrode, Tungsten Wire Microdissection, Straight, Ultra Sharp, Single Use 5/box

Code	Length (L) in	Length (L) cm
97130BX	1 1/4in	3.2cm
97131BX	1 5/8in	4.1cm
97132BX	2in	5.1cm
97133BX	2 1/2in	6.4cm
97134BX	4 5/8in	11.7cm
97135BX	5 1/2 in	14cm
97136BX	6 5/8in	16.8cm

Needle Electrode, Tungsten Wire Microdissection, Angled Tip, Ultra Sharp Single Use 5/box

Code	Length (L) in	Length (L) cm
97141BX	1 5/8in	4.1cm
97143BX	2in	5.1cm
97144BX	4in	10.2cm
97145BX	4 5/8in	11.7cm
97147BX	6 5/8in	16.8cm
97148BX	7in	17.8cm

Needle Electrode, Tungsten Wire Microdissection, Angled Shaft, Ultra Sharp, Single Use 5/box

Code	Length (L) in	Length (L) cm
97125BX	6in	15.2cm
97150BX	2in	5.1cm
97152BX	4 7/8in	12.4cm
97154BX	6 7/8in	17.4cm

Needle Electrode, Tungsten Wire Microdissection, Curved Shaft, Ultra Sharp, Single Use 5/box

Code	Length (L) in	Length (L) cm
97170BX	4in	10.2cm
97171BX	5 7/8in	14.9cm

Needle Electrode, Tungsten Wire Microdissection, Straight, SuperCut Single Use 5/box

Code	Length (L) cm
ES60	3cm
ES61	3cm
ES61HS	3cm
ES62	4.5cm
ES63	4.5cm

Needle Electrode, Tungsten Wire Microdissection, Straight, SuperCut Single Use 5/box

Code	Length (L) in	Length (L) cm
97110BX	1 1/2in	3.8cm
97111BX	2in	5.1cm
97112BX	4 7/8in	12.4cm
97116BX	3 1/2in	8.9cm

Needle Electrode, Tungsten Wire Microdissection, Straight, SuperCut Reusable

Code	Length (L) cm
ES60R	3cm
ES62R	4.5cm