


## Introducing Reddick® Cholangiogram Catheter


- Scoop tip curved introducer and balloon design for stable catheter placement – without clamping
- Depth marks to visualize catheter location
- Dual lumen, non-radiopaque catheter
- Needle introducer option

*Simplifies dye-enhanced imaging for more complete stone detection and anatomy identification for gallbladder removal*



## Reddick® Cholangiogram Catheter

# Reddick® Scoop Tip Cholangiogram Catheter



## Introducer Options:

Scoop Tip Introducer and Needle Introducer



## Complementary Laparoscopic Instrument

Reddick-Saye Screw Retractor Kit

Ideal when grasping instruments fail. Screws into distended gall bladder or cystic ovary wall for positive fixation and excellent retraction.



Visit [symmetrysurgical.com](http://symmetrysurgical.com) or contact your local Sales Representative for more information.

### Reddick® Scoop Tip Cholangiogram Catheter

CODE	DESCRIPTION	DIAMETER	LENGTH
e2401-50	Reddick® Cholangiogram Catheter	4F	50 cm
e2401-51	Reddick® Cholangiogram Stiffer Catheter	4F	50 cm
e2401-52	Reddick® Cholangiogram Catheter with Needle Introducer (Needle Introducer: 13 gauge by 8 inches)	4F	50 cm
e2401-53	Reddick® Cholangiogram Stiffer Catheter with Needle Introducer (Needle Introducer: 13 gauge by 8 inches)	4F	50 cm

### Reddick® Cholangiogram Catheter

CODE	DESCRIPTION	DIAMETER	LENGTH
e2400-50	Reddick® Cholangiogram Catheter	4F	50 cm
e2400-51	Reddick® Cholangiogram Stiffer Catheter	4F	50 cm
e2400-52	Reddick® Cholangiogram Catheter with Needle Introducer (Needle Introducer: 13 gauge by 8 inches)	4F	50 cm
e2400-53	Reddick® Cholangiogram Stiffer Catheter with Needle Introducer (Needle Introducer: 13 gauge by 8 inches)	4F	50 cm

### Laparoscopic Instruments (4 per box)

CODE	DESCRIPTION	DIAMETER	LENGTH
e2504-01	Reddick-Saye Screw Retractor Kit	5 mm	30 cm

1 Reddick®-Saye Screw Retractor Kit and 1 13G Aspiration Needle

