



- **SIX INCH FLEXIBLE SHAFT**
- **8, 10 OR 12 FRENCH DIAMETERS**
- **FOOT OR HAND ACTIVATION**
- **FLUID EVACUATION AND/OR COAGULATION**

Continuing to expand its electrosurgical offering, Bovie Medical Corporation is pleased to introduce six Suction Coagulators. They are compatible with Bovie® and Aaron® generators, as well as, those of most other companies. Featuring three different diameters, flexible shaft and hand or foot activation, these suction coagulators allow suction with or without simultaneous coagulation. These single use devices are individually packaged sterile, ten per shelf pack.

HANDSWITCH SUCTION COAGULATORS

Bovie® Catalog Number	Description	Tube Diameter (French Size)	Cable Length
SCH08	Disposable	08 Fr	3m
SCH10	Disposable	10 Fr	3m
SCH12	Disposable	12 Fr	3m

FOOTSWITCH SUCTION COAGULATORS

Bovie® Catalog Number	Description	Tube Diameter (French Size)	Cable Length
SCF08	Disposable	08 Fr	3m
SCF10	Disposable	10 Fr	3m
SCF12	Disposable	12 Fr	3m

SUCTION COAGULATORS CONVERSION CHART

Bovie®	Valley Lab™	ConMed®	Megadyne®
SCH08	E2608-6		0042-10
SCH10	E2610-6	130188	0041-10
SCH12	E2612-6		0043-10
SCF08			
SCF10	E2505-10FR	130187	0040-10
SCF12			



	●	●	●
French Gauge	12	10	8
Diameter (mm)	4	3.3	2.7
Diameter (inches)	0.158	0.131	0.105



Bovie is a registered trademark of Bovie Medical Corporation.
 Valleylab is a trademark of Covidien.
 ConMed is a registered trademark of ConMed Electrosurgery.
 Megadyne is a registered trademark of Megadyne Medical Products, Inc.