




Introducing the Bovie® SharkSkin

Ergonomic smoke evacuation pencil adaptor


- Transforms electrosurgical pencil into a device that captures smoke plume at the surgical site
- Eliminates the need to hold smoke evacuation tubing near the surgical site
- Attaches quickly and easily
- Does not significantly increase the overall diameter of the assembled handpiece
- Compatible with many electrosurgical pencils



(Electrosurgical pencil not included)

Bovie® SharkSkin

Bovie® SharkSkin



Shown with ESPR2 electrocautery pencil (not included)

Designed to accommodate the reusable handpieces of the Bovie® DERM 101, 102, A942 and Bantam|PRO (A952) high frequency desiccators

Molded plastic housing with lightweight 3/8" tubing provides cost effective option for surgical smoke evacuation

Transparent Capture Port improves visibility of the surgical site

Visit boviemedical.com or contact your local Sales Representative for more information.

Bovie® SharkSkin

CODE	DESCRIPTION	QUANTITY
SS95	SharkSkin Electrocautery Pencil Adapter With 3/8" x 10' Tubing (packaged sterile individually)	Box of 20

