


...FROM A NAME YOU CAN TRUST

All Aaron® blade, ball, needle, loop, and square disposable electrodes use the standard 3/32" diameter stainless steel shafts. These electrodes, like all Aaron® electrodes, feature safety grip insulators, combining patient and user safety with easy insertion into and removal from the surgical pencil.




Many Electrode Sizes Available



The loop and square electrodes feature tungsten wire for superior shape and wire integrity throughout the excision procedure, while minimizing thermal damage and yielding excellent histopathology specimens. Each molded loop and square electrode includes our patented tray to protect the loops from bending during shipping and handling.

Ball electrodes in 3mm and 5mm sizes are available for fulguration and desiccation, while needle electrodes are available for precise, pinpoint cutting and coagulation.



List of Electrode Styles and Conversion Chart...

| | | | STRYKER/ COLORADO | PREMIER | WALLACH | CONMED | VALLEY | ACMI/ CIRCON | UTAH | MILTEX |
|---|-------------------------|---|----------------------|-----------|-------------|---------------|---------|-----------------|---------|---------|
| SuperCut™ Tungsten Needles | ES60 | 3CM Super Fine Needle, 5/box | N101A | | | | | | | |
| | ES61 | 3CM Modified Super Fine Needle, 5/box | N103A | | | | E1651 | | | |
| | ES62 | 4.5CM Super Fine Needle, 5/box | N102A | | | | | | | |
| | ES63 | 4.5CM Modified Super Fine Needle, 5/box | N104A | | | | | | | |
| | | | COOPER | PREMIER | WALLACH | CONMED | VALLEY | ACMI/ CIRCON | UTAH | MILTEX |
| Needle Electrodes | ES02 | Standard Needle, 25/box | | | 909135 | 275-1502/9113 | E1552 | | | |
| | ES03 | Extended Needle, 25/box | | | | 275-1504/9116 | E1552-6 | | | |
| | ES38 | Modified Needle, 25/box | | | | | | | | |
| | ES40 | Extended Modified Needle, 25/box | N2512 | 9006240 | 909005 | 138025 | | 6162901 | | |
| Dermal Loops | ES22-8 | 1/4" Short Shaft Loop, 5/box | | 9006406 | | 138011 | | | | |
| | ES23-8 | 5/16" Short Shaft Loop, 5/box | | | | | | | | |
| | ES24-8 | 3/8" Short Shaft Loop, 5/box | | 9006407 | | 138010 | | | | |
| | ES25-8 | 1/2" Short Shaft Loop, 5/box | | 9006400 | | | | | | |
| | ES20-8 | 3/4" Short Shaft Loop, 5/box | | | | 138012 | | | | |
| Dermal Tips | A804 | Sharp Dermal Tip, Non-Sterile, 100/box | | | | | | | | |
| | A805 | Sharp Dermal Tip, Sterile, 50/box | | | | | | | | |
| | A806 | Blunt Dermal Tip, Non-Sterile, 100/box | | | | | | | | |
| | A807 | Blunt Dermal Tip, Sterile, 50/box | | | | | | | | |
| Blade Electrodes | ES01 | Standard Blade, 50/box | | | | 275-1500/9111 | E1551G | | | |
| | ES04 | Extended Blade, 25/box | | | | 275-1503/9114 | E1551-6 | | | |
| | ES37 | Modified Blade, 25/box | | | | 275-1524 | E1553 | | | |
| | ES39 | Extended Modified Blade 25/box | | | | 275-1525 | | | | |
| | ES16 | Angled Blade, 5/box | | | | | | | | |
| Ball Electrodes | ES06 | Extended 3mm Ball, 5/box | B0312 | 9006225 | 909038 | 60-5181-102 | E1563 | 6156901 | DBL-311 | 31-5000 |
| | ES07 | Extended 5mm Ball, 5/box | B0512 | 9006230 | 909003 | 60-5181-103 | E1564 | 6157901 | DBL-511 | 31-5002 |
| | ES20 | 3/16" Ball, 25/box | | | | 138103 | E1550 | | | |
| | ES21 | 5/32" Ball, 25/box | | | 909136 | | | | | |
| | ES60 | 1/8" Ball, 25/box | | | | | | | | |
| LLETZ Loop Electrodes | ES08 | 5mm x 5mm Loop, 5/box | | | | | | 6133901 | DLP-T11 | 31-3056 |
| | ES09 | 10mm x 10mm Loop, 5/box | R1010 | 9006200 | 909007 | | E1559 | 6134901 | DLP-S11 | 31-5072 |
| | ES10 | 15mm x 8mm Loop, 5/box | R1507 * | 9006205 | 909011 | | | 6135901 | | |
| | ES11 | 15mm x 10mm Loop, 5/box | | | | 60-5180-102 | | 6136901 | DLP-M11 | 31-5074 |
| | ES12 | 20mm x 8mm Loop, 5/box | R2008 | 9006215 | 909013 | | | 6137901 | | 31-5090 |
| | ES13 | 20mm x 15mm Loop, 5/box | R2008 | 9006215 | 909009 | | | 6137901 | | 31-5096 |
| | ES31 | 20mm x 20mm Loop, 5/box | R2020 | | | | | | | 31-5098 |
| | ES42 | 20mm x 12mm Loop, 5/box | | 9006220 * | | | E1561 | | DLP-W11 | 31-5094 |
| | ES43 | 15mm x 15mm Loop, 5/box | R1515 | 9006210 | | | | | | |
| | ES44 | 15mm x 5mm Loop, 5/box | R1505 | | 909134 | | | | | |
| | ES45 | 13mm x 13mm Loop, 5/box | | | | | | | | |
| | ES46 | 10mm x 5mm Loop, 5/box | | | | | | | | |
| | ES47 | 25mm x 10mm Loop, 5/box | R2510 | | | | | | | |
| | ES48 | 10mm x 3mm Loop, 5/box | | | | | | | | |
| ES49 | 20mm x 10mm Loop, 5/box | R2010 | | 909132 | 60-5180-104 | | | | | |

| | | | COOPER | PREMIER | WALLACH | CONMED | VALLEY | ACMI/CIRCON | UTAH | MILTEX |
|--------------------------------|------|--------------------------------|--------|---------|----------|---------------|--------|-------------|---------|---------|
| LLETZ Loop Electrodes | ES51 | 13mm x 8mm Loop, 5/box | | | | 60-5180-101 * | | | | |
| | ES52 | 10mm x 8mm Loop, 5/box | | | 909133* | | | | | 31-5054 |
| | ES53 | 15mm x 12mm Loop, 5/box | | | | | E1560 | | | |
| LLETZ Square Electrodes | ES14 | 5mm x 5 mm Square, 5/box | | | | | E1565 | 6153901 | | |
| | ES15 | 10mm x 4mm Square, 5/box | S1004 | 9006235 | | | E1566 | 6154901 | DLP-SQ2 | |
| | ES16 | 10mm x 8mm Square, 5/box | | | 909017 | | E1567 | 6155901 | DLP-SQ1 | |
| | ES41 | 5mm x 10mm Square, 5/box | | | | | | | | |
| Arthroscopic Electrodes | AR00 | Electrode Hook 90, 5/box | | | 88-1510 | | E1510 | | | |
| | AR01 | Electrode Hook 45, 5/box | | | 88-1511 | | E1511 | | | |
| | AR02 | Electrode Acromioplasty, 5/box | | | 88-1512 | CS012A | E1512 | | | |
| | AR03 | Electrode Meniscectomy, 5/box | | | 88-15103 | CS010A | E1513 | | | |

* shows 1mm difference in size

ARTHROSCOPIC ELECTRODES

The Aaron Arthroscopic electrodes are available in four tip styles and designed to fit all standard electro-surgical pencils.



AR00*



AR01*



AR02*



AR03*



ES63*


ES62*

ES61*

ES60*

ES62R

ES60R



SUPERCUT™ NEEDLES

SuperCut™ Needles (ES 60-63) are made of Tungsten and keep their shape throughout the procedure.

*STERILE: Individually packaged, sterile

Short Tungsten Loops...



ES22-8

ES23-8

ES24-8

ES25-8

ES26-8

DISPOSABLE SHORT TUNGSTEN .008 WIRE LOOPS

ES22-8 — ¼" short shaft loop, sterile, .008" wire - 5/box

ES23-8 — ⅝" short shaft loop, sterile, .008" wire - 5/box

ES24-8 — ⅜" short shaft loop, sterile, .008" wire - 5/box

ES25-8 — ½" short shaft loop, sterile, .008" wire - 5/box

ES26-8 — ¾" short shaft loop, sterile, .008" wire - 5/box



ES22-R

ES23-R

ES24-R

ES25-R

ES26-R

REUSABLE SHORT TUNGSTEN LOOPS

ES22R — ¼" short shaft loop, non-sterile - 1/box

ES23R — ⅝" short shaft loop, non-sterile - 1/box

ES24R — ⅜" short shaft loop, non-sterile - 1/box

ES25R — ½" short shaft loop, non-sterile - 1/box

ES26R — ¾" short shaft loop, non-sterile - 1/box

Distributed By:



3034 Owen Drive • Antioch, TN 37013
U.S. Phone 1-800-251-3000 • Int'l Phone +1-615-964-5532
www.symmetrysurgical.com • customerservice@symmetrysurgical.com