

Bard-Parker® Protected Blades and Scalpels

Protect Yourself.
Protect Your Staff.

Choose from blade systems with reusable stainless steel handles or fully disposable scalpels.

Designed with a similar weight and feel as conventional products, while using the same high-quality stainless steel blade.

Help provide protection from injuries that may occur before use, during steps of a procedure, or after use.

Features a locking retractable shield that boasts easy one-handed activation.

Protected blades are quick and easy to load and unload.


Bard-Parker® Conventional Blades and Scalpels

In the Operating Room, Emergency Room, Labor & Delivery and other Healthcare Departments, Quality Matters.

High quality blades, including the exclusive Carbon Rib-Back® design, allow for enhanced strength and cutting control.

Our blades are packaged in a foil pouch with a unique inner liner for maximum protection of blade and for the end user.

Choose from blade systems with reusable handles or fully disposable scalpels.



Bard-Parker® and Rib-Back® are registered trademarks of Aspen Surgical Products, Inc. SafetyLock™ is a trademark of Aspen Surgical Products, Inc. Aspen Surgical reserves the right to make changes without notice in design, specifications and models.


©2020 Aspen Surgical Products, Inc. ALL RIGHTS RESERVED.
M-SM-CAT-1005176 rev 7 15-AUG-2020 ENG - US

For further information about this product or a service, please contact your local Aspen Surgical representative or visit our webpage:

Phone 616-698-7100

Toll-free 888-364-7004

www.aspensurgical.com



Bard-Parker®

Surgical Blades, Handles and Disposable Scalpels

For over 100 years, the surgeons' choice for a high level of quality, sharpness, and precision.

 aspen surgical®

The Bard-Parker® Brand

Surgery Requires Precision.
Demand Bard-Parker® Brand Products.


For over 100 years, the Bard-Parker® brand has remained a leader in surgical blades and scalpels by delivering constant sharpness and a high level of quality.

Bard-Parker® blades and scalpels are proudly made in the United States to ensure tight control over product quality and lead-times.

Choose from a wide selection of surgical blades, handles and disposable scalpels, including both conventional and protected styles.

Our manufacturing capacity provides supply reliability. With this our Bard-Parker® blades and scalpels are manufactured with a proprietary steel composition and process to provide high-performing surgical instruments.

Contact our committed team of dedicated sales consultants to help evaluate staff needs and find the appropriate product fit.




Disposable Scalpels

Stainless Steel Blade

Surgical Blades

Carbon and Stainless Steel Blades

Surgical Blade Handles



Sterile Conventional	Sterile Protected	Sterile Protected (Long Handle)	Non-Sterile Conventional	Non-Sterile Protected	Blade Image	Size	Sterile Carbon Steel Rib-Back®	Sterile Carbon Steel SafetyLock™ Rib-Back®	Sterile Stainless Steel Conventional	Sterile Stainless Steel Protected	Non-Sterile Carbon Steel Rib-Back®
371610	372610	372710	371630	373810		10	371110	371150	371210	373910	371310
371611	372611	372711	371631	373811		11	371111	371151	371211	373911	371311
371612	372612					12	371112	371152	371212		371312
371615	372615	372715	371633	373815		15	371115	371153	371215	373915	371315
371620	372620		371634	373820		20	371120	371154	371220	373920	371320
371621	372621		371635	373821		21	371121	371155	371221	373921	371321
	372622			373822		22	371122	371156	371222	373922	371322
371623	372623			373823		23	371123	371157	371223	373923	371323
						24					371324
						25					371325
						60					371340
						10A			371700		
371613	372613					12B			371712		
						15C			371716	374305	
						11					

371641*

Conventional Handles – Non-Sterile

Reorder No.	Handle Size	Handle Image	Blade Sizes
371030	3		10-15
371031	3L		10-15
371040	4		20-25
371050	5		10-15
371060	6		20-25
371070	7		10-15
371080	8		60
371090	9		10-15

Protected Handles – Non-Sterile

374030	3		10-15
374031	3L		10-15
374040	4		20-23
374041	4L		20-23

Disposable Scalpel Handle & Cover – Sterile

ST-935D	3		10-15
ST-830D	3		10-15C

10 per box, 10 boxes per case
 * Mini Scalpels 100 per bag, 10 bags per case
 Protected Non-Sterile Scalpels: 10 per Strip, 50 strips per case

Individually packaged, 50 blades per box, 3 boxes per case
 Non-Sterile, Carbon Steel, Rib-Back®: Individually packaged,
 25 strips of 6 blades, 150 blades per case

5 per case

