

Colpo-Master™ II

C205T-LED Camera-Ready
Trinocular head
with camera



COLPO-MASTER™ II CENTER-POST MOUNT COLPOSCOPES

Colpo-Master™ II is a compact Center-Post Colposcope with unique, continuous zoom magnification and sharp focusing for crisp 3-D imaging.

Features:

- **HEAD:** Full continuous zoom magnification range: 4 – 27x with 20x eyepieces. Binocular or video-ready Trinocular heads, with 298mm working distance. The field of view (FOV) ranges from 58.5mm at 4x to 8.5mm at 27x
- 3-D Stereoscopic Viewing
- Dual diopters allow user to independently focus each eyepiece to adjust for the natural vision differential between eyes for viewing comfort. Camera-ready Trinocular models allow the user to view through both eyepieces while using the camera.
- Sturdy, 5-legged, easy-to-maneuver base slides in and under exam tables for optimal positioning
- Dual light intensity controls for both white and green LEDs
- Adjustable height center post with smooth suspension system
- S-Video and digital camera options available, trinocular version with built-in C-Mount Camera adapter

Complete with dust cover, power cord, manual, tool kit and 7 Year Limited Warranty, 1 year on electrical system (excludes bulbs and fuses). Optional accessories are available - please contact your Bovie Representative for additional information.

Item No.	Description
CS-205LED	Colpo-Master™ II Binocular Zoom Center Post Colposcope
CS-205LED-E	Colpo-Master™ II Binocular Zoom Center Post Colposcope, 220V
CS-205TLED	Colpo-Master™ II Camera-Ready Zoom Center Post Colposcope
CS-205TLED-E	Colpo-Master™ II Camera-Ready Zoom Center Post Colposcope, 220V

Colpo-Master™ is a registered trademark of United Products and Instruments, Inc; DBA UNICO